

INTOXILYZER 8000 INSPECTION

Serial No.: 80-003066

Location: WARD CO

- | | Check When Done |
|---|-----------------|
| A. Is the warm-up time less than 20 minutes? | <u>✓</u> |
| B. Is a three-pronged grounded outlet used? | <u>✓</u> |
| C. Is the breath tube heated? | <u>✓</u> |
| D. Is the diagnostic check complete? | <u>✓</u> |
| E. Is the time, date, and year correct (re-set if necessary)? | <u>✓</u> |
| F. Print test. (Attach test record.) | <u>✓</u> |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)
Sim SN <u>DR5132</u> Lot No. <u> </u> Rep AC <u>0.021</u> | <u>✓</u> |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN <u>DR5134</u> Lot No. <u> </u> Rep AC <u>0.294</u> | <u>✓</u> |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.
(Attach test record.)
Sim SN <u>DR3846</u> Lot No. <u> </u> Rep AC <u>0.10AC + acetone</u> | <u>✓</u> |
| J. RFI check. Run CMS mode. Key radio on first room air.
(Attach test record.) | <u>✓</u> |
| K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.) | |

Lot No. <u>659358</u>	Cyl No. <u>32</u>	Exp. Date <u>09/03/2011</u>
Test 1 <u>0.081</u>	Test 4 <u>0.081</u>	Test 7 <u>0.081</u>
Test 2 <u>0.081</u>	Test 5 <u>0.081</u>	Test 8 <u>0.082</u>
Test 3 <u>0.080</u>	Test 6 <u>0.080</u>	Test 9 <u>0.082</u>
Average <u>0.081</u>		

L. Tank pressure, Level 3 Func D Sub G P: 283 psi Reg 270 psi

M. Remarks/Maintenance record: Annual inspection

03 Apr 2011
Date

Bob Kashur
Field Inspector's Signature

4-14-11
Date

[Signature]
Reviewed by

80-003066
P. 129 of 137
FM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003066
Location = WARD 8164.13.00 06/09
04/03/2011 12:35

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:35
02 Std. Sol.	0.018	12:36
03 Room Air	0.000	12:37
04 Std. Sol.	0.018	12:37
05 Room Air	0.000	12:38
06 Std. Sol.	0.018	12:39
07 Room Air	0.000	12:39

08 Sim Temp = 34.0°C

Simul Ser No = 0

Std Sol No = 0

County = 51

Oper No. = 777777

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

low AC 0.021AC

Form 106-I8000

80-003066
P. 131 of 137
JKM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003066
Location = WARD 8164.13.00 06/09
04/03/2011 12:57

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:58
02 Std. Sol.	0.291	12:58
03 Room Air	0.000	12:59
04 Std. Sol.	0.292	13:00
05 Room Air	0.000	13:00
06 Std. Sol.	0.293	13:01
07 Room Air	0.000	13:02

08 Sim Temp = 34.0°C

Simul Ser No = 0

Std Sol No = 0

County = 51

Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Linearity test 0.294 AC

Form 106-I8000

80-003066
P. 132 of 137
AL

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003066
Location = WARD 8164.13.00 06/09
04/03/2011 12:23

Test	AC	Time
01 Room Air	0.000	12:24
02 *Subject Test	INT*	12:24
03 Room Air	0.000	12:25

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 659358
Cyl No = 32
Expiration Date = 09/03/2011
County = 51 Oper No. = 777777

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

Interferent check 0.10AC + acetone

Form 106-I8000

80-003066
P. 133 of 137
JAM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003066
Location = WARD 8164.13.00 06/09
04/03/2011 12:26

Test	AC	Time
01 Diagnostic	OK	12:27
02 Room Air	RFI*	12:27
03 Room Air	0.000	12:28

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 659358
Cyl No = 32
Expiration Date = 09/03/2011
County = 51 Oper No. = 777777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
DEB KASHUR

Remarks:

RFI check

Form 106-I8000

80-003066
P.134 of 137
FM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003066
Location = WARD 8164.13.00 06/09
04/03/2011 12:09

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:10
02 Std. Gas	0.081	12:10
03 Room Air	0.000	12:11
04 Std. Gas	0.081	12:11
05 Room Air	0.000	12:12
06 Std. Gas	0.080	12:12
07 Room Air	0.000	12:12

Lot No = 659358
Cyl No = 32
Exp Date = 09/03/2011
County = ~~08-DK~~ 51 Oper No. = 777777

Deb Kashur
Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080AC

Form 106-I8000

80-003066
P.135 of 137
FM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003066
Location = WARD 8164.13.00 06/09
04/03/2011 12:13

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:13
02 Std. Gas	0.081	12:14
03 Room Air	0.000	12:14
04 Std. Gas	0.081	12:14
05 Room Air	0.000	12:15
06 Std. Gas	0.080	12:15
07 Room Air	0.000	12:16

Lot No = 659358
Cyl No = 32
Exp Date = 09/03/2011
County = 51 Oper No. = 777777

Deb Kasher
Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080 AC

Form 106-I8000

80-003066
P.136 of 137
AM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003066
Location = WARD 8164.13.00 06/09
04/03/2011 12:16

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:17
02 Std. Gas	0.081	12:17
03 Room Air	0.000	12:18
04 Std. Gas	0.082	12:18
05 Room Air	0.000	12:18
06 Std. Gas	0.082	12:19
07 Room Air	0.000	12:19

Lot No = 659358
Cyl No = 32
Exp Date = 09/03/2011
County = 51 Oper No. = 777777

Deb Kashur
Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080 AC

Form 106-I8000

80-003066
P. 137 of 137
RN