

INTOXILYZER 8000 INSPECTION

Serial No.: 80-004953

Location: TOXL

- | | Check When Done |
|--|-----------------|
| A. Is the warm-up time less than 20 minutes? | <u>✓</u> |
| B. Is a three-pronged grounded outlet used? | <u>✓</u> |
| C. Is the breath tube heated? | <u>✓</u> |
| D. Is the diagnostic check complete? | <u>✓</u> |
| E. Is the time, date, and year correct (re-set if necessary)? | <u>✓</u> |
| F. Print test. (Attach test record.) | <u>✓</u> |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)
Sim SN <u>DR5144</u> Lot No. <u>—</u> Rep AC <u>0.021</u> | <u>✓</u> |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN <u>DR3847</u> Lot No. <u>—</u> Rep AC <u>0.294</u> | <u>✓</u> |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.
(Attach test record.)
Sim SN <u>DR3452</u> Lot No. <u>—</u> Rep AC <u>0.10</u> | <u>✓</u> |
| J. RFI check. Run CMS mode. Key radio on first room air.
(Attach test record.) | <u>✓</u> |
| K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.) | |

Lot No. <u>15410080A1</u>	Cyl No. <u>12</u>	Exp. Date <u>7-1-12</u>
Test 1 <u>0.079</u>	Test 4 <u>0.079</u>	Test 7 <u>0.080</u>
Test 2 <u>0.080</u>	Test 5 <u>0.079</u>	Test 8 <u>0.079</u>
Test 3 <u>0.079</u>	Test 6 <u>0.080</u>	Test 9 <u>0.080</u>
Average <u>0.079</u>		

L. Tank pressure, Level 3 Func D Sub G P: 335 psi Reg 330 psi

M. Remarks/Maintenance record: Annual Inspection

<u>4-5-11</u> Date	<u>[Signature]</u> Field Inspector's Signature
<u>05 Apr 2011</u> Date	<u>Deb Kashur</u> Reviewed by

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004953
Location = TOXL 8164.13.00 06/09
04/05/2011 14:10

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? No
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 15410080A1
Standard Cyl #? ~~40~~ 12
Standard Expiration? 07/01/2012
Oper No? 555555

Flow Cal. Date: 01/19/2010
Slope 691
Intercept -628068

IR Calibration Date: 06/08/2010
 3um 9um

0th Coef(*100): -10343 -18245
1st Coef(*100): 249276 131289
2nd Coef(*100): 4014 1439
H2O adj(mg/l*10k): 697 655

***** Printer Test End *****



Operator Signature
ASHLEY AYD

Remarks: *Print test*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004953
Location = TOXL 8164.13.00 06/09
04/05/2011 14:11

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:12
02 Std. Sol.	0.021	14:12
03 Room Air	0.000	14:13
04 Std. Sol.	0.022	14:14
05 Room Air	0.000	14:14
06 Std. Sol.	0.020	14:15
07 Room Air	0.000	14:15


08 Sim Temp = 34.0°C

Simul Ser No = 0

Std Sol No = 0

County = 08

Oper No. = 555555



Operator Signature
ASHLEY AYD

Remarks: low AC 0.021 AC

Form 106-I8000

80-004953
P. 47 of 53
JAM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004953
Location = TOXL 8164.13.00 06/09
04/05/2011 14:16

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:17
02 Std. Sol.	0.301	14:18
03 Room Air	0.000	14:18
04 Std. Sol.	0.303	14:19
05 Room Air	0.000	14:19
06 Std. Sol.	0.303	14:20
07 Room Air	0.000	14:21


08 Sim Temp = 34.0°C

Simul Ser No = 0

Std Sol No = 0

County = 08

Oper No. = 555555



Operator Signature
ASHLEY AYD

Remarks: High AC 0.294 AC

Form 106-I8000

80-004953
p. 48 of 53
SM


Intoxilyzer Test Record and Checklist
 NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
 North Dakota Model 8000 SN 80-004953
 Location = TOXL 8164.13.00 06/09
 04/05/2011 14:22

Test	AC	Time
01 Room Air	0.000	14:23
02 *Subject Test	INT*	14:23
03 Room Air	0.000	14:24

*Invalid Test
 Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
 Sub DOB = 01/01/1982
 Sub Sex = Female Weight = 150
 Test = DUI Cit = NA
 Dr. Lic. = ND/DIS821456
 Lot No = 15410080A1
 Cyl No = 12
 Expiration Date = 07/01/2012
 County = 08 Oper No. = 555555



Operator Signature
 ASHLEY AYD

Remarks: *Interferent check*

Form 106-I8000

80-004953
 p. 49 of 53
8m

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004953
Location = TOXL 8164.13.00 06/09
04/05/2011 14:24

Test	AC	Time
01 Diagnostic	OK	14:25
02 Room Air	RFI*	14:25
03 Room Air	0.000	14:26

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 15410080A1
Cyl No = 12
Expiration Date = 07/01/2012
County = 08 Oper No. = 555555

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
ASHLEY AYD

Remarks: RFI

Form 106-I8000

80-004953
P. 50 of 53
vkm


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004953
Location = TOXL 8164.13.00 06/09
04/05/2011 14:27

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:27
02 Std. Gas	0.079	14:28
03 Room Air	0.000	14:28
04 Std. Gas	0.080	14:28
05 Room Air	0.000	14:29
06 Std. Gas	0.079	14:29
07 Room Air	0.000	14:30

Lot No = 15410080A1
Cyl No = 12
Exp Date = 07/01/2012
County = 08 Oper No. = 555555



Operator Signature
ASHLEY AYD

Remarks: Calibration Check 0.080 AC

Form 106-I8000

80-004953
P. 51453
YAM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004953
Location = TOXL 8164.13.00 06/09
04/05/2011 14:59

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:00
02 Std. Gas	0.079	15:00
03 Room Air	0.000	15:01
04 Std. Gas	0.079	15:01
05 Room Air	0.000	15:02
06 Std. Gas	0.080	15:02
07 Room Air	0.000	15:02

Lot No = 15410080A1
Cyl No = 12
Exp Date = 07/01/2012
County = 08 Oper No. = 555555



Operator Signature
ASHLEY AYD

Remarks: Calibration check 0.080 AC

Form 106-I8000

80-004953
P. 52 of 53
FM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004953
Location = TOXL 8164.13.00 06/09
04/05/2011 15:04

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:04
02 Std. Gas	0.080	15:05
03 Room Air	0.000	15:05
04 Std. Gas	0.079	15:06
05 Room Air	0.000	15:06
06 Std. Gas	0.080	15:06
07 Room Air	0.000	15:07

Lot No = 15410080A1
Cyl No = 12
Exp Date = 07/01/2012
County = 08 Oper No. = 555555



Operator Signature
ASHLEY AYD

Remarks: Calibration check 0.080 AC

Form 106-I8000

80-004953
p. 53 of 53
FM