

INTOXILYZER 8000 INSPECTION

Serial No.: 80-003069

Location: TOXL

- | | Check When Done |
|--|-----------------|
| A. Is the warm-up time less than 20 minutes? | <u>✓</u> |
| B. Is a three-pronged grounded outlet used? | <u>✓</u> |
| C. Is the breath tube heated? | <u>✓</u> |
| D. Is the diagnostic check complete? | <u>✓</u> |
| E. Is the time, date, and year correct (re-set if necessary)? | <u>✓</u> |
| F. Print test. (Attach test record.) | <u>✓</u> |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)
Sim SN <u>DR5144</u> Lot No. <u>—</u> Rep AC <u>0.021</u> | <u>✓</u> |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN <u>DR3847</u> Lot No. <u>—</u> Rep AC <u>0.294</u> | <u>✓</u> |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.
(Attach test record.)
Sim SN <u>DR3452</u> Lot No. <u>—</u> Rep AC <u>0.10</u> | <u>✓</u> |
| J. RFI check. Run CMS mode. Key radio on first room air.
(Attach test record.) | <u>✓</u> |
| K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.) | |

Lot No. <u>15410080A1</u>	Cyl No. <u>12</u>	Exp. Date <u>7-1-12</u>
Test 1 <u>0.078</u>	Test 4 <u>0.078</u>	Test 7 <u>0.079</u>
Test 2 <u>0.077</u>	Test 5 <u>0.078</u>	Test 8 <u>0.079</u>
Test 3 <u>0.079</u>	Test 6 <u>0.079</u>	Test 9 <u>0.079</u>
Average <u>0.078</u>		

L. Tank pressure, Level 3 Func D Sub G P: 303 psi Reg 300 psi

M. Remarks/Maintenance record: Annual Inspection

4-6-11
Date

Deb Kashur
Field Inspector's Signature

06 Apr 11
Date

Deb Kashur
Reviewed by

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003069
Location = TOXL 8164.13.00 06/09
04/06/2011 08:24

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 15410080A1
Standard Cyl #? 12
Standard Expiration? 07/01/2012
Oper No? 555555

Flow Cal. Date: 04/08/2008
Slope 691
Intercept -599264

IR Calibration Date: 06/04/2008
 3um 9um

0th Coef(*100): -14879 -13696
1st Coef(*100): 280678 136270
2nd Coef(*100): 177 615
H2O adj(mg/l*10k): 391 185

***** Printer Test End *****



Operator Signature
ASHLEY AYD

Remarks: *Print test*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003069
Location = TOXL 8164.13.00 06/09
04/06/2011 08:29

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	08:30
02 Std. Sol.	0.019	08:30
03 Room Air	0.000	08:31
04 Std. Sol.	0.019	08:32
05 Room Air	0.000	08:32
06 Std. Sol.	0.020	08:33
07 Room Air	0.000	08:33


08 Sim Temp = 34.0°C

Simul Ser No = 0

Std Sol No = 0

County = 08

Oper No. = 555555



Operator Signature
ASHLEY AYD

Remarks: Low AC 0.021 AC

Form 106-I8000

80-003069
p. 98 of 104
KTM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003069
Location = TOXL 8164.13.00 06/09
04/06/2011 08:35

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	08:35
02 Std. Sol.	0.295	08:36
03 Room Air	0.000	08:36
04 Std. Sol.	0.297	08:37
05 Room Air	0.000	08:38
06 Std. Sol.	0.296	08:38
07 Room Air	0.000	08:39


08 Sim Temp = 34.0°C

Simul Ser No = 0

Std Sol No = 0

County = 08

Oper No. = 555555



Operator Signature
ASHLEY AYD

Remarks: High AC 0.294 AC

Form 106-I8000

80-003069
P. 99 of 104
KTM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003069
Location = TOXL 8164.13.00 06/09
04/06/2011 08:41

Test	AC	Time
01 Room Air	0.000	08:42
02 *Subject Test	INT*	08:42
03 Room Air	0.000	08:43

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 15410080A1
Cyl No = 12
Expiration Date = 07/01/2012
County = 08 Oper No. = 555555



Operator Signature
ASHLEY AYD

Remarks: *Interferent 0.10 AC + Acetone*

Form 106-I8000

80-003069
P. 100 of 104
RAM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003069
Location = TOXL 8164.13.00 06/09
04/06/2011 08:58

Test	AC	Time
01 Diagnostic	OK	08:59
02 Room Air	RFI*	08:59
03 Room Air	0.000	08:59

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 15410080A1
Cyl No = 12
Expiration Date = 07/01/2012
County = 08 Oper No. = 555555

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
ASHLEY AYD

Remarks: RFI

Form 106-I8000

80-003069
P. 101 of 104
vkm


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003069
Location = TOXL 8164.13.00 06/09
04/06/2011 09:00

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:00
02 Std. Gas	0.078	09:01
03 Room Air	0.000	09:01
04 Std. Gas	0.077	09:02
05 Room Air	0.000	09:02
06 Std. Gas	0.079	09:02
07 Room Air	0.000	09:03

Lot No = 15410080A1
Cyl No = 12
Exp Date = 07/01/2012
County = 08 Oper No. = 555555



Operator Signature
ASHLEY AYD

Remarks: Calibration check 0.080AC

Form 106-I8000

80-003069
P. 102 of 104
VM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003069
Location = TOXL 8164.13.00 06/09
04/06/2011 09:03

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:04
02 Std. Gas	0.078	09:04
03 Room Air	0.000	09:05
04 Std. Gas	0.078	09:05
05 Room Air	0.000	09:06
06 Std. Gas	0.079	09:06
07 Room Air	0.000	09:06

Lot No = 15410080A1
Cyl No = 12
Exp Date = 07/01/2012
County = 08 Oper No. = 555555



Operator Signature
ASHLEY AYD

Remarks: Calibration check 0.080 AC

Form 106-I8000

80-003069
P. 103 of 104
KM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003069
Location = TOXL 8164.13.00 06/09
04/06/2011 09:07

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:07
02 Std. Gas	0.079	09:08
03 Room Air	0.000	09:08
04 Std. Gas	0.079	09:08
05 Room Air	0.000	09:09
06 Std. Gas	0.079	09:09
07 Room Air	0.000	09:10

Lot No = 15410080A1
Cyl No = 12
Exp Date = 07/01/2012
County = 08 Oper No. = 555555



Operator Signature
ASHLEY AYD

Remarks: *Calibration check 0.080 AC*

Form 106-I8000

*80-003069
P. 104 of 104
FM*