

INTOXILYZER 8000 INSPECTION

Serial No.: 80-004937 Location: TOXL

- |   | Check When Done |
|---|-----------------|
| A. Is the warm-up time less than 20 minutes?  | <u>✓</u>        |
| B. Is a three-pronged grounded outlet used?   | <u>✓</u>        |
| C. Is the breath tube heated?   | <u>✓</u>        |
| D. Is the diagnostic check complete?  | <u>✓</u>        |
| E. Is the time, date, and year correct (re-set if necessary)?   | <u>✓</u>        |
| F. Print test. (Attach test record.)  | <u>✓</u>        |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR3455</u> Lot No. <u>10290</u> Rep AC <u>0.020</u>   | <u>✓</u>        |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR3847</u> Lot No. <u>11110</u> Rep AC <u>0.300</u>  | <u>✓</u>        |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.<br>(Attach test record.)<br>Sim SN <u>DR3377</u> Lot No. <u>    </u> Rep AC <u>0.100 AC + ACETONE</u> | <u>✓</u>        |
| J. RFI check. Run CMS mode. Key radio on first room air.<br>(Attach test record.)   | <u>✓</u>        |
| K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.<br>(Attach test record.)   |                 |

Lot No. <u>23411080A1</u>	Cyl No. <u>65</u>	Exp. Date <u>10/01/2013</u>
Test 1 <u>0.081</u>	Test 4 <u>0.080</u>	Test 7 <u>0.081</u>
Test 2 <u>0.081</u>	Test 5 <u>0.080</u>	Test 8 <u>0.081</u>
Test 3 <u>0.081</u>	Test 6 <u>0.080</u>	Test 9 <u>0.081</u>
Average <u>0.081</u>		

L. Tank pressure, Level 3 Func D Sub G P: 532 psi Reg 525 psi

M. Remarks/Maintenance record: ANNUAL INSPECTION  
19 APR 12 Date Charles E. Du Field Inspector's Signature  
19 Apr 2012 Date Deb Kashur Reviewed by

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004937  
Location = TOXL                    8164.13.00 06/09  
04/19/2012                            14:23

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?         1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            23411080A1  
Standard Cyl #?            65  
Standard Expiration?       10/01/2013  
Oper No?                    666666

Flow Cal. Date:            01/18/2010  
Slope                       730  
Intercept                  -779991

IR Calibration Date:        05/24/2011  
                                    3um                    9um

-----  
0th Coef(\*100):            -12537                    -15884  
1st Coef(\*100):            262613                    132554  
2nd Coef(\*100):            1600                      951  
H2O adj(mg/l\*10k):        529                        537

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:    19 APR 12    PRINT TEST

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004937  
Location = TOXL      8164.13.00 06/09  
04/19/2012      14:34

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:35
02 Std. Sol.	0.020	14:36
03 Room Air	0.000	14:36
04 Std. Sol.	0.022	14:37
05 Room Air	0.000	14:38
06 Std. Sol.	0.021	14:38
07 Room Air	0.000	14:39

08 Sim Temp = 34.0°C

Simul Ser No = 000000

Std Sol No = 0000

County = 08

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

Low AC  
0.020AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004937  
Location = TOXL      8164.13.00 06/09  
04/19/2012      14:51

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:51
02 Std. Sol.	0.295	14:52
03 Room Air	0.000	14:53
04 Std. Sol.	0.295	14:53
05 Room Air	0.000	14:54
06 Std. Sol.	0.297	14:55
07 Room Air	0.000	14:55

08 Sim Temp = 34.0°C

Simul Ser No = 000000

Std Sol No = 0000

County = 08

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

LINEARITY TEST  
0.300 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004937  
Location = TOXL                8164.13.00 06/09  
04/19/2012                        14:57

Test	AC	Time
01 Room Air	0.000	14:58
02 *Subject Test	INT*	14:59
03 Room Air	0.000	14:59

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = QWERTYUIOP  
Dr. Lic. = ND/DIS821456  
Simul Ser No = 000000  
Std Sol No = 0000  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *INTERFERENT CHECK*  
*0.100 AC + ACETONE*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004937  
Location = TOXL                    8164.13.00 06/09  
04/19/2012                            15:01

Test	AC	Time
01 Diagnostic	OK	15:02
02 Room Air	RFI*	15:02
03 Room Air	0.000	15:02

\*Invalid Test  
  Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = ASDFGHJKL  
Dr. Lic. = ND/DIS821456  
Simul Ser No = 000000  
Std Sol No = 0000  
County = 08                            Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks: *RFI CHECK*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004937  
Location = TOXL                    8164.13.00 06/09  
04/19/2012                            15:04

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:05
02 Std. Gas	0.081	15:05
03 Room Air	0.000	15:05
04 Std. Gas	0.081	15:06
05 Room Air	0.000	15:06
06 Std. Gas	0.081	15:07
07 Room Air	0.000	15:07

Lot No = 23411080A1  
Cyl No = 65  
Exp Date = 10/01/2013  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK

~~0.0~~ 0.080 AC  
CEC

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004937  
Location = TOXL                    8164.13.00 06/09  
04/19/2012                            15:09

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:10
02 Std. Gas	0.080	15:10
03 Room Air	0.000	15:11
04 Std. Gas	0.080	15:11
05 Room Air	0.000	15:12
06 Std. Gas	0.080	15:12
07 Room Air	0.000	15:13

Lot No = 23411080A1  
Cyl No = 65  
Exp Date = 10/01/2013  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK

0.080 AC

Form 106-I8000




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004937  
Location = TOXL      8164.13.00 06/09  
04/19/2012      15:14

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:14
02 Std. Gas	0.081	15:15
03 Room Air	0.000	15:15
04 Std. Gas	0.081	15:16
05 Room Air	0.000	15:16
06 Std. Gas	0.081	15:16
07 Room Air	0.000	15:17

Lot No = 23411080A1  
Cyl No = 65  
Exp Date = 10/01/2013  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC