

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-006496

Location: TOXL

- Check When Done
- A. Is the breath tube heated? ✓
  - B. Is the diagnostic check complete? ✓
  - C. Is the time, date, and year correct (re-set if necessary)? ✓
  - D. Print test. (Attach test record.) ✓
  - E. Low AC. Use  $\leq 0.03$  AC in ACA mode. (Attach test record.)  
Sim SN DR5114 Lot No. ~~201505B-6~~ Rep AC 0.030  
201505B
  - F. High AC. Use  $\geq 0.25$  AC in ACA mode. (Attach test record.)  
Sim SN MP3068 Lot No. 14102 Rep AC 0.400
  - G. Interferent check. Use 0.05% acetone in ethanol solution in ABA mode.  
(Attach test record.)  
Sim SN DR5135 Lot No. ICS4 Solution 0.10AC + 0.05% ACETONE
  - H. RFI check. Run CMS mode. Key radio during CMS test.  
(Attach test record.) ✓
  - I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.  
(Attach test record.) 0.080 AC

Lot No. 13314080A3 Cyl No. 028 Exp. Date 6/5/2016

Test 1 <u>0.080</u>	Test 4 <u>0.081</u>	Test 7 <u>0.081</u>
Test 2 <u>0.080</u>	Test 5 <u>0.081</u>	Test 8 <u>0.082</u>
Test 3 <u>0.081</u>	Test 6 <u>0.081</u>	Test 9 <u>0.081</u>

Average 0.081

J. Tank pressure, Level 3 Func D Sub G P: 480 psi Reg ~500 psi

K. Remarks/Maintenance record: ANNUAL INSPECTION

4.12.2016  
Date  
12 Apr 2016  
Date

Charles E. Edwards  
Field Inspector's Signature  
Deb Shanover  
Reviewed by

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer           Alcohol Analyzer  
North Dakota Model 8000         SN 80-006496  
Location = TOXL                 8164.13.00 06/09  
04/12/2016                         15:03

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=  
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=  
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

Current Instrument Setup

Data Entry Mode:           Enabled  
Start Test Sequence:       DABACABA  
Display Prelim Rslt?       Yes  
Display Third Digit?       Yes  
Inhib Printer(Y/N)?       No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?         1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            13314080A3  
Standard Cyl #?            28  
Standard Expiration?       06/05/2016  
Oper No?                    666666

Flow Cal. Date:            08/18/2015  
Slope                       673  
Intercept                  -654167

IR Calibration Date:       08/18/2015  
                              3um                   9um

	3um	9um
0th Coef(*100):	-22081	-15963
1st Coef(*100):	278490	136818
2nd Coef(*100):	2086	1282
H2O adj(mg/l*10k):	598	386

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:

PRINT TEST

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006496  
Location = TOXL      8164.13.00 06/09  
04/12/2016      15:06

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:07
02 Std. Sol.	0.027	15:08
03 Room Air	0.000	15:08
04 Std. Sol.	0.028	15:09
05 Room Air	0.000	15:09
06 Std. Sol.	0.028	15:10
07 Room Air	0.000	15:11

08 Sim Temp = 34.0°C

Simul Ser No = DR5114  
Std Sol No = ~~14102~~ CEE 201505B  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:      Low AC      0.030 AC

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006496  
Location = TOXL      8164.13.00 06/09  
04/12/2016      15:13

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:14
02 Std. Sol.	0.402	15:15
03 Room Air	0.000	15:16
04 Std. Sol.	0.402	15:16
05 Room Air	0.000	15:17
06 Std. Sol.	0.402	15:18
07 Room Air	0.000	15:18

08 Sim Temp = 34.0°C

Simul Ser No = MP3068  
Std Sol No = 14102  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

HIGH AC  
0.400 AC


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006496  
Location = TOXL      8164.13.00 06/09  
04/12/2016      15:20

Test	AC	Time
01 Room Air	0.000	15:21
02 *Subject Test	INT*	15:22
03 Room Air	0.000	15:23

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = INTERFERENT  
Dr. Lic. = ND/DIS821456  
Lot No = 13314080A3  
Cyl No = 28  
Expiration Date = 06/05/2016  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*INTERFERENT CHECK*  
*0.10AC + 0.05% ACETONE*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-006496  
Location = TOXL                    8164.13.00 06/09  
04/12/2016                            15:24

Test	AC	Time
01 Diagnostic	OK	15:25
02 Room Air	0.000	15:25
03 *Subject Test	RFI*	15:25
04 Room Air	0.000	15:26

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = RFI CHECK  
Dr. Lic. = ND/DIS821456  
Lot No = 13314080A3  
Cyl No = 28  
Expiration Date = 06/05/2016  
County = 08                            Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

RFI CHECK

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006496  
Location = TOXL      8164.13.00 06/09  
04/12/2016      15:28

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:29
02 Std. Gas	0.080	15:29
03 Room Air	0.000	15:30
04 Std. Gas	0.080	15:30
05 Room Air	0.000	15:31
06 Std. Gas	0.081	15:31
07 Room Air	0.000	15:32

Lot No = 13314080A3  
Cyl No = 28  
Exp Date = 06/05/2016  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006496  
Location = TOXL                    8164.13.00 06/09  
04/12/2016                            15:32

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:32
02 Std. Gas	0.081	15:33
03 Room Air	0.000	15:33
04 Std. Gas	0.081	15:34
05 Room Air	0.000	15:34
06 Std. Gas	0.081	15:35
07 Room Air	0.000	15:35

Lot No = 13314080A3  
Cyl No = 28  
Exp Date = 06/05/2016  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080AC




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006496  
Location = TOXL                    8164.13.00 06/09  
04/12/2016                            15:36

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:36
02 Std. Gas	0.081	15:37
03 Room Air	0.000	15:37
04 Std. Gas	0.082	15:38
05 Room Air	0.000	15:38
06 Std. Gas	0.081	15:39
07 Room Air	0.000	15:39

Lot No = 13314080A3  
Cyl No = 28  
Exp Date = 06/05/2016  
County = 08                            Oper No. = 666666

  
Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
          0.080 AC