

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-002669

Location: TOXL

- |  | Check When Done |
|--|-----------------|
| A. Is the breath tube heated?  | <u>✓</u>        |
| B. Is the diagnostic check complete?   | <u>✓</u>        |
| C. Is the time, date, and year correct (re-set if necessary)?  | <u>✓</u>        |
| D. Print test. (Attach test record.)   | <u>✓</u>        |
| E. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR3803</u> Lot No. <u>13070</u> Rep AC <u>0.020</u>  | <u>✓</u>        |
| F. High AC. Use ≥ 0.25 AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR3374</u> Lot No. <u>13080</u> Rep AC <u>0.300</u>  | <u>✓</u>        |
| G. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.<br>(Attach test record.)<br>Sim SN <u>DR3375</u> Lot No. <u>ICS 2</u> Rep AC <u>0.10AC + 0.05% acetone</u> | <u>✓</u>        |
| H. RFI check. Run CMS mode. Key radio on first room air.<br>(Attach test record.)  | <u>✓</u>        |
| I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.<br>(Attach test record.)  |                 |

Lot No. <u>165/2080A2</u>	Cyl No. <u>22</u>	Exp. Date <u>08/01/2014</u>
Test 1 <u>0.080</u>	Test 4 <u>0.079</u>	Test 7 <u>0.079</u>
Test 2 <u>0.080</u>	Test 5 <u>0.079</u>	Test 8 <u>0.080</u>
Test 3 <u>0.079</u>	Test 6 <u>0.080</u>	Test 9 <u>0.080</u>
Average <u>0.080</u>		

J. Tank pressure, Level 3 Func D Sub G P: 439 psi Reg 450 psi

K. Remarks/Maintenance record: Annual inspection, replaced simulator return o-ring.

11 Apr 2014  
Date

14 APR 2014 ✓  
Date

Deb Kashner  
Field Inspector's Signature

[Signature]  
Reviewed by

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-002669  
Location = TOXL                8164.13.00 06/09  
04/11/2014                      09:29

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?         1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            16512080A2  
Standard Cyl #?            22  
Standard Expiration?       08/01/2014  
Oper No?                    777777

Flow Cal. Date:            09/09/2011  
Slope                      662  
Intercept                  -510765

IR Calibration Date:       05/21/2013

	3um	9um
0th Coef(*100):	-24386	-12187
1st Coef(*100):	269217	134678
2nd Coef(*100):	3135	1513
H2O adj(mg/l*10k):	637	486

\*\*\*\*\* Printer Test End \*\*\*\*\*

Deb Kashur

Operator Signature  
DEB KASHUR

Remarks:

Print test

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-002669  
Location = ~~MNAF~~ *TOXL*      8164.13.00 06/09  
04/11/2014      *OK*      09:17

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:18
02 Std. Sol.	0.019	09:19
03 Room Air	0.000	09:19
04 Std. Sol.	0.019	09:20
05 Room Air	0.000	09:20
06 Std. Sol.	0.020	09:21
07 Room Air	0.000	09:22

08 Sim Temp = 34.0°C

Simul Ser No = DR3803

Std Sol No = 13070

County = 08

Oper No. = 777777

*Deb Kasher*

Operator Signature  
DEB KASHUR

Remarks:

*Low AC      0.020 Ae*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-002669  
Location = ~~MNAF TOXL~~      8164.13.00 06/09  
04/11/2014      DK      09:05

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:05
02 Std. Sol.	0.305	09:06
03 Room Air	0.000	09:07
04 Std. Sol.	0.304	09:07
05 Room Air	0.000	09:08
06 Std. Sol.	0.304	09:09
07 Room Air	0.000	09:09

08 Sim Temp = 34.0°C

Simul Ser No = DR3374

Std Sol No = 13080

County = 08

Oper No. = 777777

*Deb Kashur*

Operator Signature  
DEB KASHUR

Remarks:

*High AC      0.300 AC*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-002669  
Location = TOXL      8164.13.00 06/09  
04/11/2014      09:30

Test	AC	Time
01 Room Air	0.000	09:30
02 *Subject Test	INT*	09:31
03 Room Air	0.000	09:31

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = NA  
Dr. Lic. = ND/DIS921456  
Lot No = 16512080A2  
Cyl No = 22  
Expiration Date = 08/01/2014  
County = 08      Oper No. = 777777

Deb Kashur

Operator Signature  
DEB KASHUR

Remarks:

*Interferent check 0.10AC + 0.05x acetone*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-002669  
Location = TOXL      8164.13.00 06/09  
04/11/2014      09:32

Test	AC	Time
01 Diagnostic	OK	09:32
02 Room Air	RFI*	09:32
03 Room Air	0.000	09:33

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = NA  
Dr. Lic. = ND/DIS921456  
Lot No = 16512080A2  
Cyl No = 22  
Expiration Date = 08/01/2014  
County = 08      Oper No. = 777777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Deb Kasher

Operator Signature  
DEB KASHUR

Remarks:

RFI check

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-002669  
Location = ~~MNAF Tox~~      8164.13.00 06/09  
04/11/2014      ~~DK~~      09:00

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:01
02 Std. Gas	0.080	09:01
03 Room Air	0.000	09:02
04 Std. Gas	0.080	09:02
05 Room Air	0.000	09:03
06 Std. Gas	0.079	09:03
07 Room Air	0.000	09:04

Lot No = 16512080A2  
Cyl No = 22  
Exp Date = 08/01/2014  
County = 08      Oper No. = 777777

Deb Kasher  
Operator Signature  
DEB KASHUR

Remarks:

Calibration check 0080AC

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-002669  
Location = TOXL      8164.13.00 06/09  
04/11/2014      09:34

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:34
02 Std. Gas	0.079	09:34
03 Room Air	0.000	09:35
04 Std. Gas	0.079	09:35
05 Room Air	0.000	09:36
06 Std. Gas	0.080	09:36
07 Room Air	0.000	09:37

Lot No = 16512080A2  
Cyl No = 22  
Exp Date = 08/01/2014  
County = 08      Oper No. = 777777

*Deb Kasher*

Operator Signature  
DEB KASHUR

Remarks:

*Calibration check 0.080 AC*

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-002669  
Location = TOXL      8164.13.00 06/09  
04/11/2014      09:37

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:38
02 Std. Gas	0.079	09:38
03 Room Air	0.000	09:39
04 Std. Gas	0.080	09:39
05 Room Air	0.000	09:39
06 Std. Gas	0.080	09:40
07 Room Air	0.000	09:40

Lot No = 16512080A2

Cyl No = 22

Exp Date = 08/01/2014

County = 08

Oper No. = 777777

Deb Kasher

Operator Signature  
DEB KASHUR

Remarks:

Calibration check      0.080 AC

Form 106-I8000