

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-004199

Location: TOXL

- |  | Check When Done |
|--|-----------------|
| A. Is the breath tube heated?  | <u>✓</u>        |
| B. Is the diagnostic check complete?   | <u>✓</u>        |
| C. Is the time, date, and year correct (re-set if necessary)?  | <u>✓</u>        |
| D. Print test. (Attach test record.)   | <u>✓</u>        |
| E. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR3803</u> Lot No. <u>13070</u> Rep AC <u>0.020</u>  | <u>✓</u>        |
| F. High AC. Use ≥ 0.25 AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR3374</u> Lot No. <u>13080</u> Rep AC <u>0.300</u>  | <u>✓</u>        |
| G. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.<br>(Attach test record.)<br>Sim SN <u>DR 3375</u> Lot No. <u>ICS 2</u> Rep <u>AC 0.10 AC + 0.05% acetone</u> | <u>✓</u>        |
| H. RFI check. Run CMS mode. Key radio on first room air.<br>(Attach test record.)  | <u>✓</u>        |
| I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.<br>(Attach test record.)  |                 |

Lot No. <u>16512080AA</u>	Cyl No. <u>2a</u>	Exp. Date <u>08/01/2014</u>
Test 1 <u>0.080</u>	Test 4 <u>0.081</u>	Test 7 <u>0.080</u>
Test 2 <u>0.081</u>	Test 5 <u>0.081</u>	Test 8 <u>0.081</u>
Test 3 <u>0.080</u>	Test 6 <u>0.080</u>	Test 9 <u>0.080</u>
Average <u>0.080</u>		

J. Tank pressure, Level 3 Func D Sub G P: 445 psi Reg 450 psi

K. Remarks/Maintenance record: Annual inspection

10 Apr 2014  
Date

14 APR 14 ✓  
Date

Deb Kashner  
Field Inspector's Signature

[Signature]  
Reviewed by



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004199  
Location = TOXL      8164.13.00 06/09  
04/10/2014      16:45

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:46
02 Std. Sol.	0.019	16:46
03 Room Air	0.000	16:47
04 Std. Sol.	0.020	16:48
05 Room Air	0.000	16:48
06 Std. Sol.	0.021	16:49
07 Room Air	0.000	16:49

08 Sim Temp = 34.0°C

Simul Ser No = DR3803

Std Sol No = 13070

County = 08

Oper No. = 777777

*Deb Kasher*

Operator Signature  
DEB KASHUR

Remarks:

*Low AC      0.020 AC*

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004199  
Location = TOXL      8164.13.00 06/09  
04/10/2014      16:30

Test	AC	Time
01 Room Air	0.000	16:31
02 *Subject Test	INT*	16:31
03 Room Air	0.000	16:32

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = NA  
Dr. Lic. = ND/DIS921456  
Lot No = 16512080A2  
Cyl No = 22  
Expiration Date = 08/01/2014  
County = 08      Oper No. = 777777

Deb Kashur

Operator Signature  
DEB KASHUR

Remarks:

*Interferent check 0.100AE + 0.05% acetone*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004199  
Location = TOXL                    8164.13.00 06/09  
04/10/2014                            16:55

Test	AC	Time
01 Diagnostic	OK	16:56
02 Room Air	RFI*	16:56
03 Room Air	0.000	16:56

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = NA  
Dr. Lic. = ND/DIS921456  
Lot No = 16512080A2  
Cyl No = 22  
Expiration Date = 08/01/2014  
County = 08                            Oper No. = 777777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
DEB KASHUR

Remarks:

RFI check

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004199  
Location = TOXL      8164.13.00 06/09  
04/10/2014      16:16

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:17
02 Std. Gas	0.080	16:17
03 Room Air	0.000	16:17
04 Std. Gas	0.081	16:18
05 Room Air	0.000	16:18
06 Std. Gas	0.080	16:19
07 Room Air	0.000	16:19

Lot No = 16512080A2  
Cyl No = 22  
Exp Date = 08/01/2014  
County = 08      Oper No. = 777777

Deb Kasher  
Operator Signature  
DEB KASHUR

Remarks:

Calibration check      0.080AC

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004199  
Location = TOXL      8164.13.00 06/09  
04/10/2014      16:20

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:20
02 Std. Gas	0.081	16:21
03 Room Air	0.000	16:21
04 Std. Gas	0.081	16:21
05 Room Air	0.000	16:22
06 Std. Gas	0.080	16:22
07 Room Air	0.000	16:23

Lot No = 16512080A2  
Cyl No = 22  
Exp Date = 08/01/2014  
County = 08      Oper No. = 777777

*Deb Kasher*

Operator Signature  
DEB KASHUR

Remarks:

*Calibration check 0.080 AC*

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004199  
Location = TOXL      8164.13.00 06/09  
04/10/2014      16:25

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:25
02 Std. Gas	0.080	16:26
03 Room Air	0.000	16:26
04 Std. Gas	0.081	16:27
05 Room Air	0.000	16:27
06 Std. Gas	0.080	16:28
07 Room Air	0.000	16:28

Lot No = 16512080A2  
Cyl No = 22  
Exp Date = 08/01/2014  
County = 08      Oper No. = 777777

Deb Kasher  
Operator Signature  
DEB KASHUR

Remarks:  
Calibration check 0.080AC