

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-004953

Location: TOXL

- |  | Check When Done |
|--|-----------------|
| A. Is the breath tube heated?  | <u>✓</u>        |
| B. Is the diagnostic check complete?   | <u>✓</u>        |
| C. Is the time, date, and year correct (re-set if necessary)?  | <u>✓</u>        |
| D. Print test. (Attach test record.)   | <u>✓</u>        |
| E. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR 3803</u> Lot No. <u>13070</u> Rep AC <u>0.020</u>   | <u>✓</u>        |
| F. High AC. Use ≥ 0.25 AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR 3374</u> Lot No. <u>13080</u> Rep AC <u>0.300</u>   | <u>✓</u>        |
| G. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.<br>(Attach test record.)<br>Sim SN <u>DR 3375</u> Lot No. <u>ICS2</u> Rep AC <u>0.10Ac + 0.05% acetone</u> | <u>✓</u>        |
| H. RFI check. Run CMS mode. Key radio on first room air.<br>(Attach test record.)  | <u>✓</u>        |
| I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.<br>(Attach test record.)  |                 |

Lot No. <u>16512080A2</u>	Cyl No. <u>22</u>	Exp. Date <u>08/01/2014</u>
Test 1 <u>0.079</u>	Test 4 <u>0.080</u>	Test 7 <u>0.081</u>
Test 2 <u>0.080</u>	Test 5 <u>0.080</u>	Test 8 <u>0.081</u>
Test 3 <u>0.080</u>	Test 6 <u>0.080</u>	Test 9 <u>0.081</u>
Average <u>0.080</u>		

J. Tank pressure, Level 3 Func D Sub G P: 407 psi Reg 400 psi

K. Remarks/Maintenance record: Annual inspection, replaced o-ring on simulator return

14 Apr 2014  
Date

14 APR 2014 ✓  
Date

Deb Kashur  
Field Inspector's Signature

Charles E. Egan  
Reviewed by

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004953  
Location = TOXL      8164.13.00 06/09  
04/14/2014      13:49

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%^&\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%^&\*()\_+?

Current Instrument Setup

Data Entry Mode:      Enabled  
Start Test Sequence:      DABACABA  
Display Prelim Rslt?      Yes  
Display Third Digit?      Yes  
Inhib Printer(Y/N)?      No  
Display Volume?      No  
Disable On Memfull?      Yes  
# of Print Copies?      1  
Select Std (D/W/I)?      Dry  
Standard Value?      0.080  
Standard Lot #?      16512080A2  
Standard Cyl #?      22  
Standard Expiration?      08/01/2014  
Oper No?      777777

Flow Cal. Date:      01/19/2010  
Slope      691  
Intercept      -628068

IR Calibration Date:      06/08/2010

	3um	9um
0th Coef(*100):	-10343	-18245
1st Coef(*100):	249276	131289
2nd Coef(*100):	4014	1439
H2O adj(mg/l*10k):	697	655

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
DEB KASHUR

Remarks:

*Print Test*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004953  
Location = TOXL      8164.13.00 06/09  
04/14/2014      10:55

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:55
02 Std. Sol.	0.020	10:56
03 Room Air	0.000	10:57
04 Std. Sol.	0.022	10:57
05 Room Air	0.000	10:58
06 Std. Sol.	0.021	10:59
07 Room Air	0.000	10:59

08 Sim Temp = 34.0°C

Simul Ser No = DR3803

Std Sol No = 13070

County = 08

Oper No. = 777777



Operator Signature  
DEB KASHUR

Remarks:

Low Ac 0.020 AC

Form 106-18000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004953  
Location = TOXL      8164.13.00 06/09  
04/14/2014      10:45

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:46
02 Std. Sol.	0.301	10:47
03 Room Air	0.000	10:47
04 Std. Sol.	0.303	10:48
05 Room Air	0.000	10:48
06 Std. Sol.	0.303	10:49
07 Room Air	0.000	10:50

08 Sim Temp = 34.0°C

Simul Ser No = DR3374

Std Sol No = 13080

County = 08

Oper No. = 777777

*Deb Kashur*

Operator Signature  
DEB KASHUR

Remarks:

*High AC 0.300AC*

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004953  
Location = TOXL      8164.13.00 06/09  
04/14/2014      13:46

Test	AC	Time
01 Room Air	0.000	13:47
02 *Subject Test	INT*	13:47
03 Room Air	0.000	13:48

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = NA  
Dr. Lic. = ND/DIS921456  
Lot No = 16512080A2  
Cyl No = 22  
Expiration Date = 08/01/2014  
County = 08      Oper No. = 777777

  
\_\_\_\_\_  
Operator Signature  
DEB KASHUR

Remarks:

*Interferent check 0.05% acetone + 0.10 AC*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004953  
Location = TOXL                8164.13.00 06/09  
04/14/2014                        13:50

Test	AC	Time
01 Diagnostic	OK	13:51
02 Room Air	RFI*	13:51
03 Room Air	0.000	13:51

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                Weight = 150  
Test = DUI                        Cit = NA  
Dr. Lic. = ND/DIS921456  
Lot No = 16512080A2  
Cyl No = 22  
Expiration Date = 08/01/2014  
County = 08                        Oper No. = 777777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
DEB KASHUR

Remarks:

*RFI check*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004953  
Location = TOXL      8164.13.00 06/09  
04/14/2014      10:41

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:41
02 Std. Gas	0.079	10:41
03 Room Air	0.000	10:42
04 Std. Gas	0.080	10:42
05 Room Air	0.000	10:43
06 Std. Gas	0.080	10:43
07 Room Air	0.000	10:44

Lot No = 16512080A2  
Cyl No = 22  
Exp Date = 08/01/2014  
County = 08      Oper No. = 777777

*Deb Kashur*

Operator Signature  
DEB KASHUR

Remarks:

*Calibration check 0.080 AC*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004953  
Location = TOXL      8164.13.00 06/09  
04/14/2014      13:52

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:52
02 Std. Gas	0.080	13:53
03 Room Air	0.000	13:53
04 Std. Gas	0.080	13:53
05 Room Air	0.000	13:54
06 Std. Gas	0.080	13:54
07 Room Air	0.000	13:55

Lot No = 16512080A2  
Cyl No = 22  
Exp Date = 08/01/2014  
County = 08      Oper No. = 777777

Deb Kashur

Operator Signature  
DEB KASHUR

Remarks:

Calibration check 0.080AC

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004953  
Location = TOXL      8164.13.00 06/09  
04/14/2014      13:57

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:57
02 Std. Gas	0.081	13:58
03 Room Air	0.000	13:58
04 Std. Gas	0.081	13:59
05 Room Air	0.000	13:59
06 Std. Gas	0.081	13:59
07 Room Air	0.000	14:00

Lot No = 16512080A2

Cyl No = 22

Exp Date = 08/01/2014

County = 08

Oper No. = 777777

*Deb Kashur*

Operator Signature  
DEB KASHUR

Remarks:

*Calibration check 0.080AC*

Form 106-I8000