

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-003069

Location: TOXL

- |  | Check When Done |
|--|-----------------|
| A. Is the breath tube heated?  | <u>✓</u>        |
| B. Is the diagnostic check complete?   | <u>✓</u>        |
| C. Is the time, date, and year correct (re-set if necessary)?  | <u>✓</u>        |
| D. Print test. (Attach test record.)   | <u>✓</u>        |
| E. Low AC. Use $\leq 0.03$ AC in ACA mode. (Attach test record.)<br>Sim SN <u>MP3066</u> Lot No. <u>15160</u> Rep AC <u>0.020</u>  | <u>✓</u>        |
| F. High AC. Use $\geq 0.25$ AC in ACA mode. (Attach test record.)<br>Sim SN <u>MP3068</u> Lot No. <u>14102</u> Rep AC <u>0.400</u>   | <u>✓</u>        |
| G. Interferent check. Use 0.05% acetone in ethanol solution in ABA mode.<br>(Attach test record.)<br>Sim SN <u>DR5135</u> Lot No. <u>ICS4</u> Solution <u>0.10AC + 0.05% ACETONE</u> | <u>✓</u>        |
| H. RFI check. Run CMS mode. Key radio during CMS test.<br>(Attach test record.)  | <u>✓</u>        |
| I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.<br>(Attach test record.) <u>0.080</u> AC  |                 |

Lot No. 13314080A3 Cyl No. 028 Exp. Date 6/5/2016

Test 1 0.079 Test 4 0.079 Test 7 0.079

Test 2 0.077 Test 5 0.078 Test 8 0.078

Test 3 0.079 Test 6 0.078 Test 9 0.078

Average 0.078

J. Tank pressure, Level 3 Func D Sub G P: 290 psi Reg ~275 psi

K. Remarks/Maintenance record: ANNUAL INSPECTION, REPLACED BREATH HOSE O-RING

4.14.2016  
Date

[Signature]  
Field Inspector's Signature

14 Apr 2016  
Date

[Signature]  
Reviewed by

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer          Alcohol Analyzer  
North Dakota Model 8000          SN 80-003069  
Location = TOXL                  8164.13.00 06/09  
04/14/2016                                  09:21

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890- = |  
ABCDEFGHIJKLMNOpqrstuvwxyz!@#%&^\* () \_ + ?

abcdefghijklmnopqrstuvwxyz1234567890- = |  
ABCDEFGHIJKLMNOpqrstuvwxyz!@#%&^\* () \_ + ?

Current Instrument Setup


Data Entry Mode:                  Enabled  
Start Test Sequence:              DABACABA  
Display Prelim Rslt?              Yes  
Display Third Digit?              Yes  
Inhib Printer(Y/N)?              No  
Display Volume?                  No  
Disable On Memfull?              Yes  
# of Print Copies?                1  
Select Std (D/W/I)?              Dry  
Standard Value?                  0.080  
Standard Lot #?                  14214080A6  
Standard Cyl #?                  72  
Standard Expiration?              07/05/2016  
Oper No?                          666666

Flow Cal. Date:                  09/17/2013  
Slope                              678  
Intercept                         -734495

IR Calibration Date:              04/08/2015  
    3um                          9um

-----  
0th Coef(\*100):                  -13291                          -10923  
1st Coef(\*100):                  270590                         132326  
2nd Coef(\*100):                  4128                             1532  
H2O adj(mg/l\*10k):              348                              226

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:

PRINT TEST

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003069  
Location = TOXL      8164.13.00 06/09  
04/14/2016      09:24

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:24
02 Std. Sol.	0.019	09:25
03 Room Air	0.000	09:26
04 Std. Sol.	0.019	09:26
05 Room Air	0.000	09:27
06 Std. Sol.	0.019	09:28
07 Room Air	0.000	09:28

08 Sim Temp = 34.0°C

Simul Ser No = MP3066  
Std Sol No = 15160  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *Low AC 0.020 AC*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003069  
Location = TOXL      8164.13.00 06/09  
04/14/2016      09:30

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:31
02 Std. Sol.	0.397	09:32
03 Room Air	0.000	09:32
04 Std. Sol.	0.398	09:33
05 Room Air	0.000	09:34
06 Std. Sol.	0.399	09:34
07 Room Air	0.000	09:35

08 Sim Temp = 34.0°C

Simul Ser No = MP3068  
Std Sol No = 14102  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *HIGH AC 0.400*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003069  
Location = TOXL      8164.13.00 06/09  
04/14/2016      09:37

Test	AC	Time
01 Room Air	0.000	09:38
02 *Subject Test	INT*	09:39
03 Room Air	0.000	09:39

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = INTERFERENT  
Dr. Lic. = ND/DIS821456  
Lot No = 13314080A3  
Cyl No = 28  
Expiration Date = 06/05/2016  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *Interferent Check*  
*0.10 AC + 0.05% ACETONE*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-003069  
Location = TOXL                    8164.13.00 06/09  
04/14/2016                            09:44

Test	AC	Time
01 Diagnostic	OK	09:45
02 Room Air	0.000	09:45
03 *Subject Test	RFI*	09:46
04 Room Air	0.000	09:46

\*Invalid Test  
  Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = RFI CHECK  
Dr. Lic. = ND/DIS821456  
Lot No = 13314080A3  
Cyl No = 28  
Expiration Date = 06/05/2016  
County = 08                            Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

*RFI CHECK*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003069  
Location = TOXL      8164.13.00 06/09  
04/14/2016      09:40

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:41
02 Std. Gas	0.079	09:41
03 Room Air	0.000	09:42
04 Std. Gas	0.077	09:42
05 Room Air	0.000	09:43
06 Std. Gas	0.079	09:43
07 Room Air	0.000	09:43

Lot No = 13314080A3  
Cyl No = 28  
Exp Date = 06/05/2016  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003069  
Location = TOXL      8164.13.00 06/09  
04/14/2016      09:47

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:47
02 Std. Gas	0.079	09:48
03 Room Air	0.000	09:48
04 Std. Gas	0.078	09:49
05 Room Air	0.000	09:49
06 Std. Gas	0.078	09:49
07 Room Air	0.000	09:50

Lot No = 13314080A3  
Cyl No = 28  
Exp Date = 06/05/2016  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK

0.080 AC

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003069  
Location = TOXL      8164.13.00 06/09  
04/14/2016      09:50

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:51
02 Std. Gas	0.079	09:51
03 Room Air	0.000	09:52
04 Std. Gas	0.078	09:52
05 Room Air	0.000	09:52
06 Std. Gas	0.078	09:53
07 Room Air	0.000	09:53

Lot No = 13314080A3  
Cyl No = 28  
Exp Date = 06/05/2016  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000