

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-004949

Location: TOXL

- | | Check When Done |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|
| A. Is the breath tube heated? | <u>✓</u> |
| B. Is the diagnostic check complete? | <u>✓</u> |
| C. Is the time, date, and year correct (re-set if necessary)? | <u>✓</u> |
| D. Print test. (Attach test record.) | <u>✓</u> |
| E. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)
Sim SN <u>DR3377</u> Lot No. <u>13270</u> Rep AC <u>0.020</u> | <u>✓</u> |
| F. High AC. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN <u>DR5189</u> Lot No. <u>13300</u> Rep AC <u>0.300</u> | <u>✓</u> |
| G. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.
(Attach test record.)
Sim SN <u>DR5135</u> Lot No. <u>ICS#2</u> Rep AC <u>0.100AC + 0.05% ACETONE</u> | <u>✓</u> |
| H. RFI check. Run CMS mode. Key radio on first room air.
(Attach test record.) | <u>✓</u> |
| I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.) | |

Lot No. <u>16512080A2</u>	Cyl No. <u>034</u>	Exp. Date <u>08/01/2014</u>
Test 1 <u>0.079</u>	Test 4 <u>0.080</u>	Test 7 <u>0.080</u>
Test 2 <u>0.079</u>	Test 5 <u>0.080</u>	Test 8 <u>0.081</u>
Test 3 <u>0.079</u>	Test 6 <u>0.080</u>	Test 9 <u>0.080</u>
Average <u>0.080</u>		

J. Tank pressure, Level 3 Func D Sub G P: 92 psi Reg 100 psi

K. Remarks/Maintenance record: REPLACED O-RING ON SIMULATOR RETURN,

ANNUAL INSPECTION

15 APR 2014
Date

15 Apr 2014
Date

Charles E. Eden
Field Inspector's Signature

Deb Kashur
Reviewed by

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004949
Location = TOXL 8164.13.00 06/09
04/15/2014 10:32

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-|=|
ABCDEFGHIJKLMNPNOPQRSTUVWXYZ!@#%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-|=|
ABCDEFGHIJKLMNPNOPQRSTUVWXYZ!@#%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 16512080A2
Standard Cyl #? 34
Standard Expiration? 08/01/2014
Oper No? 666666

Flow Cal. Date: 09/24/2010
Slope 784
Intercept -862571

IR Calibration Date: 09/24/2010
 3um 9um

0th Coef(*100): -23764 -22121
1st Coef(*100): 267615 132337
2nd Coef(*100): 1069 992
H2O adj(mg/l*10k): 531 480

***** Printer Test End *****



Operator Signature
CHARLES EDER

Remarks:

PRINT TEST

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004949
Location = TOXL 8164.13.00 06/09
04/15/2014 10:35

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:35
02 Std. Sol.	0.020	10:36
03 Room Air	0.000	10:37
04 Std. Sol.	0.019	10:37
05 Room Air	0.000	10:38
06 Std. Sol.	0.019	10:39
07 Room Air	0.000	10:39

08 Sim Temp = 34.0°C

Simul Ser No = DR3377
Std Sol No = 13270
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

LOW AC
0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004949
Location = TOXL 8164.13.00 06/09
04/15/2014 10:40

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:41
02 Std. Sol.	0.297	10:42
03 Room Air	0.000	10:43
04 Std. Sol.	0.297	10:43
05 Room Air	0.000	10:44
06 Std. Sol.	0.295	10:45
07 Room Air	0.000	10:45


08 Sim Temp = 34.0°C

Simul Ser No = DR5189

Std Sol No = 13300

County = 08

Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

HIGH AC
0.300 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004949
Location = TOXL 8164.13.00 06/09
04/15/2014 10:46

Test	AC	Time
01 Room Air	0.000	10:47
02 *Subject Test	INT*	10:48
03 Room Air	0.000	10:49

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = INTERFERENT
Dr. Lic. = ND/DIS821456
Lot No = 16512080A2
Cyl No = 34
Expiration Date = 08/01/2014
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

INTERFERENT CHECK
0.100AC + 0.05% ACETONE

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004949
Location = TOXL 8164.13.00 06/09
04/15/2014 10:50

Test	AC	Time
01 Diagnostic	OK	10:51
02 Room Air	RFI*	10:51
03 Room Air	0.000	10:51

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = RFI CHECK
Dr. Lic. = ND/DIS821456
Lot No = 16512080A2
Cyl No = 34
Expiration Date = 08/01/2014
County = 08 Oper No. = 666666

I followed the Approved Method and the
instructions displayed by the Intoxilyzer
in conducting this test.



Operator Signature
CHARLES EDER

Remarks: RFI CHECK


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004949
Location = TOXL 8164.13.00 06/09
04/15/2014 10:52

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:53
02 Std. Gas	0.079	10:53
03 Room Air	0.000	10:54
04 Std. Gas	0.079	10:54
05 Room Air	0.000	10:55
06 Std. Gas	0.079	10:55
07 Room Air	0.000	10:56

Lot No = 16512080A2
Cyl No = 34
Exp Date = 08/01/2014
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004949
Location = TOXL 8164.13.00 06/09
04/15/2014 10:56

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:56
02 Std. Gas	0.080	10:57
03 Room Air	0.000	10:57
04 Std. Gas	0.080	10:58
05 Room Air	0.000	10:58
06 Std. Gas	0.080	10:58
07 Room Air	0.000	10:59

Lot No = 16512080A2
Cyl No = 34
Exp Date = 08/01/2014
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004949
Location = TOXL 8164.13.00 06/09
04/15/2014 10:59

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:00
02 Std. Gas	0.080	11:00
03 Room Air	0.000	11:01
04 Std. Gas	0.081	11:01
05 Room Air	0.000	11:02
06 Std. Gas	0.080	11:02
07 Room Air	0.000	11:02

Lot No = 16512080A2
Cyl No = 34
Exp Date = 08/01/2014
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC