

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-004201

Location: TOXL

- Check When Done
- A. Is the breath tube heated? ✓
  - B. Is the diagnostic check complete? ✓
  - C. Is the time, date, and year correct (re-set if necessary)? ✓
  - D. Print test. (Attach test record.) ✓
  - E. Low AC. Use  $\leq 0.03$  AC in ACA mode. (Attach test record.)  
Sim SN MP3066 Lot No. 15160 Rep AC 0.020 ✓
  - F. High AC. Use  $\geq 0.25$  AC in ACA mode. (Attach test record.)  
Sim SN MP3068 Lot No. 14102 Rep AC 0.400 ✓
  - G. Interferent check. Use 0.05% acetone in ethanol solution in ABA mode.  
(Attach test record.)  
Sim SN DR5135 Lot No. ICS 4 Solution 0.10AC+0.05% ACETONE ✓
  - H. RFI check. Run CMS mode. Key radio during CMS test.  
(Attach test record.) ✓
  - I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.  
(Attach test record.) 0.080 AC

Lot No. 13314080A3 Cyl No. 028 Exp. Date 6/5/2016

Test 1 <u>0.078</u>	Test 4 <u>0.078</u>	Test 7 <u>0.078</u>
Test 2 <u>0.078</u>	Test 5 <u>0.078</u>	Test 8 <u>0.077</u>
Test 3 <u>0.078</u>	Test 6 <u>0.077</u>	Test 9 <u>0.078</u>

Average 0.078

J. Tank pressure, Level 3 Func D Sub G P: 90 psi Reg ~100 psi

K. Remarks/Maintenance record: ANNUAL INSPECTION.

4.15.2016  
Date  
15 Apr 2016  
Date

Charles E. Ed...  
Field Inspector's Signature  
Deb Shanaver  
Reviewed by



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004201  
Location = TOXL                    8164.13.00 06/09  
04/15/2016                            14:50

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:51
02 Std. Sol.	0.017	14:51
03 Room Air	0.000	14:52
04 Std. Sol.	0.017	14:53
05 Room Air	0.000	14:53
06 Std. Sol.	0.017	14:54
07 Room Air	0.000	14:54

08 Sim Temp = 34.0°C

Simul Ser No = MP3066  
Std Sol No = 15160  
County = 08                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:            *Low AC*  
                          *0.020 AC*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004201  
Location = TOXL      8164.13.00 06/09  
04/15/2016      14:56

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:56
02 Std. Sol.	0.396	14:57
03 Room Air	0.000	14:58
04 Std. Sol.	0.398	14:59
05 Room Air	0.000	14:59
06 Std. Sol.	0.399	15:00
07 Room Air	0.000	15:00

08 Sim Temp = 34.0°C

Simul Ser No = MP3068  
Std Sol No = 14102  
County = 08      Oper No. = 666666

  
Operator Signature  
CHARLES EDER

Remarks:      HIGH AC  
                    0.400 AC


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004201  
Location = TOXL      8164.13.00 06/09  
04/15/2016      15:04

Test	AC	Time
01 Room Air	0.000	15:05
02 *Subject Test	INT*	15:06
03 Room Air	0.000	15:06

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = INTERFERENT  
Dr. Lic. = ND/DIS821456  
Lot No = 13314080A3  
Cyl No = 28  
Expiration Date = 06/05/2016  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Interferent Check*  
*0.10AC + 0.05% ACETONE*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004201  
Location = TOXL                    8164.13.00 06/09  
04/15/2016                            15:07

Test	AC	Time
01 Diagnostic	OK	15:07
02 Room Air	0.000	15:08
03 *Subject Test	RFI*	15:08
04 Room Air	0.000	15:08

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = RFI CHECK  
Dr. Lic. = ND/DIS821456  
Lot No = 13314080A3  
Cyl No = 28  
Expiration Date = 06/05/2016  
County = 08                            Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks: *RFI CHECK*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004201  
Location = TOXL      8164.13.00 06/09  
04/15/2016      15:10

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:10
02 Std. Gas	0.078	15:10
03 Room Air	0.000	15:11
04 Std. Gas	0.078	15:11
05 Room Air	0.000	15:12
06 Std. Gas	0.078	15:12
07 Room Air	0.000	15:13

Lot No = 13314080A3  
Cyl No = 28  
Exp Date = 06/05/2016  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004201  
Location = TOXL      8164.13.00 06/09  
04/15/2016      15:17

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:17
02 Std. Gas	0.078	15:17
03 Room Air	0.000	15:18
04 Std. Gas	0.078	15:18
05 Room Air	0.000	15:19
06 Std. Gas	0.077	15:19
07 Room Air	0.000	15:20

Lot No = 13314080A3  
Cyl No = 28  
Exp Date = 06/05/2016  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

CALIBRATION CHECK  
0.080 AC

Form 106-I8000




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004201  
Location = TOXL      8164.13.00 06/09  
04/15/2016      15:20

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:21
02 Std. Gas	0.078	15:21
03 Room Air	0.000	15:22
04 Std. Gas	0.077	15:22
05 Room Air	0.000	15:23
06 Std. Gas	0.078	15:23
07 Room Air	0.000	15:24

Lot No = 13314080A3  
Cyl No = 28  
Exp Date = 06/05/2016  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC