

INTOXILYZER 8000 INSPECTION

Serial No.: 80-003067

Location: TOX LAB

- | | Check When Done |
|---|-----------------|
| A. Is the warm-up time less than 20 minutes? | <u>✓</u> |
| B. Is a three-pronged grounded outlet used? | <u>✓</u> |
| C. Is the breath tube heated? | <u>✓</u> |
| D. Is the diagnostic check complete? | <u>✓</u> |
| E. Is the time, date, and year correct (re-set if necessary)? | <u>✓</u> |
| F. Print test. (Attach test record.) | <u>✓</u> |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)
Sim SN <u>DR5190</u> Lot No. <u>—</u> Rep AC <u>0.019 AC</u> | <u>✓</u> |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN <u>DR5113</u> Lot No. <u>—</u> Rep AC <u>0.285</u> | <u>✓</u> |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.
(Attach test record.)
Sim SN <u>DR5135</u> Lot No. <u>—</u> Rep AC <u>0.10 AC + 0.05% Acetone</u> | <u>✓</u> |
| J. RFI check. Run CMS mode. Key radio on first room air.
(Attach test record.) | <u>✓</u> |
| K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.) | |

Lot No. <u>709352</u>	Cyl No. <u>39</u>	Exp. Date <u>12/23/2011</u>
Test 1 <u>0.079</u>	Test 4 <u>0.079</u>	Test 7 <u>0.079</u>
Test 2 <u>0.080</u>	Test 5 <u>0.079</u>	Test 8 <u>0.079</u>
Test 3 <u>0.079</u>	Test 6 <u>0.079</u>	Test 9 <u>0.079</u>
Average <u>0.079</u>		

L. Tank pressure, Level 3 Func D Sub G P: 415 psi Reg 400 psi

M. Remarks/Maintenance record: Annual Inspection

<u>16 April 2010</u> Date	<u>Margy Pearson</u> Field Inspector's Signature
<u>19 Apr 2010</u> Date	<u>Deb Keshur</u> Reviewed by

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003067
Location = TOXL 8164.13.00 06/09
04/16/2010 14:54

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:55
02 Std. Sol.	0.017	14:56
03 Room Air	0.000	14:57
04 Std. Sol.	0.018	14:57
05 Room Air	0.000	14:58
06 Std. Sol.	0.017	14:58
07 Room Air	0.000	14:59


08 Sim Temp = 34.0°C

Simul Ser No = 000000

Std Sol No = 0000

County = 08

Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Low AC 0.019AC

Form 106-I8000

*80-003067
7/1/80 KM*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003067
Location = TOXL 8164.13.00 06/09
04/16/2010 15:09

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:09
02 Std. Sol.	0.271	15:11
03 Room Air	0.000	15:11
04 Std. Sol.	0.270	15:12
05 Room Air	0.000	15:12
06 Std. Sol.	0.268	15:13
07 Room Air	0.000	15:13

08 Sim Temp = 34.0°C

Simul Ser No = 000000

Std Sol No = 0000

County = 08

Oper No. = 107501

Margy Pearson

Operator Signature
MARGY PEARSON

Remarks:

Linearity 0.285 AC

Form 106-I8000

80-003067
75/80 AH


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003067
Location = TOXL 8164.13.00 06/09
04/16/2010 15:04

Test	AC	Time
01 Room Air	0.000	15:04
02 *Subject Test	INT*	15:05
03 Room Air	0.000	15:05

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Simul Ser No = 000000
Std Sol No = 0000
County = 08 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

0.10AC + 0.05% Acetone

Form 106-I8000

80-003067
7/6/80 *km*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003067
Location = TOXL 8164.13.00 06/09
04/16/2010 15:20

Test	AC	Time
01 Diagnostic	OK	15:21
02 Room Air	RFI*	15:21
03 Room Air	0.000	15:22

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Simul Ser No = 000000
Std Sol No = 0000
County = 08 Oper No. = 107501

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
MARGY PEARSON

Remarks:

RFI Test

Form 106-I8000

80-003067
7/80 15A


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003067
Location = TOXL 8164.13.00 06/09
04/16/2010 14:41

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:42
02 Std. Gas	0.079	14:42
03 Room Air	0.000	14:43
04 Std. Gas	0.080	14:43
05 Room Air	0.000	14:44
06 Std. Gas	0.079	14:44
07 Room Air	0.000	14:44

Lot No = 709352
Cyl No = 39
Exp Date = 12/23/2011
County = 08 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Cal Check

Form 106-I8000

80-003067
7/1/10 *PL*


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003067
Location = TOXL 8164.13.00 06/09
04/16/2010 14:46

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:46
02 Std. Gas	0.079	14:47
03 Room Air	0.000	14:47
04 Std. Gas	0.079	14:47
05 Room Air	0.000	14:48
06 Std. Gas	0.079	14:48
07 Room Air	0.000	14:49

Lot No = 709352
Cyl No = 39
Exp Date = 12/23/2011
County = 08 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Cal Check

Form 106-I8000

80-003067
7/9/10 KH


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003067
Location = TOXL 8164.13.00 06/09
04/16/2010 14:50

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:50
02 Std. Gas	0.079	14:50
03 Room Air	0.000	14:51
04 Std. Gas	0.079	14:51
05 Room Air	0.000	14:52
06 Std. Gas	0.079	14:52
07 Room Air	0.000	14:53

Lot No = 709352
Cyl No = 39
Exp Date = 12/23/2011
County = 08 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Cal Check

Form 106-I8000

80-003067
80/80 192