

INTOXILYZER 8000 INSPECTION

Serial No.: 80-004191

Location: TOX LAB

- | | Check When Done |
|--|-----------------|
| A. Is the warm-up time less than 20 minutes? | <u>✓</u> |
| B. Is a three-pronged grounded outlet used? | <u>✓</u> |
| C. Is the breath tube heated? | <u>✓</u> |
| D. Is the diagnostic check complete? | <u>✓</u> |
| E. Is the time, date, and year correct (re-set if necessary)? | <u>✓</u> |
| F. Print test. (Attach test record.) | <u>✓</u> |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)
Sim SN <u>DR 5190</u> Lot No. <u>—</u> Rep AC <u>0.019</u> | <u>✓</u> |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN <u>DR 5113</u> Lot No. <u>—</u> Rep AC <u>0.285 AC</u> | <u>✓</u> |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.
(Attach test record.)
Sim SN <u>DR 5135</u> Lot No. <u>—</u> Rep AC <u>0.10 AC + 0.05% Acetone</u> | <u>✓</u> |
| J. RFI check. Run CMS mode. Key radio on first room air.
(Attach test record.) | <u>✓</u> |
| K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.) | |

Lot No. <u>659357</u>	Cyl No. <u>42</u>	Exp. Date <u>9/03/2011</u>
Test 1 <u>0.081</u>	Test 4 <u>0.081</u>	Test 7 <u>0.081</u>
Test 2 <u>0.080</u>	Test 5 <u>0.080</u>	Test 8 <u>0.081</u>
Test 3 <u>0.081</u>	Test 6 <u>0.080</u>	Test 9 <u>0.081</u>
Average <u>0.081 AC</u>		

L. Tank pressure, Level 3 Func D Sub G P: 445 psi Reg 480 psi

M. Remarks/Maintenance record: — Annual Inspection

<u>4/16/2010</u> Date	<u>Mary Pearson</u> Field Inspector's Signature
<u>19 Apr 2010</u> Date	<u>Deb Keshner</u> Reviewed by

80-004191
3/6/14 *JK*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004191
Location = TOXL 8164.13.00 06/09
04/16/2010 08:47

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 659358
Standard Cyl #? Y
Standard Expiration? 09/03/2011
Oper No? 107501

Flow Cal. Date: 01/28/2010
Slope 738
Intercept -777574

IR Calibration Date: 01/28/2010

	3um	9um
0th Coef(*100):	-23128	-33609
1st Coef(*100):	272485	131724
2nd Coef(*100):	1490	1056
H2O adj(mg/l*10k):	432	614

***** Printer Test End *****



Operator Signature
MARGY PEARSON

Remarks:

Print Test

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004191
Location = TOXL 8164.13.00 06/09
04/16/2010 09:21

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:22
02 Std. Sol.	0.016	09:23
03 Room Air	0.000	09:23
04 Std. Sol.	0.016	09:24
05 Room Air	0.000	09:25
06 Std. Sol.	0.016	09:25
07 Room Air	0.000	09:26

08 Sim Temp = 34.0°C

Simul Ser No = 00009=000
Std Sol No = 00000
County = 08 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:
Low 0.019 AC

Form 106-I8000

*80-004191
38/44 KM*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004191
Location = TOXL 8164.13.00 06/09
04/16/2010 09:32

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:32
02 Std. Sol.	0.285	09:33
03 Room Air	0.000	09:34
04 Std. Sol.	0.286	09:35
05 Room Air	0.000	09:35
06 Std. Sol.	0.286	09:36
07 Room Air	0.000	09:37

08 Sim Temp = 34.0°C

Simul Ser No = 00009=000

Std Sol No = 00000

County = 08

Oper No. = 107501

Margy Pearson

Operator Signature
MARGY PEARSON

Remarks:

Linearity 0.285 AC

Form 106-I8000

*80-004191
39/44 KA*


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004191
Location = TOXL 8164.13.00 06/09
04/16/2010 09:17

Test	AC	Time
01 Room Air	0.000	09:18
02 *Subject Test	INT*	09:19
03 Room Air	0.000	09:19

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Simul Ser No = 00009=000
Std Sol No = 00000
County = 08 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Alc + Acetone

Form 106-I8000

80-004191
4/10/44 192

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004191
Location = TOXL 8164.13.00 06/09
04/16/2010 09:07

Test	AC	Time
01 Diagnostic	OK	09:08
02 Room Air	RFI*	09:08
03 Room Air	0.000	09:08

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 659358
Cyl No = Y
Expiration Date = 09/03/2011
County = 08 Oper No. = 107501

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
MARGY PEARSON

Remarks:

RFI

Form 106-I8000

80-004191
4/14/12


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004191
Location = TOXL 8164.13.00 06/09
04/16/2010 09:44

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:45
02 Std. Gas	0.081	09:45
03 Room Air	0.000	09:46
04 Std. Gas	0.080	09:46
05 Room Air	0.000	09:47
06 Std. Gas	0.081	09:47
07 Room Air	0.000	09:48

Lot No = 659357
Cyl No = 42
Exp Date = 09/03/2011
County = 08 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Cal Check

Form 106-I8000

80-004191
4/16/10 PM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004191
Location = TOXL 8164.13.00 06/09
04/16/2010 09:48

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:49
02 Std. Gas	0.081	09:49
03 Room Air	0.000	09:49
04 Std. Gas	0.080	09:50
05 Room Air	0.000	09:50
06 Std. Gas	0.080	09:51
07 Room Air	0.000	09:51

Lot No = 659357
Cyl No = 42
Exp Date = 09/03/2011
County = 08 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Cal Check

Form 106-I8000

80-004191
4/3/11 PA


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004191
Location = TOXL 8164.13.00 06/09
04/16/2010 09:54

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:55
02 Std. Gas	0.081	09:55
03 Room Air	0.000	09:56
04 Std. Gas	0.081	09:56
05 Room Air	0.000	09:56
06 Std. Gas	0.081	09:57
07 Room Air	0.000	09:57

Lot No = 659357
Cyl No = 42
Exp Date = 09/03/2011
County = 08 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Cal Check

Form 106-18000

*80-004191
4/16/11 PM*