

INTOXILYZER 8000 INSPECTION

Serial No.: 80-004934

Location: TOX LAB

- | | Check When Done |
|--|-----------------|
| A. Is the warm-up time less than 20 minutes? | <u>✓</u> |
| B. Is a three-pronged grounded outlet used? | <u>✓</u> |
| C. Is the breath tube heated? | <u>✓</u> |
| D. Is the diagnostic check complete? | <u>✓</u> |
| E. Is the time, date, and year correct (re-set if necessary)? | <u>✓</u> |
| F. Print test. (Attach test record.) | <u>✓</u> |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)
Sim SN <u>DRS190</u> Lot No. <u>—</u> Rep AC <u>0.019</u> | <u>✓</u> |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN <u>DRS113</u> Lot No. <u>—</u> Rep AC <u>0.285 AC</u> | <u>✓</u> |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.
(Attach test record.)
Sim SN <u>DR 5135</u> Lot No. <u>—</u> Rep AC <u>0.10 AC + 0.05% Acetone</u> | <u>✓</u> |
| J. RFI check. Run CMS mode. Key radio on first room air.
(Attach test record.) | <u>✓</u> |

K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.)

Lot No. ^{MAP} 704659357 Cyl No. 11 Exp. Date 09/03/2011

Test 1 <u>0.080</u>	Test 4 <u>0.080</u>	Test 7 <u>0.080</u>
Test 2 <u>0.080</u>	Test 5 <u>0.081</u>	Test 8 <u>0.080</u>
Test 3 <u>0.080</u>	Test 6 <u>0.081</u>	Test 9 <u>0.080</u>

Average 0.080

L. Tank pressure, Level 3 Func D Sub G P: 216 psi Reg 200 psi

M. Remarks/Maintenance record: Annual Inspection

16 April 2010
Date

Mary Pearson
Field Inspector's Signature

19 Apr 2010
Date

Deb Kashner
Reviewed by

80-004934
3/31/11 *MP*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004934
Location = TOXL 8164.13.00 06/09
04/16/2010 12:17

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 659357
Standard Cyl #? 11
Standard Expiration? 09/03/2011
Oper No? 107501

Flow Cal. Date: 01/20/2010
Slope 698
Intercept -694859

IR Calibration Date: 03/12/2010
 3um 9um

0th Coef(*100): -18922 -21325
1st Coef(*100): 252748 128359
2nd Coef(*100): 509 713
H2O adj(mg/l*10k): 732 642

***** Printer Test End *****



Operator Signature
MARGY PEARSON

Remarks:

Print

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004934
Location = TOXL 8164.13.00 06/09
04/16/2010 13:58

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:58
02 Std. Sol.	0.016	14:01
03 Room Air	0.000	14:01
04 Std. Sol.	0.017	14:02
05 Room Air	0.000	14:03
06 Std. Sol.	0.018	14:03
07 Room Air	0.000	14:04


08 Sim Temp = 34.0°C

Simul Ser No = 0000000

Std Sol No = 000000

County = 08

Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Low AC 0.019AC

Form 106-I8000

80-004934
3/5/11 *AK*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004934
Location = TOXL 8164.13.00 06/09
04/16/2010 13:52

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:53
02 Std. Sol.	0.276	13:54
03 Room Air	0.000	13:54
04 Std. Sol.	0.276	13:55
05 Room Air	0.000	13:55
06 Std. Sol.	0.276	13:56
07 Room Air	0.000	13:57


08 Sim Temp = 34.0°C

Simul Ser No = 0000000

Std Sol No = 000000

County = 08

Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Linearity 0.285 AC

Form 106-18000

80-004934
3/6/11 pm


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004934
Location = TOXL 8164.13.00 06/09
04/16/2010 14:10

Test	AC	Time
01 Room Air	0.000	14:11
02 *Subject Test	INT*	14:12
03 Room Air	0.000	14:12

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Simul Ser No = 0000000
Std Sol No = 000000
County = 08 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

0.10 AC + 0.05% A/E tone

Form 106-I8000

80-004934
3/4/11 FM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004934
Location = TOXL 8164.13.00 06/09
04/16/2010 13:44

Test	AC	Time
01 Diagnostic	OK	13:45
02 Room Air	RFI*	13:45
03 Room Air	0.000	13:46

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 659357
Cyl No = 11
Expiration Date = 09/03/2011
County = 08 Oper No. = 107501

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
MARGY PEARSON

Remarks:

RFI

Form 106-I8000

80-004934
3/11 PM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004934
Location = TOXL 8164.13.00 06/09
04/16/2010 12:04

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:05
02 Std. Gas	0.080	12:05
03 Room Air	0.000	12:06
04 Std. Gas	0.080	12:06
05 Room Air	0.000	12:06
06 Std. Gas	0.080	12:07
07 Room Air	0.000	12:07

Lot No = 659357
Cyl No = 11
Exp Date = 09/03/2011
County = 08 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Cal Check

Form 106-I8000

80-004934
3/1/11 JH


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004934
Location = TOXL 8164.13.00 06/09
04/16/2010 12:09

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:10
02 Std. Gas	0.080	12:10
03 Room Air	0.000	12:11
04 Std. Gas	0.081	12:11
05 Room Air	0.000	12:11
06 Std. Gas	0.081	12:12
07 Room Air	0.000	12:12

Lot No = 659357
Cyl No = 11
Exp Date = 09/03/2011
County = 08 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Cal Check

Form 106-I8000

*80-004934
4/16/11 FM*


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004934
Location = TOXL 8164.13.00 06/09
04/16/2010 12:14

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:14
02 Std. Gas	0.080	12:15
03 Room Air	0.000	12:15
04 Std. Gas	0.080	12:16
05 Room Air	0.000	12:16
06 Std. Gas	0.080	12:16
07 Room Air	0.000	12:17

Lot No = 659357
Cyl No = 11
Exp Date = 09/03/2011
County = 08 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Cal Check

Form 106-I8000

80-004934
4/16/11 PM