

INTOXILYZER 8000 INSPECTION

Serial No.: 80-004942

Location: TOX LAB

- | | Check When Done |
|---|-----------------|
| A. Is the warm-up time less than 20 minutes? | <u>✓</u> |
| B. Is a three-pronged grounded outlet used? | <u>✓</u> |
| C. Is the breath tube heated? | <u>✓</u> |
| D. Is the diagnostic check complete? | <u>✓</u> |
| E. Is the time, date, and year correct (re-set if necessary)? | <u>✓</u> |
| F. Print test. (Attach test record.) | <u>✓</u> |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)
Sim SN <u>DR5190</u> Lot No. <u>—</u> Rep AC <u>0.019</u> | <u>✓</u> |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN <u>DR51</u> Lot No. <u>—</u> Rep AC <u>0.285</u> | <u>✓</u> |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.
(Attach test record.)
Sim SN <u>DR5135</u> Lot No. <u>—</u> Rep AC <u>0.10 AC + 0.05% Acetone</u> | <u>✓</u> |
| J. RFI check. Run CMS mode. Key radio on first room air.
(Attach test record.) | <u>✓</u> |
| K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.) | |

Lot No. <u>709352</u>	Cyl No. <u>26</u>	Exp. Date <u>12/23/2011</u>
Test 1 <u>0.080</u>	Test 4 <u>0.079</u>	Test 7 <u>0.080</u>
Test 2 <u>0.080</u>	Test 5 <u>0.080</u>	Test 8 <u>0.080</u>
Test 3 <u>0.080</u>	Test 6 <u>0.080</u>	Test 9 <u>0.081</u>
Average <u>0.080 AC</u>		

L. Tank pressure, Level 3 Func D Sub G P: 599 psi Reg 650 psi

M. Remarks/Maintenance record: Annual Inspection

16 April 2010
Date

Margy Pleuro
Field Inspector's Signature

19 Apr 2010
Date

Deb Kashner
Reviewed by

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004942
Location = TOXL 8164.13.00 06/09
04/16/2010 13:49

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:50
02 Std. Sol.	0.018	13:51
03 Room Air	0.000	13:51
04 Std. Sol.	0.017	13:52
05 Room Air	0.000	13:53
06 Std. Sol.	0.018	13:53
07 Room Air	0.000	13:54


08 Sim Temp = 34.0°C

Simul Ser No = 0000000

Std Sol No = 000000

County = 08

Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Low AC 0.019AC

Form 106-I8000

80-004942
3/10 FM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004942
Location = TOXL 8164.13.00 06/09
04/16/2010 13:58

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:58
02 Std. Sol.	0.283	13:59
03 Room Air	0.000	14:00
04 Std. Sol.	0.283	14:00
05 Room Air	0.000	14:01
06 Std. Sol.	0.283	14:02
07 Room Air	0.000	14:02

08 Sim Temp = 34.0°C

Simul Ser No = 0000000
Std Sol No = 000000
County = 08 Oper No. = 107501


Operator Signature
MARGY PEARSON

Remarks:

Linearity 0.285 AC

Form 106-I8000

80-004942
35/40 KAL


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004942
Location = TOXL 8164.13.00 06/09
04/16/2010 13:54

Test	AC	Time
01 Room Air	0.000	13:55
02 *Subject Test	INT*	13:56
03 Room Air	0.000	13:56

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Simul Ser No = 0000000
Std Sol No = 000000
County = 08 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

0.10 AC + 0.05% Acetone

Form 106-I8000

80-004942
3/10/10

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004942
Location = TOXL 8164.13.00 06/09
04/16/2010 13:42

Test	AC	Time
01 Diagnostic	OK	13:43
02 Room Air	RFI*	13:43
03 Room Air	0.000	13:43

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 709352
Cyl No = 26
Expiration Date = 12/23/2011
County = 08 Oper No. = 107501

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
MARGY PEARSON

Remarks:

RFI

Form 106-I8000

80-004942
37/10 KM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004942
Location = TOXL 8164.13.00 06/09
04/16/2010 12:03

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:04
02 Std. Gas	0.080	12:04
03 Room Air	0.000	12:05
04 Std. Gas	0.080	12:05
05 Room Air	0.000	12:05
06 Std. Gas	0.080	12:06
07 Room Air	0.000	12:06

Lot No = 709352
Cyl No = 26
Exp Date = 12/23/2011
County = 08 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Cal Check

Form 106-I8000

80-004942
38/40 FM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004942
Location = TOXL 8164.13.00 06/09
04/16/2010 12:06

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:07
02 Std. Gas	0.079	12:07
03 Room Air	0.000	12:08
04 Std. Gas	0.080	12:08
05 Room Air	0.000	12:09
06 Std. Gas	0.080	12:09
07 Room Air	0.000	12:10

Lot No = 709352
Cyl No = 26
Exp Date = 12/23/2011
County = 08 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Cal Check

Form 106-I8000

*80-004942
39/40 PK*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004942
Location = TOXL 8164.13.00 06/09
04/16/2010 12:11

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:11
02 Std. Gas	0.080	12:12
03 Room Air	0.000	12:12
04 Std. Gas	0.080	12:13
05 Room Air	0.000	12:13
06 Std. Gas	0.081	12:14
07 Room Air	0.000	12:14

Lot No = 709352
Cyl No = 26
Exp Date = 12/23/2011
County = 08 Oper No. = 107501

Margy Pearson

Operator Signature
MARGY PEARSON

Remarks:

Cal Check

Form 106-I8000

80-004942
4/16/10 FM