

INTOXILYZER 8000 INSPECTION

Serial No.: 80-005360

Location: ToxL

Check When Done

- A. Is the warm-up time less than 20 minutes? ✓
- B. Is a three-pronged grounded outlet used? ✓
- C. Is the breath tube heated? ✓
- D. Is the diagnostic check complete? ✓
- E. Is the time, date, and year correct (re-set if necessary)? ✓
- F. Print test. (Attach test record.) ✓
- G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)  
Sim SN DR3455 Lot No. 10290 Rep AC 0.020 ✓
- H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)  
Sim SN DR3847 Lot No. 11110 Rep AC 0.300 ✓
- I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.  
(Attach test record.)  
Sim SN DR3377 Lot No.      Rep AC 0.100 ✓
- J. RFI check. Run CMS mode. Key radio on first room air.  
(Attach test record.) ✓
- K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.  
(Attach test record.)

Lot No. 2341108041 Cyl No. 55 Exp. Date 10/01/2013  
Test 1 0.083 Test 4 0.081 Test 7 0.082  
Test 2 0.083<sup>CEE</sup>2 Test 5 0.082 Test 8 0.082  
Test 3 0.083 Test 6 0.081 Test 9 0.081  
Average 0.082

L. Tank pressure, Level 3 Func D Sub G P: 840 psi Reg 850 psi

M. Remarks/Maintenance record: ANNUAL INSPECTION  
16 APR 12 Date Charles E. Ed Field Inspector's Signature  
17 Apr 2012 Date Deb Kashner Reviewed by

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-005360  
Location = TOXL                    8164.13.00 06/09  
04/16/2012    16:36

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?          1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            ~~659358~~ 23411080A1  
Standard Cyl #?            ~~2755~~  
Standard Expiration?        ~~09/03/2011~~ 10/01/2013  
Oper No?                    666666

Flow Cal. Date:            08/17/2011  
Slope                      639  
Intercept                  -343114

IR Calibration Date:        08/17/2011  
                                    3um                      9um  
-----  
0th Coef(\*100):            -17266                      -25747  
1st Coef(\*100):            263773                      142804  
2nd Coef(\*100):            2606                        1406  
H2O adj(mg/l\*10k):        499                         577

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:

PRINT TEST

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005360  
Location = TOXL      8164.13.00 06/09  
04/16/2012      16:55

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:55
02 Std. Sol.	0.018	16:56
03 Room Air	0.000	16:57
04 Std. Sol.	0.018	16:57
05 Room Air	0.000	16:58
06 Std. Sol.	0.019	16:59
07 Room Air	0.000	16:59

08 Sim Temp = 34.0°C

Simul Ser No = 000000  
Std Sol No = 0000  
County = 08      Oper No. = 666666

  
\_\_\_\_\_  
Operator Signature  
CHARLES EDER

Remarks:      *LOW AC*  
  
   *0.020 AC*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005360  
Location = TOXL      8164.13.00 06/09  
04/16/2012      16:49

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:50
02 Std. Sol.	0.302	16:50
03 Room Air	0.000	16:51
04 Std. Sol.	0.303	16:52
05 Room Air	0.000	16:52
06 Std. Sol.	0.303	16:53
07 Room Air	0.000	16:54

08 Sim Temp = 34.0°C

Simul Ser No = 000000

Std Sol No = 0000

County = 08

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

LINEARITY TEST  
0.300 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005360  
Location = TOXL      8164.13.00 06/09  
04/16/2012      17:00

Test	AC	Time
01 Room Air	0.000	17:02
02 *Subject Test	INT*	17:02
03 Room Air	0.000	17:03

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = 123456  
Dr. Lic. = ND/DIS821456  
Simul Ser No = 000000  
Std Sol No = 0000  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *INTERFERENT CHECK*  
*0.100 AC + ACETONE*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005360  
Location = TOXL      8164.13.00 06/09  
04/16/2012      17:04

Test	AC	Time
01 Diagnostic	OK	17:05
02 Room Air	RFI*	17:05
03 Room Air	0.000	17:05

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = 123456  
Dr. Lic. = ND/DIS821456  
Simul Ser No = 000000  
Std Sol No = 0000  
County = 08      Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks: *RFI CHECK*


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005360  
Location = TOXL      8164.13.00 06/09  
04/16/2012      17:07

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:08
02 Std. Gas	0.083	17:08
03 Room Air	0.000	17:09
04 Std. Gas	0.082	17:09
05 Room Air	0.000	17:09
06 Std. Gas	0.083	17:10
07 Room Air	0.000	17:10

Lot No = 23411080A1  
Cyl No = 55  
Exp Date = 10/01/2013  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005360  
Location = TOXL      8164.13.00 06/09  
04/16/2012      17:11

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:12
02 Std. Gas	0.081	17:12
03 Room Air	0.000	17:12
04 Std. Gas	0.082	17:13
05 Room Air	0.000	17:13
06 Std. Gas	0.081	17:14
07 Room Air	0.000	17:14

Lot No = 23411080A1  
Cyl No = 55  
Exp Date = 10/01/2013  
County = 08      Oper No. = 666666

  
Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005360  
Location = TOXL      8164.13.00 06/09  
04/16/2012      17:14

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:15
02 Std. Gas	0.082	17:15
03 Room Air	0.000	17:16
04 Std. Gas	0.082	17:16
05 Room Air	0.000	17:17
06 Std. Gas	0.081	17:17
07 Room Air	0.000	17:18

Lot No = 23411080A1  
Cyl No = 55  
Exp Date = 10/01/2013  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000