

INTOXILYZER 8000 INSPECTION

Serial No.: 80-005363

Location: ToxL

- | | Check When Done |
|--|-----------------|
| A. Is the warm-up time less than 20 minutes? | <u>✓</u> |
| B. Is a three-pronged grounded outlet used? | <u>✓</u> |
| C. Is the breath tube heated? | <u>✓</u> |
| D. Is the diagnostic check complete? | <u>✓</u> |
| E. Is the time, date, and year correct (re-set if necessary)? | <u>✓</u> |
| F. Print test. (Attach test record.) | <u>✓</u> |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.) Sim SN <u>DR3455</u> Lot No. <u>10290</u> Rep AC <u>0.020</u> | <u>✓</u> |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.) Sim SN <u>DR3847</u> Lot No. <u>1110</u> Rep AC <u>0.300</u> | <u>✓</u> |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode. (Attach test record.) Sim SN <u>DR3377</u> Lot No. <u> </u> Rep AC <u>0.100</u> | <u>✓</u> |
| J. RFI check. Run CMS mode. Key radio on first room air. (Attach test record.) | <u>✓</u> |

K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.)

| | | |
|---------------------------|---------------------|-----------------------------|
| Lot No. <u>23411080A1</u> | Cyl No. <u>55</u> | Exp. Date <u>10/01/2013</u> |
| Test 1 <u>0.082</u> | Test 4 <u>0.082</u> | Test 7 <u>0.082</u> |
| Test 2 <u>0.081</u> | Test 5 <u>0.081</u> | Test 8 <u>0.082</u> |
| Test 3 <u>0.081</u> | Test 6 <u>0.081</u> | Test 9 <u>0.082</u> |
| Average <u>0.082</u> | | |

L. Tank pressure, Level 3 Func D Sub G P: 851 psi Reg 850 psi

M. Remarks/Maintenance record: ANNUAL INSPECTION

16 APR 12
Date
17 Apr 2012
Date

[Signature]
Field Inspector's Signature
Deb Kashner
Reviewed by

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005363
Location = TOXL 8164.13.00 06/09
04/16/2012 10:45

***** Printer Test *****

abcdefghijklmnopqrstuvwxy1234567890-=
ABCDEFGHIJKLMNORSTUVWXYZ!@#%^&*()_+?

abcdefghijklmnopqrstuvwxy1234567890-=
ABCDEFGHIJKLMNORSTUVWXYZ!@#%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 23411080A1
Standard Cyl #? 55
Standard Expiration? 10/01/2013
Oper No? 666666

Flow Cal. Date: 08/19/2011
Slope 660
Intercept -560671

IR Calibration Date: 08/19/2011
3um 9um

0th Coef(*100): -22812 -27781
1st Coef(*100): 271435 139375
2nd Coef(*100): 3844 1471
H2O adj(mg/l*10k): 585 583

***** Printer Test End *****



Operator Signature
CHARLES EDER

Remarks:

PRINT TEST

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005363
Location = TOXL 8164.13.00 06/09
04/16/2012 11:06

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 11:07 |
| 02 Std. Sol. | 0.018 | 11:08 |
| 03 Room Air | 0.000 | 11:08 |
| 04 Std. Sol. | 0.018 | 11:09 |
| 05 Room Air | 0.000 | 11:09 |
| 06 Std. Sol. | 0.019 | 11:10 |
| 07 Room Air | 0.000 | 11:11 |

08 Sim Temp = 34.0°C

Simul Ser No = 000000

Std Sol No = 0000

County = 08

Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

Low AC 0.020 AC

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005363
Location = TOXL 8164.13.00 06/09
04/16/2012 11:18

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 11:19 |
| 02 Std. Sol. | 0.304 | 11:20 |
| 03 Room Air | 0.000 | 11:20 |
| 04 Std. Sol. | 0.304 | 11:21 |
| 05 Room Air | 0.000 | 11:22 |
| 06 Std. Sol. | 0.303 | 11:22 |
| 07 Room Air | 0.000 | 11:23 |

08 Sim Temp = 34.0°C

Simul Ser No = 000000

Std Sol No = 0000

County = 08

Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

LINEARITY TEST
0.300 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005363
Location = TOXL 8164.13.00 06/09
04/16/2012 11:01

| Test | AC | Time |
|------------------|-------|-------|
| 01 Room Air | 0.000 | 11:02 |
| 02 *Subject Test | INT* | 11:03 |
| 03 Room Air | 0.000 | 11:03 |

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = 1223456
Dr. Lic. = ND/DIS821456
Lot No = 23411080A1
Cyl No = 55
Expiration Date = 10/01/2013
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *INTERFERENT CHECK*
0.100 AC + ACETONE

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005363
Location = TOXL 8164.13.00 06/09
04/16/2012 10:46

| Test | AC | Time |
|---------------|-------|-------|
| 01 Diagnostic | OK | 10:47 |
| 02 Room Air | RFI* | 10:47 |
| 03 Room Air | 0.000 | 10:47 |

*Invalid Test
 Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = 123456
Dr. Lic. = ND/DIS821456
Lot No = 23411080A1
Cyl No = 55
Expiration Date = 10/01/2013
County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks: *RFI CHECK*


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005363
Location = TOXL 8164.13.00 06/09
04/16/2012 10:48

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 10:49 |
| 02 Std. Gas | 0.082 | 10:49 |
| 03 Room Air | 0.000 | 10:50 |
| 04 Std. Gas | 0.081 | 10:50 |
| 05 Room Air | 0.000 | 10:51 |
| 06 Std. Gas | 0.081 | 10:51 |
| 07 Room Air | 0.000 | 10:51 |

Lot No = 23411080A1
Cyl No = 55
Exp Date = 10/01/2013
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005363
Location = TOXL 8164.13.00 06/09
04/16/2012 10:52

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 10:52 |
| 02 Std. Gas | 0.082 | 10:53 |
| 03 Room Air | 0.000 | 10:53 |
| 04 Std. Gas | 0.081 | 10:54 |
| 05 Room Air | 0.000 | 10:54 |
| 06 Std. Gas | 0.081 | 10:55 |
| 07 Room Air | 0.000 | 10:55 |

Lot No = 23411080A1

Cyl No = 55

Exp Date = 10/01/2013

County = 08

Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

CALIBRATION CHECK

0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005363
Location = TOXL 8164.13.00 06/09
04/16/2012 10:55

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 10:56 |
| 02 Std. Gas | 0.082 | 10:56 |
| 03 Room Air | 0.000 | 10:57 |
| 04 Std. Gas | 0.082 | 10:57 |
| 05 Room Air | 0.000 | 10:58 |
| 06 Std. Gas | 0.082 | 10:58 |
| 07 Room Air | 0.000 | 10:59 |

Lot No = 23411080A1
Cyl No = 55
Exp Date = 10/01/2013
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

CALIBRATION CHECK
0.080 AC

Form 106-I8000