

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-004944

Location: West Fargo PD

Check When Done

- A. Is the breath tube heated? ✓
- B. Is the diagnostic check complete? ✓
- C. Is the time, date, and year correct (re-set if necessary)? ✓
- D. Print test. (Attach test record.) ✓
- E. Low AC. Use  $\leq 0.03$  AC in ACA mode. (Attach test record.)  
Sim SN DR5143 Lot No. 201501 B Rep AC 0.020 ✓
- F. High AC. Use  $\geq 0.25$  AC in ACA mode. (Attach test record.)  
Sim SN DR3455 Lot No. 14270 Rep AC 0.300 ✓
- G. Interferent check. Use 0.05% acetone in ethanol solution in ABA mode.  
(Attach test record.)  
Sim SN DR5132 Lot No. ICS#2 Solution 0.10 AC + 0.05% acetone ✓
- H. RFI check. Run CMS mode. Key radio during CMS test.  
(Attach test record.) ✓
- I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.  
(Attach test record.) 0.080 AC

Lot No. 14214080A6 Cyl No. 59 Exp. Date 07/05/2016  
Test 1 0.082 Test 4 0.082 Test 7 0.082  
Test 2 0.082 Test 5 0.082 Test 8 0.082  
Test 3 0.082 Test 6 0.083 Test 9 0.082  
Average 0.082

J. Tank pressure, Level 3 Func D Sub G P: 609 psi Reg 600 psi

K. Remarks/Maintenance record: Annual inspection, replaced simulator return  
O-ring, checked breath hose screen - good.

15 Apr 15  
Date

17 APR 15  
Date

Deb Kashur  
Field Inspector's Signature

Charles E. [Signature] ✓  
Reviewed by

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004944  
Location = WFPD                    8164.13.00 06/09  
04/15/2015                            15:06

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%&^\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%&^\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?          1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            14214080A6  
Standard Cyl #?            59  
Standard Expiration?        07/05/2016  
Oper No?                    777777

Flow Cal. Date:            07/01/2011  
Slope                        685  
Intercept                   -634740

IR Calibration Date:        07/01/2011  
                                  3um                            9um

-----  
0th Coef(\*100):            -22591                        -17234  
1st Coef(\*100):            265915                        132691  
2nd Coef(\*100):            2841                            1562  
H2O adj(mg/l\*10k):        559                            425

\*\*\*\*\* Printer Test End \*\*\*\*\*

  
-----  
Operator Signature  
DEB KASHUR

Remarks:

Print test

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004944  
Location = WFPD                    8164.13.00 06/09  
04/15/2015                            15:24

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:24
02 Std. Sol.	0.021	15:25
03 Room Air	0.000	15:26
04 Std. Sol.	0.021	15:27
05 Room Air	0.000	15:27
06 Std. Sol.	0.021	15:28
07 Room Air	0.000	15:28

08 Sim Temp = 34.0°C

Simul Ser No = DR5143  
Std Sol No = 201501B  
County = 09                            Oper No. = 777777

Deb Kasher  
Operator Signature  
DEB KASHUR

Remarks: <sup>DIC</sup>  
Calib Low AC <sup>pk 0.020</sup> ~~0.021~~ AC

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004944  
Location = WFPD                    8164.13.00 06/09  
04/15/2015                            15:38

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:38
02 Std. Sol.	0.311	15:39
03 Room Air	0.000	15:40
04 Std. Sol.	0.310	15:41
05 Room Air	0.000	15:41
06 Std. Sol.	0.309	15:42
07 Room Air	0.000	15:42

08 Sim Temp = 34.0°C

Simul Ser No = DR3455  
Std Sol No = 14270  
County = 09                            Oper No. = 777777

Deb Kasher  
Operator Signature  
DEB KASHUR

Remarks:

High AC 0.300 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004944  
Location = WFPD                    8164.13.00 06/09  
04/15/2015                            15:44

Test	AC	Time
01 Room Air	0.000	15:45
02 *Subject Test	INT*	15:45
03 Room Air	0.000	15:46

\*Invalid Test  
Interferent Detected

Sub Name = TEST, NA NA  
Sub DOB = 08/13/1980  
Sub Sex = Unknown                    Weight = NA  
Test = DUI                              Cit = NA  
Dr. Lic. = ND/NA  
Lot No = 14214080A6  
Cyl No = 59  
Expiration Date = 07/05/2016  
County = 09                              Oper No. = 777777



Operator Signature  
DEB KASHUR

Remarks:

*Interferent check 0.10AC + 0.05% acetone*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004944  
Location = WFPD                    8164.13.00 06/09  
04/15/2015                            15:46

Test	AC	Time
01 Diagnostic	OK	15:47
02 Room Air	RFI*	15:47
03 Room Air	0.000	15:48

\*Invalid Test  
Inhibited - RFI

Sub Name = RFI, NA NA  
Sub DOB = 08/13/1980  
Sub Sex = Unknown                    Weight = NA  
Test = DUI                              Cit = NA  
Dr. Lic. = ND/NA  
Lot No = 14214080A6  
Cyl No = 59  
Expiration Date = 07/05/2016  
County = 09                              Oper No. = 777777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
DEB KASHUR

Remarks:

RFI check

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004944  
Location = WFPD      8164.13.00 06/09  
04/15/2015      15:11

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:11
02 Std. Gas	0.082	15:12
03 Room Air	0.000	15:12
04 Std. Gas	0.082	15:13
05 Room Air	0.000	15:13
06 Std. Gas	0.082	15:14
07 Room Air	0.000	15:14

Lot No = 14214080A6  
Cyl No = 59  
Exp Date = 07/05/2016  
County = 09      Oper No. = 777777

Deb Kashur  
Operator Signature  
DEB KASHUR

Remarks:

Calibration check 0.080 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004944  
Location = WFPD      8164.13.00 06/09  
04/15/2015      15:15

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:16
02 Std. Gas	0.082	15:16
03 Room Air	0.000	15:17
04 Std. Gas	0.082	15:17
05 Room Air	0.000	15:17
06 Std. Gas	0.083	15:18
07 Room Air	0.000	15:18

Lot No = 14214080A6  
Cyl No = 59  
Exp Date = 07/05/2016  
County = 09      Oper No. = 777777

  
\_\_\_\_\_  
Operator Signature  
DEB KASHUR

Remarks:

Calibration check 0.080 AC



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004944  
Location = WFPD      8164.13.00 06/09  
04/15/2015      15:20

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:20
02 Std. Gas	0.082	15:21
03 Room Air	0.000	15:21
04 Std. Gas	0.082	15:22
05 Room Air	0.000	15:22
06 Std. Gas	0.082	15:23
07 Room Air	0.000	15:23

Lot No = 14214080A6  
Cyl No = 59  
Exp Date = 07/05/2016  
County = 09      Oper No. = 777777

Deb Kashur  
Operator Signature  
DEB KASHUR

Remarks:

Calibration check 0.080 AC