



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59281 (5-2011)

Serial Number <i>80-005946</i>	Instrument Location <i>Watford City Police Dept.</i>
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level 1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. RFI Test (Use CMS Mode; Key Radio During First Room Air).
- 7. Repair and/or Maintenance Performed (if Any): *None Needed*
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>[Signature]</i>	Date <i>04.24.2014</i>
Reviewed By (Crime Laboratory Use Only) <i>Deb Kashner</i>	Date <i>25 Apr 2014 ✓</i>

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005946
Location = WCPD 8164.13.00 06/09
04/24/2014 11:25

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%&^*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%&^*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 16512080A2
Standard Cyl #? 64
Standard Expiration? 08/01/2014
Oper No? 120701

Flow Cal. Date: 11/01/2013
Slope 628
Intercept -526243

IR Calibration Date: 04/01/2014
 3um 9um

0th Coef(*100): -43472 -36988
1st Coef(*100): 280525 139494
2nd Coef(*100): 2586 750
H2O adj(mg/l*10k): 846 658

***** Printer Test End *****



Operator Signature
DANIEL BLOOD

Remarks: Print Test


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005946
Location = WCPD 8164.13.00 06/09
04/24/2014 11:26

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:27
02 Std. Gas	0.077	11:27
03 Room Air	0.000	11:28
04 Std. Gas	0.077	11:28
05 Room Air	0.000	11:29
06 Std. Gas	0.078	11:29
07 Room Air	0.000	11:30

Lot No = 16512080A2
Cyl No = 64
Exp Date = 08/01/2014
County = 27 Oper No. = 120701



Operator Signature
DANIEL BLOOD

Remarks: ACA Test

Intoxilyzer Test Record and Checklist
 NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
 North Dakota Model 8000 SN 80-005946
 Location = WCPD 8164.13.00 06/09
 04/24/2014 11:30

Test	AC	Time
01 Diagnostic	OK	11:32
02 Room Air	RFI*	11:33
03 Room Air	0.000	11:33

*Invalid Test
 Inhibited - RFI

Sub Name = TEST, NA NA
 Sub DOB = 04/23/2014
 Sub Sex = Unknown Weight = 000
 Test = DUI Cit = NA
 Dr. Lic. = ND/NA
 Lot No = 16512080A2
 Cyl No = 64
 Expiration Date = 08/01/2014
 County = 27 Oper No. = 120701

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



 Operator Signature
 DANIEL BLOOD

Remarks: RFI Test (CMS mode)