

INTOXILYZER 8000 INSPECTION

Serial No.: 80-003061

Location: TOXL

- | | Check When Done |
|---|-----------------|
| A. Is the warm-up time less than 20 minutes? | <u>✓</u> |
| B. Is a three-pronged grounded outlet used? | <u>✓</u> |
| C. Is the breath tube heated? | <u>✓</u> |
| D. Is the diagnostic check complete? | <u>✓</u> |
| E. Is the time, date, and year correct (re-set if necessary)? | <u>✓</u> |
| F. Print test. (Attach test record.) | <u>✓</u> |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)
Sim SN <u>DR5132</u> Lot No. <u>—</u> Rep AC <u>0.021</u> | <u>✓</u> |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN <u>DR5134</u> Lot No. <u>—</u> Rep AC <u>0.294</u> | <u>✓</u> |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.
(Attach test record.)
Sim SN <u>DR3846</u> Lot No. <u>—</u> Rep AC <u>0.10 AC + acetone</u> | <u>✓</u> |
| J. RFI check. Run CMS mode. Key radio on first room air.
(Attach test record.) | <u>✓</u> |
| K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.) | |

Lot No. <u>15410080A1</u>	Cyl No. <u>19</u>	Exp. Date <u>07/01/2012</u>
Test 1 <u>0.079</u>	Test 4 <u>0.080</u>	Test 7 <u>0.081</u>
Test 2 <u>0.080</u>	Test 5 <u>0.081</u>	Test 8 <u>0.080</u>
Test 3 <u>0.080</u>	Test 6 <u>0.081</u>	Test 9 <u>0.081</u>
Average <u>0.080</u>		

L. Tank pressure, Level 3 Func D Sub G P: 210 psi Reg 220 psi

M. Remarks/Maintenance record: Annual inspection & calibration

27 Apr 2011
Date

Deb Kashur
Field Inspector's Signature

4-27-11
Date

[Signature]
Reviewed by

80-003061
P. 119 of 129
DFM

TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-003061
 04/26/2011 08:54:32

Auto Calibration

pg 1 of 2

```

  <<<<<      3um      >>>>>      <<<<<      9um      >>>>>
  -----
  Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1
  Sample      % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
  Sample #1   0.0290      (-0.0290)        0.0170      (0.0020)
  Sample #2   0.0150      (-0.0150)       -0.0020      (-0.0020)
  Sample #3   0.0180      (-0.0160)        0.0040      (0.0030)
  Sample #4   0.0170      (-0.0130)        0.0210      (0.0010)
  Avg % Abs   0.0167      (-0.0147)        0.0077      (0.0007)
  STD DEV     0.0015      (0.0015)         0.0119      (0.0025)
  REL STD DEV 9.165       (10.415)         155.613     (377.492)
  
```

```

  -----
  Solution = 0.021 g/210L or 0.1000 mg/l, Samples = 4, Discarded = 1
  Sample      % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
  Sample #1   0.4260      (-0.0020)        0.7810      (0.0120)
  Sample #2   0.4530      (0.0110)         0.8070      (0.0040)
  Sample #3   0.4500      (0.0260)         0.8050      (0.0040)
  Sample #4   0.4110      (0.0370)         0.7880      (0.0180)
  Avg % Abs   0.4380      (0.0247)         0.8000      (0.0087)
  STD DEV     0.0234      (0.0131)         0.0104      (0.0081)
  REL STD DEV 5.349       (52.910)         1.305       (93.264)
  
```

```

  -----
  Solution = 0.040 g/210L or 0.1905 mg/l, Samples = 4, Discarded = 1
  Sample      % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
  Sample #1   0.7390      (0.0000)         1.3600      (0.0170)
  Sample #2   0.7240      (0.0150)         1.4050      (0.0030)
  Sample #3   0.7640      (0.0150)         1.4130      (0.0120)
  Sample #4   0.7710      (0.0150)         1.4080      (0.0250)
  Avg % Abs   0.7530      (0.0150)         1.4087      (0.0133)
  STD DEV     0.0254      (0.0000)         0.0040      (0.0111)
  REL STD DEV 3.368       (0.000)          0.287       (82.953)
  
```

```

  -----
  Solution = 0.120 g/210L or 0.5714 mg/l, Samples = 4, Discarded = 1
  Sample      % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
  Sample #1   2.1800      (-0.0270)        4.0490      (-0.0220)
  Sample #2   2.1450      (-0.0040)        4.0300      (0.0000)
  Sample #3   2.1180      (0.0230)         4.0670      (-0.0160)
  Sample #4   2.1280      (0.0200)         4.0630      (0.0000)
  Avg % Abs   2.1303      (0.0130)         4.0533      (-0.0053)
  STD DEV     0.0137      (0.0148)         0.0203      (0.0092)
  REL STD DEV 0.641       (113.836)        0.501       (173.205)
  
```

```

  -----
  Solution = 0.294 g/210L or 1.4000 mg/l, Samples = 4, Discarded = 1
  Sample      % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
  Sample #1   5.0890      (-0.0430)        9.4310      (-0.0130)
  Sample #2   5.0700      (0.0090)         9.4550      (0.0320)
  Sample #3   5.0490      (0.0280)         9.4540      (0.0410)
  Sample #4   5.0210      (0.0400)         9.4410      (0.0540)
  Avg % Abs   5.0467      (0.0257)         9.4500      (0.0423)
  STD DEV     0.0246      (0.0156)         0.0078      (0.0111)
  REL STD DEV 0.487       (60.901)         0.083       (26.127)
  
```

*Calibration
 Deb Kashur*

80-003061
 P.120-4/29
 FM

TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-003061
 04/26/2011 08:54:32

Auto Calibration

pg 2 of 2

<<<<< 3um >>>>>

 Zero Order Coef -101.39
 First Order Coef 2659.84
 Second Order Coef 26.78

<<<<< 9um >>>>>

 Zero Order Coef -47.93
 First Order Coef 1366.08
 Second Order Coef 12.77

Act (g/210L)	Fit (g/210L)	Residual (g/210L)
0.000	-0.001	0.0012
0.021	0.022	-0.0014
0.040	0.040	-0.0002
0.120	0.119	0.0006
0.294	0.294	-0.0001

Act (g/210L)	Fit (g/210L)	Residual (g/210L)
0.000	-0.001	0.0008
0.021	0.022	-0.0011
0.040	0.040	0.0001
0.120	0.120	0.0003
0.294	0.294	-0.0000

<<<<< 3um >>>>>

 Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1
 Sample

<<<<< 9um >>>>>

 Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1
 Sample

Sample #1	3304.00
Sample #2	3309.00
Sample #3	3327.00
Sample #4	3219.00
Avg	3285.0000
STD DEV	57.8619
REL STD DEV	1.761
H2O adjust (mg/l*10k)	524

Sample #1	3389.00
Sample #2	3369.00
Sample #3	3365.00
Sample #4	3394.00
Avg	3376.0000
STD DEV	15.7162
REL STD DEV	0.466
H2O adjust (mg/l*10k)	433

Atmospheric Pressure = 1004

*****CALIBRATION SUCCESSFUL*****

*Calibration
 Deb Kashur*

80-003061
 P. 12/12/29
 KM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003061
Location = TOXL 8164.13.00 06/09
04/27/2011 09:13

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 15410080A1
Standard Cyl #? 19
Standard Expiration? 07/01/2012
Oper No? 777777

Flow Cal. Date: 03/04/2011
Slope 700
Intercept -623726

IR Calibration Date: 04/26/2011
 3um 9um

0th Coef(*100): -10139 -4793
1st Coef(*100): 265984 136608
2nd Coef(*100): 2678 1277
H2O adj(mg/l*10k): 524 433

***** Printer Test End *****



Operator Signature
DEB KASHUR

Remarks:

Print test

Form 106-I8000

80-003061
P.122 of 129
KFM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003061
Location = TOXL 8164.13.00 06/09
04/27/2011 08:44

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	08:45
02 Std. Sol.	0.022	08:45
03 Room Air	0.000	08:46
04 Std. Sol.	0.021	08:47
05 Room Air	0.000	08:47
06 Std. Sol.	0.022	08:48
07 Room Air	0.000	08:48

08 Sim Temp = 34.0°C

Simul Ser No = 0

Std Sol No = 0

County = 08

Oper No. = 777777

Deb Kasher
Operator Signature
DEB KASHUR

Remarks:

Low AC 0.021 AC

Form 106-I8000

80-003061
P. 123 of 129
AM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003061
Location = TOXL 8164.13.00 06/09
04/27/2011 08:58

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	08:59
02 Std. Sol.	0.292	09:00
03 Room Air	0.000	09:00
04 Std. Sol.	0.294	09:01
05 Room Air	0.000	09:01
06 Std. Sol.	0.294	09:02
07 Room Air	0.000	09:03

08 Sim Temp = 34.0°C

Simul Ser No = 0

Std Sol No = 0

County = 08

Oper No. = 777777

Deb Kasher
Operator Signature
DEB KASHUR

Remarks:

Linearity test 0.294 AC

Form 106-I8000

80-003061
P.124 of 129
KML

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003061
Location = TOXL 8164.13.00 06/09
04/27/2011 09:04

Test	AC	Time
01 Room Air	0.000	09:05
02 *Subject Test	INT*	09:05
03 Room Air	0.000	09:06

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 15410080A1
Cyl No = 19
Expiration Date = 07/01/2012
County = 08 Oper No. = 777777

Deb Kasher
Operator Signature
DEB KASHUR

Remarks:

Interferent check 0.10AC + acetone

Form 106-I8000

80-003061
p. 125 of 129
JFM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003061
Location = TOXL 8164.13.00 06/09
04/27/2011 09:14

Test	AC	Time
01 Diagnostic	OK	09:15
02 Room Air	RFI*	09:15
03 Room Air	0.000	09:15

*Invalid Test
 Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 15410080A1
Cyl No = 19
Expiration Date = 07/01/2012
County = 08 Oper No. = 777777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

RFI test

Form 106-I8000

80-003061
P. 126 of 129
FM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003061
Location = TOXL 8164.13.00 06/09
04/26/2011 09:41

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:42
02 Std. Gas	0.079	09:42
03 Room Air	0.000	09:43
04 Std. Gas	0.080	09:43
05 Room Air	0.000	09:44
06 Std. Gas	0.080	09:44
07 Room Air	0.000	09:44

Lot No = 15410080A1
Cyl No = 19
Exp Date = 07/01/2012
County = 08 Oper No. = 777777

Deb Kashur
Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080AC

Form 106-I8000

80-003061
P.127 4/129
JFM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003061
Location = TOXL 8164.13.00 06/09
04/26/2011 10:11

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:11
02 Std. Gas	0.080	10:12
03 Room Air	0.000	10:12
04 Std. Gas	0.081	10:13
05 Room Air	0.000	10:13
06 Std. Gas	0.081	10:13
07 Room Air	0.000	10:14

Lot No = 15410080A1
Cyl No = 19
Exp Date = 07/01/2012
County = 08 Oper No. = 777777

Deb Kasher
Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080 AC

Form 106-I8000

80-003061
P. 128 of 129
YML

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003061
Location = TOXL 8164.13.00 06/09
04/27/2011 08:39

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	08:40
02 Std. Gas	0.081	08:40
03 Room Air	0.000	08:40
04 Std. Gas	0.080	08:41
05 Room Air	0.000	08:41
06 Std. Gas	0.081	08:42
07 Room Air	0.000	08:42

Lot No = 15410080A1
Cyl No = 19
Exp Date = 07/01/2012
County = 08 Oper No. = 777777

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080 AC

Form 106-I8000

80-003061
P. 129 of 129
KM