



# INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

OFFICE OF ATTORNEY GENERAL  
CRIME LABORATORY DIVISION  
Toxicology Section/Breath Alcohol Program  
SFN 59281 (5-2011)

Serial Number <i>80-004952</i>	Instrument Location <i>Hett - Adams County SO</i>
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level 1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. RFI Test (Use CMS Mode; Key Radio During First Room Air).
- 7. Repair and/or Maintenance Performed (if Any): \_\_\_\_\_
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>J Collins - 4100</i>	Date <i>04/27/2012</i>
Reviewed By (Crime Laboratory Use Only) <i>Deb Kashur</i>	Date <i>03 May 2012</i>

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004952  
Location = HETT                8164.13.00 06/09  
04/27/2012                      17:40

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^&\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^&\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?         1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            23411080A1  
Standard Cyl #?            59  
Standard Expiration?       10/01/2013  
Oper No?                    020533

Flow Cal. Date:            01/19/2010  
Slope                      694  
Intercept                  -638616

IR Calibration Date:        01/19/2010

	3um	9um
0th Coef(*100):	-32832	-34603
1st Coef(*100):	262566	131011
2nd Coef(*100):	1418	975
H2O adj(mg/l*10k):	686	550

\*\*\*\*\* Printer Test End \*\*\*\*\*

\_\_\_\_\_

Operator Signature  
TRAVIS COLLINS

Remarks:


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004952  
Location = HETT      8164.13.00 06/09  
04/27/2012      17:45

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:45
02 Std. Gas	0.078	17:46
03 Room Air	0.000	17:46
04 Std. Gas	0.078	17:47
05 Room Air	0.000	17:47
06 Std. Gas	0.077	17:48
07 Room Air	0.000	17:48

Lot No = 23411080A1  
Cyl No = 59  
Exp Date = 10/01/2013  
County = 01      Oper No. = 020533

  
\_\_\_\_\_  
Operator Signature  
TRAVIS COLLINS

Remarks:

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004952  
Location = HETT      8164.13.00 06/09  
04/27/2012      17:52

Test	AC	Time
01 Diagnostic	OK	17:54
02 Room Air	RFI*	17:54
03 Room Air	0.000	17:54

\*Invalid Test  
Inhibited - RFI

Sub Name = TEST, TEST T  
Sub DOB = 04/27/2010  
Sub Sex = Male      Weight = 120  
Test = DUI      Cit = 0000000  
Dr. Lic. = ND/000000  
Lot No = 23411080A1  
Cyl No = 59  
Expiration Date = 10/01/2013  
County = 01      Oper No. = 020533

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Collins - 4100

Operator Signature  
TRAVIS COLLINS

Remarks: