

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-004200

Location: MOTT (HETTINGER Co. S.O.)

- Check When Done
- A. Is the breath tube heated? ✓
  - B. Is the diagnostic check complete? ✓
  - C. Is the time, date, and year correct (re-set if necessary)? ✓
  - D. Print test. (Attach test record.) ✓
  - E. Low AC. Use  $\leq 0.03$  AC in ACA mode. (Attach test record.)  
Sim SN DR5143 Lot No. 201501B Rep AC 0.020 ✓
  - F. High AC. Use  $\geq 0.25$  AC in ACA mode. (Attach test record.)  
Sim SN DR3455 Lot No. 14270 Rep AC 0.300 ✓
  - G. Interferent check. Use 0.05% acetone in ethanol solution in ABA mode.  
(Attach test record.)  
Sim SN DR5132 Lot No. ICS#2 Solution 0.10AC + 0.05% ACE ✓
  - H. RFI check. Run CMS mode. Key radio during CMS test.  
(Attach test record.) ✓
  - I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.  
(Attach test record.) 0.080 AC

Lot No. 33913080A4 Cyl No. 14 Exp. Date 1/1/2016  
Test 1 0.081 Test 4 0.081 Test 7 0.080  
Test 2 0.082 Test 5 0.082 Test 8 0.080  
Test 3 0.081 Test 6 0.081 Test 9 0.080  
Average 0.081

J. Tank pressure, Level 3 Func D Sub G P: 1011 psi Reg 1000+ psi

K. Remarks/Maintenance record: ANNUAL INSPECTION

4.22.15  
Date

28 Apr 2015 ✓  
Date

Charles E. Don  
Field Inspector's Signature

Deb Kashner  
Reviewed by

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004200  
Location = MOTT                    8164.13.00 06/09  
04/22/2015                            09:44

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?         1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            33913080A4  
Standard Cyl #?            14  
Standard Expiration?       01/01/2016  
Oper No?                    666666

Flow Cal. Date:            12/10/2012  
Slope                        637  
Intercept                   -382487

IR Calibration Date:        12/10/2012  
                                  3um                            9um

-----  
0th Coef(\*100):            -14942                        -22738  
1st Coef(\*100):            256196                        136921  
2nd Coef(\*100):            3535                            1396  
H2O adj(mg/l\*10k):        492                            515

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:

PRINT TEST

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004200  
Location = MOTT      8164.13.00 06/09  
04/22/2015      09:59

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:00
02 Std. Sol.	0.017	10:00
03 Room Air	0.000	10:01
04 Std. Sol.	0.018	10:02
05 Room Air	0.000	10:02
06 Std. Sol.	0.019	10:03
07 Room Air	0.000	10:03

08 Sim Temp = 34.0°C

Simul Ser No = DR5143  
Std Sol No = 201501B  
County = 21      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *Low AC*  
*0.020 AC*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004200  
Location = MOTT      8164.13.00 06/09  
04/22/2015      10:05

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:05
02 Std. Sol.	0.295	10:06
03 Room Air	0.000	10:07
04 Std. Sol.	0.300	10:07
05 Room Air	0.000	10:08
06 Std. Sol.	0.301	10:09
07 Room Air	0.000	10:09


08 Sim Temp = 34.0°C

Simul Ser No = DR3455

Std Sol No = 14270

County = 21

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

HIGH AC  
0.300 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004200  
Location = MOTT                    8164.13.00 06/09  
04/22/2015                            09:55

Test	AC	Time
01 Room Air	0.000	09:56
02 *Subject Test	INT*	09:56
03 Room Air	0.000	09:57

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = INTERFERENT  
Dr. Lic. = ND/DIS821456  
Lot No = 33913080A4  
Cyl No = 14  
Expiration Date = 01/01/2016  
County = 21                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *Interferent Check*  
*0.10Ac + 0.05% ACETONE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004200  
Location = MOTT      8164.13.00 06/09  
04/22/2015      09:45

Test	AC	Time
01 Diagnostic	OK	09:46
02 Room Air	RFI*	09:46
03 Room Air	0.000	09:47

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = RFI CHECK  
Dr. Lic. = ND/DIS821456  
Lot No = 33913080A4  
Cyl No = 14  
Expiration Date = 01/01/2016  
County = 21      Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks: *ACC RFI CHECK*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004200  
Location = MOTT      8164.13.00 06/09  
04/22/2015      09:47

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:47
02 Std. Gas	0.081	09:48
03 Room Air	0.000	09:48
04 Std. Gas	0.082	09:49
05 Room Air	0.000	09:49
06 Std. Gas	0.081	09:50
07 Room Air	0.000	09:50

Lot No = 33913080A4  
Cyl No = 14  
Exp Date = 01/01/2016  
County = 21      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: ~~ACA TEST CEE~~  
0.080 AC      CALIBRATION  
Check


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004200  
Location = MOTT      8164.13.00 06/09  
04/22/2015      09:50

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:51
02 Std. Gas	0.081	09:51
03 Room Air	0.000	09:52
04 Std. Gas	0.082	09:52
05 Room Air	0.000	09:53
06 Std. Gas	0.081	09:53
07 Room Air	0.000	09:54

Lot No = 33913080A4  
Cyl No = 14  
Exp Date = 01/01/2016  
County = 21      Oper No. = 666666

  
Operator Signature  
CHARLES EDER

Remarks: ~~ACA TEST~~ CEE  
0.080 AC CALIBRATION  
Check



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004200  
Location = MOTT      8164.13.00 06/09  
04/22/2015      10:11

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:12
02 Std. Gas	0.080	10:12
03 Room Air	0.000	10:13
04 Std. Gas	0.080	10:13
05 Room Air	0.000	10:13
06 Std. Gas	0.080	10:14
07 Room Air	0.000	10:14

Lot No = 33913080A4  
Cyl No = 14  
Exp Date = 01/01/2016  
County = 21      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *Calibration Check*  
*0.080 AC*