



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59281 (5-2011)

Serial Number	80-004956	Instrument Location	Grand Forks County Jail
Reason for Install/Repair			
<input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____			

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level 1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. RFI Test (Use CMS Mode; Key Radio During First Room Air).
- 7. Repair and/or Maintenance Performed (if Any): _____
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature	<i>Steven R Fischer</i>	Date	4/30/12 ^{DK} 4-20-12
Reviewed By (Crime Laboratory Use Only)	<i>Debi Kashner</i>	Date	03 Apr 12

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004956
Location = GFSO 8164.13.00 06/09
04/20/2012 21:49

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? No
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 23411080A1
Standard Cyl #? 05
Standard Expiration? 10/01/2013
Oper No? 049904

Flow Cal. Date: 01/25/2010
Slope 672
Intercept -627011

IR Calibration Date: 04/11/2012

	3um	9um
0th Coef(*100):	-8701	-10668
1st Coef(*100):	247454	124257
2nd Coef(*100):	4577	2072
H2O adj(mg/l*10k):	564	580

***** Printer Test End *****



Operator Signature
STEVEN FISCHER

Remarks:

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004956
Location = GFSO 8164.13.00 06/09
04/20/2012 21:45

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	21:45
02 Std. Gas	0.079	21:46
03 Room Air	0.000	21:46
04 Std. Gas	0.079	21:46
05 Room Air	0.000	21:47
06 Std. Gas	0.078	21:47
07 Room Air	0.000	21:48

Lot No = 23411080A1
Cyl No = 05
Exp Date = 10/01/2013
County = 18 Oper No. = 049904



Operator Signature
STEVEN FISCHER

Remarks:

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004956
Location = GFSO 8164.13.00 06/09
04/20/2012 21:50

DRY CAL CHECK

Test	AC	Time
01 Room Air	RFI*	21:50
02 Room Air	RFI*	21:50

*Invalid Test
Inhibited - RFI

Lot No = 23411080A1
Cyl No = 05
Exp Date = 10/01/2013
County = 18 Oper No. = 049904



Operator Signature
STEVEN FISCHER

Remarks:


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004956
Location = GFSO 8164.13.00 06/09
04/30/2012 14:00

Test	AC	Time
01 Diagnostic	OK	14:02
02 Room Air	0.000	14:02
03 *Subject Test	RFI*	14:03
04 Room Air	0.000	14:03

*Invalid Test
Inhibited - RFI

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
STEVEN FISCHER

Remarks: