



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59281 (5-2011)

Serial Number 80-004186	Instrument Location SHER
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level 1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. RFI Test (Use CMS Mode; Key Radio During First Room Air).
- 7. Repair and/or Maintenance Performed (if Any): N/A
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency. N/A
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>Daniel Hoffmann</i>	Date 05-01-2012
Reviewed By (Crime Laboratory Use Only) <i>Deb Kashner</i>	Date 03 May 2012

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004186
Location = SHER 8164.13.00 06/09
05/01/2012 10:00

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 23411080A1
Standard Cyl #? 1
Standard Expiration? 10/01/2013
Oper No? 048909

Flow Cal. Date: 06/03/2009
Slope 703
Intercept -668211

IR Calibration Date: 06/03/2009

	3um	9um
0th Coef(*100):	-12098	-31908
1st Coef(*100):	257966	136312
2nd Coef(*100):	2942	1009
H2O adj(mg/l*10k):	434	609

***** Printer Test End *****



Operator Signature
DELBERT HOFFMANN

Remarks:

PRINT TEST


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004186
Location = SHER 8164.13.00 06/09
05/01/2012 11:04

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:06
02 Std. Gas	0.080	11:07
03 Room Air	0.000	11:07
04 Std. Gas	0.079	11:08
05 Room Air	0.000	11:08
06 Std. Gas	0.079	11:09
07 Room Air	0.000	11:09

Lot No = 23411080A1
Cyl No = 1
Exp Date = 10/01/2013
County = 42 Oper No. = 048909


Operator Signature
DELBERT HOFFMANN

Remarks:

ACA

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004186
Location = SHER 8164.13.00 06/09
05/01/2012 10:16

Test	AC	Time
01 Diagnostic	OK	10:18
02 Room Air	RFI*	10:18
03 Room Air	0.000	10:19

*Invalid Test
 Inhibited - RFI

Sub Name = TEST, TEST NA
Sub DOB = 01/01/1990
Sub Sex = Male Weight = 150
Test = DUI Cit = 00056
Dr. Lic. = ND/555778
Lot No = 23411080A1
Cyl No = 1
Expiration Date = 10/01/2013
County = 42 Oper No. = 048909

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
DELBERT HOFFMANN

Remarks:

RFI

Form 106-I8000