

INTOXILYZER 8000 INSPECTION

Serial No.: 80-003060

Location: Mercer Co SD - Stanton

- |  | Check When Done |
|--|-----------------|
| A. Is the warm-up time less than 20 minutes?   | <u>✓</u>        |
| B. Is a three-pronged grounded outlet used?  | <u>✓</u>        |
| C. Is the breath tube heated?  | <u>✓</u>        |
| D. Is the diagnostic check complete?   | <u>✓</u>        |
| E. Is the time, date, and year correct (re-set if necessary)?  | <u>✓</u>        |
| F. Print test. (Attach test record.)   | <u>✓</u>        |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR5132</u> Lot No. <u>    </u> Rep AC <u>0.021</u>   | <u>✓</u>        |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR5115</u> Lot No. <u>    </u> Rep AC <u>0.298</u>                                | <u>✓</u>        |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.<br>(Attach test record.)<br>Sim SN <u>DR5186</u> Lot No. <u>    </u> Rep AC <u>0.100</u> | <u>✓</u>        |
| J. RFI check. Run CMS mode. Key radio on first room air.<br>(Attach test record.)  | <u>✓</u>        |
| K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.<br>(Attach test record.)  |                 |

Lot No. <u>659357</u>	Cyl No. <u>43</u>	Exp. Date <u>09/03/2011</u>
Test 1 <u>0.080</u>	Test 4 <u>0.081</u>	Test 7 <u>0.081</u>
Test 2 <u>0.080</u>	Test 5 <u>0.081</u>	Test 8 <u>0.081</u>
Test 3 <u>0.081</u>	Test 6 <u>0.081</u>	Test 9 <u>0.081</u>
Average <u>0.081</u>		

L. Tank pressure, Level 3 Func D Sub G P: 861 psi Reg 860 psi

M. Remarks/Maintenance record: Annual inspection

03 May 2010  
Date

Deb Kasheer  
Field Inspector's Signature

15 May 2010  
Date

Mary Pearson  
Reviewed by

80-003060

P. 784836  
JH

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-003060  
Location = STAT                    8164.13.00 06/09  
05/03/2010                            11:36

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^&\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?         1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            659357  
Standard Cyl #?            43  
Standard Expiration?       09/03/2011  
Oper No?                    777777

Flow Cal. Date:            04/08/2008  
Slope                       702  
Intercept                  -647610

IR Calibration Date:       06/04/2008  
                                    3um                    9um

-----  
0th Coef(\*100):            -5190                    -7737  
1st Coef(\*100):            270580                   138779  
2nd Coef(\*100):            2499                     956  
H2O adj(mg/l\*10k):        377                      188

\*\*\*\*\* Printer Test End \*\*\*\*\*

Deb Kasher

Operator Signature  
DEB KASHUR

Remarks:

Print test

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003060  
Location = STAT      8164.13.00 06/09  
05/03/2010      11:52

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:52
02 Std. Sol.	0.019	11:53
03 Room Air	0.000	11:54
04 Std. Sol.	0.019	11:55
05 Room Air	0.000	11:55
06 Std. Sol.	0.020	11:56
07 Room Air	0.000	11:56

08 Sim Temp = 34.0°C

Simul Ser No = 00000

Std Sol No = 0000

County = 29

Oper No. = 777777

*Deb Kashur*

Operator Signature  
DEB KASHUR

Remarks:

*LOW AC 0.021 AC*

Form 106-18000

80-003060  
P. 80-1886  
*yan*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003060  
Location = STAT      8164.13.00 06/09  
05/03/2010      11:57

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:58
02 Std. Sol.	0.297	11:58
03 Room Air	0.000	11:59
04 Std. Sol.	0.296	12:00
05 Room Air	0.000	12:00
06 Std. Sol.	0.295	12:01
07 Room Air	0.000	12:02

08 Sim Temp = 34.0°C

Simul Ser No = 00000

Std Sol No = 0000

County = 29

Oper No. = 777777

*Deb Kashur*

Operator Signature  
DEB KASHUR

Remarks:

*Linearity test 0.298 AC*

Form 106-I8000

80-003060  
P. 8/1/86  
JCM

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003060  
Location = STAT      8164.13.00 06/09  
05/03/2010      11:50

Test	AC	Time
01 Room Air	0.000	11:50
02 *Subject Test	INT*	11:51
03 Room Air	0.000	11:51

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = NA  
Dr. Lic. = ND/DIS921456  
Lot No = 659357  
Cyl No = 43  
Expiration Date = 09/03/2011  
County = 29      Oper No. = 777777

*Deb Kashur*

Operator Signature  
DEB KASHUR

Remarks:

*Interferent check      0.100 AC + acetone*

Form 106-I8000

80-003060  
P. 82 of 86  
JK

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003060  
Location = STAT      8164.13.00 06/09  
05/03/2010      11:48

Test	AC	Time
01 Diagnostic	OK	11:49
02 Room Air	RFI*	11:49
03 Room Air	0.000	11:50

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = NA  
Dr. Lic. = ND/DIS921456  
Lot No = 659357  
Cyl No = 43  
Expiration Date = 09/03/2011  
County = 29      Oper No. = 777777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Deb Kasher

Operator Signature  
DEB KASHUR

Remarks:

*RFI test*

Form 106-18000

80-003060  
P. 83 of 86  
AK

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003060  
Location = STAT      8164.13.00 06/09  
05/03/2010      11:37

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:38
02 Std. Gas	0.080	11:38
03 Room Air	0.000	11:39
04 Std. Gas	0.080	11:39
05 Room Air	0.000	11:40
06 Std. Gas	0.081	11:40
07 Room Air	0.000	11:40

Lot No = 659357

Cyl No = 43

Exp Date = 09/03/2011

County = 29

Oper No. = 777777

*Deb Kashur*

Operator Signature  
DEB KASHUR

Remarks:

*Calibration check 0.080 AC*

Form 106-I8000

80-003060  
P. 84 of 86  
JK

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003060  
Location = STAT      8164.13.00 06/09  
05/03/2010      11:41

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:41
02 Std. Gas	0.081	11:42
03 Room Air	0.000	11:42
04 Std. Gas	0.081	11:42
05 Room Air	0.000	11:43
06 Std. Gas	0.081	11:43
07 Room Air	0.000	11:44

Lot No = 659357

Cyl No = 43

Exp Date = 09/03/2011

County = 29

Oper No. = 777777

*Deb Kashur*

Operator Signature  
DEB KASHUR

Remarks:

*Calibration check 0.080 AC*

Form 106-I8000

80-003060  
P. 85 of 86  
JEM



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003060  
Location = STAT      8164.13.00 06/09  
05/03/2010      11:44

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:45
02 Std. Gas	0.081	11:45
03 Room Air	0.000	11:45
04 Std. Gas	0.081	11:46
05 Room Air	0.000	11:46
06 Std. Gas	0.081	11:47
07 Room Air	0.000	11:47

Lot No = 659357

Cyl No = 43

Exp Date = 09/03/2011

County = 29

Oper No. = 777777

*Deb Kashur*

Operator Signature  
DEB KASHUR

Remarks:

*Calibration check 0080 AC*

Form 106-I8000

80-003060  
p. 56 of 86  
gsm