

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-003066

Location: LAMO

- |  | Check When Done |
|--|-----------------|
| A. Is the breath tube heated?  | <u>✓</u>        |
| B. Is the diagnostic check complete?   | <u>✓</u>        |
| C. Is the time, date, and year correct (re-set if necessary)?  | <u>✓</u>        |
| D. Print test. (Attach test record.)   | <u>✓</u>        |
| E. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR5115</u> Lot No. <u>12111</u> Rep AC <u>0.020</u>  | <u>✓</u>        |
| F. High AC. Use ≥ 0.25 AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR5143</u> Lot No. <u>10-000</u> Rep AC <u>0.300</u><br><u>13080</u>   | <u>✓</u>        |
| G. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.<br>(Attach test record.)<br>Sim SN <u>DR3375</u> Lot No. <u>FCS #2</u> Rep AC <u>0.100AC + 0.05% ACETONE</u> | <u>✓</u>        |
| H. RFI check. Run CMS mode. Key radio on first room air.<br>(Attach test record.)  | <u>✓</u>        |
| I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.<br>(Attach test record.)  |                 |

Lot No. <u>16512080A2</u>	Cyl No. <u>007</u>	Exp. Date <u>08/01/2014</u>
Test 1 <u>0.078</u>	Test 4 <u>0.078</u>	Test 7 <u>0.078</u>
Test 2 <u>0.077</u>	Test 5 <u>0.078</u>	Test 8 <u>0.078</u>
Test 3 <u>0.077</u>	Test 6 <u>0.078</u>	Test 9 <u>0.078</u>
Average <u>0.078</u>		

J. Tank pressure, Level 3 Func D Sub G P: 987 psi Reg 1000 psi

K. Remarks/Maintenance record: ANNUAL INSPECTION

02 MAY 14  
Date

Charles E. Edgerly  
Field Inspector's Signature

05 May 14 ✓  
Date

Debi Kashner  
Reviewed by

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-003066  
Location = LAMO                8164.13.00 06/09  
05/02/2014                      13:00

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        No  
# of Print Copies?          1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            16512080A2  
Standard Cyl #?            7  
Standard Expiration?        08/01/2014  
Oper No?                    666666

Flow Cal. Date:            07/29/2010  
Slope                      806  
Intercept                  -895033

IR Calibration Date:        07/29/2010  
                              3um                    9um

-----  
0th Coef(\*100):            -14075                -5808  
1st Coef(\*100):            272113                138204  
2nd Coef(\*100):            1637                    986  
H2O adj(mg/l\*10k):        289                    148

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:

PRINT TEST

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003066  
Location = LAMO      8164.13.00 06/09  
05/02/2014      13:19

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:19
02 Std. Sol.	0.020	13:20
03 Room Air	0.000	13:21
04 Std. Sol.	0.020	13:21
05 Room Air	0.000	13:22
06 Std. Sol.	0.021	13:23
07 Room Air	0.000	13:23

08 Sim Temp = 34.0°C

Simul Ser No = DR5115

Std Sol No = 12111

County = 23

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Low AC*  
*0.020 AC*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003066  
Location = LAMO      8164.13.00 06/09  
05/02/2014      13:24

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:25
02 Std. Sol.	0.295	13:25
03 Room Air	0.000	13:26
04 Std. Sol.	0.296	13:27
05 Room Air	0.000	13:27
06 Std. Sol.	0.298	13:28
07 Room Air	0.000	13:29

08 Sim Temp = 34.0°C

Simul Ser No = DR5143

Std Sol No = 13080

County = 23

Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

HIGH AC  
0.300 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-003066  
Location = LAMO                8164.13.00 06/09  
05/02/2014                        13:15

Test	AC	Time
01 Room Air	0.000	13:15
02 *Subject Test	INT*	13:16
03 Room Air	0.000	13:17

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = INTERFERENT  
Dr. Lic. = ND/DIS821456  
Lot No = 16512080A2  
Cyl No = 7  
Expiration Date = 08/01/2014  
County = 23                            Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *INTERFERENT CHECK*  
*0.100 AC + 0.05% ACETONE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003066  
Location = LAMO      8164.13.00 06/09  
05/02/2014      13:12

Test	AC	Time
01 Diagnostic	OK	13:13
02 Room Air	RFI*	13:13
03 Room Air	0.000	13:14

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = RFI CHECK  
Dr. Lic. = ND/DIS821456  
Lot No = 16512080A2  
Cyl No = 7  
Expiration Date = 08/01/2014  
County = 23      Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

RFI CHECK

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003066  
Location = LAMO      8164.13.00 06/09  
05/02/2014      13:30

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:30
02 Std. Gas	0.078	13:31
03 Room Air	0.000	13:31
04 Std. Gas	0.077	13:32
05 Room Air	0.000	13:32
06 Std. Gas	0.077	13:33
07 Room Air	0.000	13:33

Lot No = 16512080A2  
Cyl No = 7  
Exp Date = 08/01/2014  
County = 23      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: 0.080 AC  
CALIBRATION CHECK

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-003066  
Location = LAMO                8164.13.00 06/09  
05/02/2014                        13:34

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:34
02 Std. Gas	0.078	13:35
03 Room Air	0.000	13:35
04 Std. Gas	0.078	13:36
05 Room Air	0.000	13:36
06 Std. Gas	0.078	13:36
07 Room Air	0.000	13:37

Lot No = 16512080A2  
Cyl No = 7  
Exp Date = 08/01/2014  
County = 23                      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: *CALIBRATION CHECK*  
*0.080 AC*



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003066  
Location = LAMO      8164.13.00 06/09  
05/02/2014      13:37

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:38
02 Std. Gas	0.078	13:38
03 Room Air	0.000	13:39
04 Std. Gas	0.078	13:39
05 Room Air	0.000	13:40
06 Std. Gas	0.078	13:40
07 Room Air	0.000	13:41

Lot No = 16512080A2  
Cyl No = 7  
Exp Date = 08/01/2014  
County = 23      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000