

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-005361

Location: JAME

- Check When Done
- A. Is the breath tube heated? ✓
 - B. Is the diagnostic check complete? ✓
 - C. Is the time, date, and year correct (re-set if necessary)? ✓
 - D. Print test. (Attach test record.) ✓
 - E. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)
Sim SN DR5115 Lot No. 12111 Rep AC 0.020 ✓
 - F. High AC. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN DR5143 Lot No. 13080 Rep AC 0.300 ✓
 - G. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.
(Attach test record.)
Sim SN DR3375 Lot No. ICS#2 Rep AC 0.100AC + 0.05% ACETONE ✓
 - H. RFI check. Run CMS mode. Key radio on first room air.
(Attach test record.) ✓
 - I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.)

Lot No. 16512080A²K09E Cyl No. 042 Exp. Date 08/01/2014
Test 1 0.079 Test 4 0.080 Test 7 0.080
Test 2 0.080 Test 5 0.079 Test 8 0.079
Test 3 0.080 Test 6 0.080 Test 9 0.079
Average 0.080

J. Tank pressure, Level 3 Func D Sub G P: 889 psi Reg 900 psi

K. Remarks/Maintenance record: ANNUAL INSPECTION

5/6/2014
Date

Charles E. Edler
Field Inspector's Signature

07 May 2014 ✓
Date

Deb Kashur
Reviewed by

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005361
Location = JAME 8164.13.00 06/09
05/06/2014 11:29

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 16512080A2
Standard Cyl #? 42
Standard Expiration? 08/01/2014
Oper No? 666666

Flow Cal. Date: 08/19/2011
Slope 645
Intercept -487053

IR Calibration Date: 01/13/2012

	3um	9um
0th Coef(*100):	-5818	-10599
1st Coef(*100):	251241	123805
2nd Coef(*100):	7259	2455
H2O adj(mg/l*10k):	715	798

***** Printer Test End *****



Operator Signature
CHARLES EDER

Remarks:

PRINT TEST

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005361
Location = JAME 8164.13.00 06/09
05/06/2014 11:42

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:43
02 Std. Sol.	0.017	11:43
03 Room Air	0.000	11:44
04 Std. Sol.	0.019	11:45
05 Room Air	0.000	11:45
06 Std. Sol.	0.019	11:46
07 Room Air	0.000	11:46

08 Sim Temp = 34.0°C

Simul Ser No = DR5115
Std Sol No = 12111
County = 47 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

LOW AC
0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005361
Location = JAME 8164.13.00 06/09
05/06/2014 11:53

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:54
02 Std. Sol.	0.301	11:54
03 Room Air	0.000	11:55
04 Std. Sol.	0.300	11:56
05 Room Air	0.000	11:56
06 Std. Sol.	0.301	11:57
07 Room Air	0.000	11:58

08 Sim Temp = 34.0°C

Simul Ser No = DR5143
Std Sol No = 13080
County = 47 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *HIGH AC*
0.300 AC

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005361
Location = JAME 8164.13.00 06/09
05/06/2014 11:37

Test	AC	Time
01 Room Air	0.000	11:38
02 *Subject Test	INT*	11:38
03 Room Air	0.000	11:39

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = INTERFERENT
Dr. Lic. = ND/DIS821456
Lot No = 16512080A2
Cyl No = 42
Expiration Date = 08/01/2014
County = 47 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

INTERFERENT CHECK
0.100AC + 0.05% ACETONE

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005361
Location = JAME 8164.13.00 06/09
05/06/2014 11:34

Test	AC	Time
01 Diagnostic	OK	11:35
02 Room Air	RFI*	11:35
03 Room Air	0.000	11:36

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = RFI CHECK
Dr. Lic. = ND/DIS821456
Lot No = 16512080A2
Cyl No = 42
Expiration Date = 08/01/2014
County = 47 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks:

RFI CHECK

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005361
Location = JAME 8164.13.00 06/09
05/06/2014 12:02

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:03
02 Std. Gas	0.080	12:03
03 Room Air	0.000	12:03
04 Std. Gas	0.079	12:04
05 Room Air	0.000	12:04
06 Std. Gas	0.080	12:05
07 Room Air	0.000	12:05

Lot No = 16512080A2
Cyl No = 42
Exp Date = 08/01/2014
County = 47 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005361
Location = JAME 8164.13.00 06/09
05/06/2014 12:05

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:06
02 Std. Gas	0.080	12:06
03 Room Air	0.000	12:07
04 Std. Gas	0.079	12:07
05 Room Air	0.000	12:08
06 Std. Gas	0.079	12:08
07 Room Air	0.000	12:09

Lot No = 16512080A2
Cyl No = 42
Exp Date = 08/01/2014
County = 47 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *CALIBRATION CHECK*
0.080 AC

Form 106-I8000