



# INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

OFFICE OF ATTORNEY GENERAL  
CRIME LABORATORY DIVISION  
Toxicology Section/Breath Alcohol Program  
SFN 59281 (5-2011)

|   |                                    |
|---|------------------------------------|
| Serial Number<br><b>80-004198</b>   | Instrument Location<br><b>GFPD</b> |
| Reason for Install/Repair<br><input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change<br><input type="checkbox"/> Other (Specify) _____ |                                    |

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level 1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. RFI Test (Use CMS Mode; Key Radio During First Room Air).
- 7. Repair and/or Maintenance Performed (if Any): \_\_\_\_\_
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

|   |                            |
|---|----------------------------|
| Field Inspector Signature<br><i>Jon Lugin</i>                 | Date<br><b>05-09-12</b>    |
| Reviewed By (Crime Laboratory Use Only)<br><i>Deb Kashner</i> | Date<br><b>14 May 2012</b> |

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004198  
Location = GFPD                 8164.13.00 06/09  
05/09/2012                        23:22

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPNOPQRSTUVWXYZ!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPNOPQRSTUVWXYZ!@#\$\$%^&\*()\_+?

Current Instrument Setup


Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?          1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.077  
Standard Lot #?            15410080A1  
Standard Cyl #?            37  
Standard Expiration?       07/01/2012  
Oper No?                    039947

Flow Cal. Date:            01/09/2012  
Slope                      637  
Intercept                  -467208

IR Calibration Date:        01/09/2012  
                                    3um                      9um

-----  
0th Coef(\*100):            -31777                      -32014  
1st Coef(\*100):            271021                      132277  
2nd Coef(\*100):            1454                        964  
H2O adj(mg/l\*10k):        558                         576

\*\*\*\*\* Printer Test End \*\*\*\*\*

  
Operator Signature  
JONATHAN R. LAMPI

Remarks:

Print test


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004198  
Location = GFPD      8164.13.00 06/09  
05/09/2012      23:23

DRY CAL CHECK

| Test        | AC    | Time  |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 23:24 |
| 02 Std. Gas | 0.080 | 23:24 |
| 03 Room Air | 0.000 | 23:25 |
| 04 Std. Gas | 0.080 | 23:25 |
| 05 Room Air | 0.000 | 23:26 |
| 06 Std. Gas | 0.080 | 23:26 |
| 07 Room Air | 0.000 | 23:26 |

Lot No = 15410080A1  
Cyl No = 37  
Exp Date = 07/01/2012  
County = 18      Oper No. = 039947

  
Operator Signature  
JONATHAN R. LAMPI

Remarks:  
ACA test

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004198  
Location = GFPD                8164.13.00 06/09  
05/09/2012                      23:31

| Test          | AC    | Time  |
|---------------|-------|-------|
| 01 Diagnostic | OK    | 23:32 |
| 02 Room Air   | RFI*  | 23:32 |
| 03 Room Air   | 0.000 | 23:33 |

\*Invalid Test  
Inhibited - RFI

Sub Name = CMSTEST, TEST T  
Sub DOB = 02/02/1987  
Sub Sex = Male                    Weight = 170  
Test = DUI                        Cit = 13244  
Dr. Lic. = ND/12312414  
Lot No = 15410080A1  
Cyl No = 37  
Expiration Date = 07/01/2012  
County = 18                        Oper No. = 039947

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

  
Operator Signature  
JONATHAN R. LAMPI

Remarks:

RFI TEST