

INTOXILYZER 8000 INSPECTION

Serial No.: 80-004200

Location: Trail CoSO

- | | Check When Done |
|--|-----------------|
| A. Is the warm-up time less than 20 minutes? | <u>✓</u> |
| B. Is a three-pronged grounded outlet used? | <u>✓</u> |
| C. Is the breath tube heated? | <u>✓</u> |
| D. Is the diagnostic check complete? | <u>✓</u> |
| E. Is the time, date, and year correct (re-set if necessary)? | <u>✓</u> |
| F. Print test. (Attach test record.) | <u>✓</u> |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)
Sim SN <u>DR 5132</u> Lot No. <u>—</u> Rep AC <u>0.021 AC</u> | <u>✓</u> |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN <u>DR 3375</u> Lot No. <u>—</u> Rep AC <u>0.292 AC</u> | <u>✓</u> |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.
(Attach test record.)
Sim SN <u>DR 5186</u> Lot No. <u>—</u> Rep AC <u>0.10 AC + 0.05% Acetone</u> | <u>✓</u> |
| J. RFI check. Run CMS mode. Key radio on first room air.
(Attach test record.) | <u>✓</u> |
| K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.) | |

Lot No. <u>659 358</u>	Cyl No. <u>45</u>	Exp. Date <u>9/03/2011</u>
Test 1 <u>0.080</u>	Test 4 <u>0.081</u>	Test 7 <u>0.081</u>
Test 2 <u>0.081</u>	Test 5 <u>0.080</u>	Test 8 <u>0.080</u>
Test 3 <u>0.081</u>	Test 6 <u>0.080</u>	Test 9 <u>0.080</u>
Average <u>0.080</u>		

L. Tank pressure, Level 3 Func D Sub G P: 428 psi Reg 475 psi

M. Remarks/Maintenance record: Annual Inspection

10 May 2010
Date

Mary Pearson
Field Inspector's Signature

15 May 2010
Date

Deb Kashner
Reviewed by

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = HILL 8164.13.00 06/09
05/10/2010 16:09

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.077
Standard Lot #? 659358
Standard Cyl #? 45
Standard Expiration? 09/03/2011
Oper No? 107501

Flow Cal. Date: 06/04/2009
Slope 715
Intercept -655593

IR Calibration Date: 12/21/2009
 3um 9um

0th Coef(*100):	-28821	-30199
1st Coef(*100):	265418	139385
2nd Coef(*100):	-207	427
H2O adj(mg/l*10k):	678	645

***** Printer Test End *****



Operator Signature
MARGY PEARSON

Remarks:

Print

Form 106-18000

80-004200
P. 77/84 192

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = HILL 8164.13.00 06/09
05/10/2010 16:33

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:34
02 Std. Sol.	0.018	16:35
03 Room Air	0.000	16:35
04 Std. Sol.	0.019	16:36
05 Room Air	0.000	16:37
06 Std. Sol.	0.018	16:37
07 Room Air	0.000	16:38

08 Sim Temp = 34.0°C

Simul Ser No = 000000

Std Sol No = 0000000

County = 49

Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Low 0.021 AC

Form 106-18000

80-004200
P. 7/8/84 KJK

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = HILL 8164.13.00 06/09
05/10/2010 16:27

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:28
02 Std. Sol.	0.278	16:29
03 Room Air	0.000	16:29
04 Std. Sol.	0.281	16:30
05 Room Air	0.000	16:30
06 Std. Sol.	0.284	16:31
07 Room Air	0.000	16:32


08 Sim Temp = 34.0°C

Simul Ser No = 000000

Std Sol No = 0000000

County = 49

Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Linearity 0.292 AC

Form 106-I8000

80-004200
P. 79/84 JAL


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = HILL 8164.13.00 06/09
05/10/2010 16:13

Test	AC	Time
01 Room Air	0.000	16:14
02 *Subject Test	INT*	16:14
03 Room Air	0.000	16:15

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Simul Ser No = 000000
Std Sol No = 0000000
County = 49 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Intror Alct Acctone

Form 106-I8000

80-004200
P. 80/84 FH

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = HILL 8164.13.00 06/09
05/10/2010 16:10

Test	AC	Time
01 Diagnostic	OK	16:12
02 Room Air	RFI*	16:12
03 Room Air	0.000	16:12

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 659358
Cyl No = 45
Expiration Date = 09/03/2011
County = 49 Oper No. = 107501

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
MARGY PEARSON

Remarks:

RFI

Form 106-I8000

80-004200
P. 8/84 VM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = HILL 8164.13.00 06/09
05/10/2010 16:40

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:40
02 Std. Gas	0.080	16:41
03 Room Air	0.000	16:41
04 Std. Gas	0.081	16:42
05 Room Air	0.000	16:42
06 Std. Gas	0.081	16:42
07 Room Air	0.000	16:43

Lot No = 659358
Cyl No = 45
Exp Date = 09/03/2011
County = 49 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Cal Check ACA

Form 106-I8000

80-004200
P. 83/84 KM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = HILL 8164.13.00 06/09
05/10/2010 16:43

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:44
02 Std. Gas	0.081	16:44
03 Room Air	0.000	16:45
04 Std. Gas	0.080	16:45
05 Room Air	0.000	16:46
06 Std. Gas	0.080	16:46
07 Room Air	0.000	16:46

Lot No = 659358
Cyl No = 45
Exp Date = 09/03/2011
County = 49 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Cal Check ACA

Form 106-I8000

*80-004200
P. 5/3/14 RAL*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004200
Location = HILL 8164.13.00 06/09
05/10/2010 16:47

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:47
02 Std. Gas	0.081	16:48
03 Room Air	0.000	16:48
04 Std. Gas	0.080	16:48
05 Room Air	0.000	16:49
06 Std. Gas	0.080	16:49
07 Room Air	0.000	16:50

Lot No = 659358

Cyl No = 45

Exp Date = 09/03/2011

County = 49

Oper No. = 107501


Operator Signature
MARGY PEARSON

Remarks:

Cal Check 0.080 AC

Form 106-I8000

80-004200
P. 8/1/84 *SM*