

INTOXILYZER 8000 INSPECTION

Serial No.: 80004944

Location: Carrington PD

- | | Check When Done |
|---|-----------------|
| A. Is the warm-up time less than 20 minutes? | <u>✓</u> |
| B. Is a three-pronged grounded outlet used? | <u>✓</u> |
| C. Is the breath tube heated? | <u>✓</u> |
| D. Is the diagnostic check complete? | <u>✓</u> |
| E. Is the time, date, and year correct (re-set if necessary)? | <u>✓</u> |
| F. Print test. (Attach test record.) | <u>✓</u> |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)
Sim SN <u>DR5132</u> Lot No. <u>—</u> Rep AC <u>0.021</u> | <u>✓</u> |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN <u>DR3375</u> Lot No. <u>—</u> Rep AC <u>0.292</u> | <u>✓</u> |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.
(Attach test record.)
Sim SN <u>DR5186</u> Lot No. <u>—</u> Rep AC <u>0.10 AC + Acetone</u> | <u>✓</u> |
| J. RFI check. Run CMS mode. Key radio on first room air.
(Attach test record.) | <u>✓</u> |
| K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.) | |

Lot No. <u>709352</u>	Cyl No. <u>12</u>	Exp. Date <u>12/23/11</u>
Test 1 <u>0.080</u>	Test 4 <u>0.078</u>	Test 7 <u>0.079</u>
Test 2 <u>0.079</u>	Test 5 <u>0.079</u>	Test 8 <u>0.080</u>
Test 3 <u>0.079</u>	Test 6 <u>0.079</u>	Test 9 <u>0.078</u>
Average <u>0.079 AC</u>		

L. Tank pressure, Level 3 Func D Sub G P: 750 psi Reg 740 psi

M. Remarks/Maintenance record: Annual Inspection

<u>10 May 2010</u> Date	<u>Mary Pearson</u> Field Inspector's Signature
<u>15 May 2010</u> Date	<u>Deb Kashner</u> Reviewed by

80-004944
P. 40/48 KM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004944
Location = CARR 8164.13.00 06/09
05/10/2010 12:16

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? No
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 709352
Standard Cyl #? 12
Standard Expiration? 12/23/2011
Oper No? 107501

Flow Cal. Date: 01/18/2010
Slope 693
Intercept -679634

IR Calibration Date: 01/18/2010
 3um 9um

0th Coef(*100): -29206 -25219
1st Coef(*100): 257889 131140
2nd Coef(*100): 2394 1225
H2O adj(mg/l*10k): 734 452

***** Printer Test End *****



Operator Signature
MARGY PEARSON

Remarks:

Print Test

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004944
Location = CARR 8164.13.00 06/09
05/10/2010 12:29

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:30
02 Std. Sol.	0.016	12:32
03 Room Air	0.000	12:32
04 Std. Sol.	0.017	12:33
05 Room Air	0.000	12:34
06 Std. Sol.	0.017	12:34
07 Room Air	0.000	12:35

08 Sim Temp = 34.0°C

Simul Ser No = 00000000

Std Sol No = 0000000

County = 16

Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Low AC 0.021 AC

Form 106-18000

80-004944
8.42/48 FK

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004944
Location = CARR 8164.13.00 06/09
05/10/2010 12:22

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:24
02 Std. Sol.	0.286	12:25
03 Room Air	0.000	12:25
04 Std. Sol.	0.290	12:26
05 Room Air	0.000	12:26
06 Std. Sol.	0.293	12:27
07 Room Air	0.000	12:28


08 Sim Temp = 34.0°C

Simul Ser No = 00000000

Std Sol No = 0000000

County = 16

Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Linearity 0.292 AC

Form 106-I8000

80-004944
p. 43/48 AH

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004944
Location = CARR 8164.13.00 06/09
05/10/2010 12:19

Test	AC	Time
01 Room Air	0.000	12:20
02 *Subject Test	INT*	12:21
03 Room Air	0.000	12:21

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Simul Ser No = 00000000
Std Sol No = 00000000
County = 16 Oper No. = 107501


Operator Signature
MARGY PEARSON

Remarks:

Interferer ALC + Acetone

Form 106-18000

80-604944
P. 44/48 AR

Intoxilyzer Test Record and Checklist
NDOAG-Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004944
Location = CARR 8164.13.00 06/09
05/10/2010 12:40

Test	AC	Time
01 Diagnostic	OK	12:41
02 Room Air	RFI*	12:41
03 Room Air	0.000	12:42

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Simul Ser No = 00000000
Std Sol No = 00000000
County = 16 Oper No. = 107501

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
MARGY PEARSON

Remarks:

RFI

Form 106-I8000

80-004944
p. 45/48 KM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004944
Location = CARR 8164.13.00 06/09
05/10/2010 12:45

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:45
02 Std. Gas	0.080	12:46
03 Room Air	0.000	12:46
04 Std. Gas	0.079	12:47
05 Room Air	0.000	12:47
06 Std. Gas	0.079	12:47
07 Room Air	0.000	12:48

Lot No = 709352
Cyl No = 12
Exp Date = 12/23/2011
County = 16 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Cal Check ACA

Form 106-I8000

*80-004944
P. 4/48 FM*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004944
Location = CARR 8164.13.00 06/09
05/10/2010 12:49

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:49
02 Std. Gas	0.078	12:50
03 Room Air	0.000	12:50
04 Std. Gas	0.079	12:50
05 Room Air	0.000	12:51
06 Std. Gas	0.079	12:51
07 Room Air	0.000	12:52

Lot No = 709352
Cyl No = 12
Exp Date = 12/23/2011
County = 16 Oper No. = 107501

Margy Pearson

Operator Signature
MARGY PEARSON

Remarks:

Cal Check ACA

Form 106-I8000

80-004944
P. 47/48 AM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004944
Location = CARR 8164.13.00 06/09
05/10/2010 12:52

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:53
02 Std. Gas	0.079	12:53
03 Room Air	0.000	12:54
04 Std. Gas	0.080	12:54
05 Room Air	0.000	12:54
06 Std. Gas	0.078	12:55
07 Room Air	0.000	12:55

Lot No = 709352
Cyl No = 12
Exp Date = 12/23/2011
County = 16 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Cal Check - ACA

Form 106-I8000

80-004944
P. 48/48 AM