

INTOXILYZER 8000 INSPECTION

Serial No.: 80-003057

Location: GF Corr Cent

- |  | Check When Done |
|--|-----------------|
| A. Is the warm-up time less than 20 minutes?   | <u>✓</u>        |
| B. Is a three-pronged grounded outlet used?  | <u>✓</u>        |
| C. Is the breath tube heated?  | <u>✓</u>        |
| D. Is the diagnostic check complete?   | <u>✓</u>        |
| E. Is the time, date, and year correct (re-set if necessary)?  | <u>✓</u>        |
| F. Print test. (Attach test record.)   | <u>✓</u>        |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR5132</u> Lot No. <u>—</u> Rep AC <u>0.021 AC</u>   | <u>✓</u>        |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR3375</u> Lot No. <u>—</u> Rep AC <u>0.292 AC</u>  | <u>✓</u>        |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.<br>(Attach test record.)<br>Sim SN <u>DR5186</u> Lot No. <u>—</u> Rep AC <u>0.10 AC + 0.05 Acetone</u> | <u>✓</u>        |
| J. RFI check. Run CMS mode. Key radio on first room air.<br>(Attach test record.)  | <u>✓</u>        |
| K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.<br>(Attach test record.)  |                 |

Lot No. <u>659357</u>	Cyl No. <u>26</u>	Exp. Date <u>09/03/2011</u>
Test 1 <u>0.081</u>	Test 4 <u>0.082</u>	Test 7 <u>0.082</u>
Test 2 <u><del>0.082</del> 0.081</u>	Test 5 <u>0.082</u>	Test 8 <u>0.082</u>
Test 3 <u>0.082</u>	Test 6 <u>0.082</u>	Test 9 <u>0.082</u>
Average <u>0.082</u>		

L. Tank pressure, Level 3 Func D Sub G P: 300 psi Reg 300 psi

M. Remarks/Maintenance record: Annual Inspection

11 May 2010  
Date  
15  
DK 17 May 2010  
Date

Mary Pearson  
Field Inspector's Signature  
Deb Kashner  
Reviewed by

80-003057  
P. 92 of 100  
JK

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003057  
Location = GFSO      8164.13.00 06/09  
05/11/2010      08:13

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-=  
ABCDEFGHIJKLMN0PQRSTUVWXYZ!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-=  
ABCDEFGHIJKLMN0PQRSTUVWXYZ!@#\$\$%^&\*()\_+?

Current Instrument Setup


Data Entry Mode:      Enabled  
Start Test Sequence:      DABACABA  
Display Prelim Rslt?      Yes  
Display Third Digit?      Yes  
Inhib Printer(Y/N)?      No  
Display Volume?      No  
Disable On Memfull?      Yes  
# of Print Copies?      1  
Select Std (D/W/I)?      Dry  
Standard Value?      0.080  
Standard Lot #?      659357  
Standard Cyl #?      26  
Standard Expiration?      09/03/2011  
Oper No?      107501

Flow Cal. Date:      04/07/2008  
Slope      718  
Intercept      -644065

IR Calibration Date:      06/05/2008

	3um	9um
0th Coef(*100):	-22884	-16302
1st Coef(*100):	265055	140885
2nd Coef(*100):	151	572
H2O adj(mg/l*10k):	508	293

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
MARGY PEARSON

Remarks:

Print Test

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003057  
Location = GFSSO      8164.13.00 06/09  
05/11/2010      08:30

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	08:31
02 Std. Sol.	0.018	08:31
03 Room Air	0.000	08:32
04 Std. Sol.	0.019	08:33
05 Room Air	0.000	08:33
06 Std. Sol.	0.019	08:34
07 Room Air	0.000	08:34

08 Sim Temp = 34.0°C

Simul Ser No = 0000000

Std Sol No = 0000

County = 18

Oper No. = 107501

*Margy Pearson*

Operator Signature  
MARGY PEARSON

Remarks:

*LOW AC 0.021 AC*

Form 106-18000

80-003057  
P. 944/100  
PM

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003057  
Location = GFSO      8164.13.00 06/09  
05/11/2010      08:24

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	08:24
02 Std. Sol.	0.280	08:25
03 Room Air	0.000	08:26
04 Std. Sol.	0.287	08:26
05 Room Air	0.000	08:27
06 Std. Sol.	0.291	08:28
07 Room Air	0.000	08:28

08 Sim Temp = 34.0°C

Simul Ser No = 0000000

Std Sol No = 0000

County = 18

Oper No. = 107501

*Margy Pearson*

Operator Signature  
MARGY PEARSON

Remarks:

*Linearity 0.292 AC*

Form 106-I8000

80-003057  
P. 95-1100  
*ph*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003057  
Location = GFSO      8164.13.00 06/09  
05/11/2010      08:19

Test	AC	Time
01 Room Air	0.000	08:20
02 *Subject Test	INT*	08:20
03 Room Air	0.000	08:21

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = NA  
Dr. Lic. = ND/DIS821456  
Lot No = 659357  
Cyl No = 26  
Expiration Date = 09/03/2011  
County = 18      Oper No. = 107501



Operator/Signature  
MARGY PEARSON

Remarks:

*Interferent Alc + Acetone*

Form 106-18000

80-003057  
P. 96 of 100  
KM

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003057  
Location = GFSC      8164.13.00 06/09  
05/11/2010      08:22

Test	AC	Time
01 Diagnostic	OK	08:23
02 Room Air	RFI*	08:23
03 Room Air	0.000	08:24

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = NA  
Dr. Lic. = ND/DIS821456  
Simul Ser No = 0000000  
Std Sol No = 0000  
County = 18      Oper No. = 107501

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
MARGY PEARSON

Remarks:

RFI

Form 106-I8000

80-003057  
P. 97 of 100  
VSM

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003057  
Location = GFSO      8164.13.00 06/09  
05/11/2010      08:37

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	08:38
02 Std. Gas	0.081	08:38
03 Room Air	0.000	08:39
04 Std. Gas	0.081	08:39
05 Room Air	0.000	08:39
06 Std. Gas	0.082	08:40
07 Room Air	0.000	08:40

Lot No = 659357  
Cyl No = 26  
Exp Date = 09/03/2011  
County = 18      Oper No. = 107501

*Margy Pearson*

Operator Signature  
MARGY PEARSON

Remarks:

*Cal Check 0.080 AC*

Form 106-I8000

80-003057  
P. 98 of 100  
JKM


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003057  
Location = GFSSO      8164.13.00 06/09  
05/11/2010      08:41

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	08:41
02 Std. Gas	0.082	08:42
03 Room Air	0.000	08:42
04 Std. Gas	0.082	08:43
05 Room Air	0.000	08:43
06 Std. Gas	0.082	08:43
07 Room Air	0.000	08:44

Lot No = 659357  
Cyl No = 26  
Exp Date = 09/03/2011  
County = 18      Oper No. = 107501



Operator Signature  
MARGY PEARSON

Remarks:

*Cal check 0.080 AC*

Form 106-I8000

*80-003057  
P. 99/100  
yem*



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003057  
Location = GFSO      8164.13.00 06/09  
05/11/2010      08:44

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	08:45
02 Std. Gas	0.082	08:45
03 Room Air	0.000	08:46
04 Std. Gas	0.082	08:46
05 Room Air	0.000	08:47
06 Std. Gas	0.082	08:47
07 Room Air	0.000	08:48

Lot No = 659357

Cyl No = 26

Exp Date = 09/03/2011

County = 18

Oper No. = 107501

  
Operator Signature  
MARGY PEARSON

Remarks:

*Cal Check 0.080 AC*

Form 106-18000

80-003057  
P. 100 of 100  
JRW