

INTOXILYZER 8000 INSPECTION

Serial No.: 80-004197

Location: BIA - BELCOURT

- | | Check When Done |
|---|-----------------|
| A. Is the warm-up time less than 20 minutes? | <u>✓</u> |
| B. Is a three-pronged grounded outlet used? | <u>✓</u> |
| C. Is the breath tube heated? | <u>✓</u> |
| D. Is the diagnostic check complete? | <u>✓</u> |
| E. Is the time, date, and year correct (re-set if necessary)? | <u>✓</u> |
| F. Print test. (Attach test record.) | <u>✓</u> |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.) Sim SN <u>DR5132</u> Lot No. <u>—</u> Rep AC <u>0.0211</u> | <u>✓</u> |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.) Sim SN <u>DR3375</u> Lot No. <u>—</u> Rep AC <u>0.292</u> | <u>✓</u> |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode. (Attach test record.) Sim SN <u>DR5186</u> Lot No. <u>—</u> Rep AC <u>0.10 AC + 0.05% Acetone</u> | <u>✓</u> |
| J. RFI check. Run CMS mode. Key radio on first room air. (Attach test record.) | <u>✓</u> |
| K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode. (Attach test record.) | |

| | | |
|-------------------------|---------------------|----------------------------|
| Lot No. <u>659357</u> | Cyl No. <u>38</u> | Exp. Date <u>9/03/2011</u> |
| Test 1 <u>0.082</u> | Test 4 <u>0.082</u> | Test 7 <u>0.081</u> |
| Test 2 <u>0.082</u> | Test 5 <u>0.081</u> | Test 8 <u>0.081</u> |
| Test 3 <u>0.081</u> | Test 6 <u>0.082</u> | Test 9 <u>0.080</u> |
| Average <u>0.081 AC</u> | | |

L. Tank pressure, Level 3 Func D Sub G P: 780 psi Reg 800 psi

M. Remarks/Maintenance record: Annual Inspection; Forms Load, Download & Clear
12 May 2010 Margy Pearson
Date Field Inspector's Signature
DK-17 May 2010 Deb Kasher
Date Reviewed by

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004197
Location = BELC 8164.13.00 06/09
05/12/2010 09:24

***** Printer Test *****

abcdefghijklmnopqrstuvwxy1234567890-=
|
ABCDEFGHIJKLMN0PQRSTUVWXYZ!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxy1234567890-=
|
ABCDEFGHIJKLMN0PQRSTUVWXYZ!@#\$\$%^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 659357
Standard Cyl #? 38
Standard Expiration? 09/03/2011
Oper No? 107501

Flow Cal. Date: 06/04/2009
Slope 697
Intercept -597078

IR Calibration Date: 06/03/2009
 3um 9um

0th Coef(*100): -18382 -32622
1st Coef(*100): 278299 134602
2nd Coef(*100): 929 921
H2O adj(mg/l*10k): 352 627

***** Printer Test End *****



Operator Signature
MARGY PEARSON

Remarks:

Print

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004197
Location = BELC 8164.13.00 06/09
05/12/2010 09:34

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 09:35 |
| 02 Std. Sol. | 0.018 | 09:36 |
| 03 Room Air | 0.000 | 09:36 |
| 04 Std. Sol. | 0.018 | 09:37 |
| 05 Room Air | 0.000 | 09:38 |
| 06 Std. Sol. | 0.019 | 09:38 |
| 07 Room Air | 0.000 | 09:39 |

08 Sim Temp = 34.0°C

Simul Ser No = 0000000
Std Sol No = 0000000
County = 63 Oper No. = 107501

Margy Pearson

Operator Signature
MARGY PEARSON

Remarks:

P. Low AC 0.021 AC

Form 106-I8000

80-004197
P. 7/16 192

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004197
Location = BELC 8164.13.00 06/09
05/12/2010 09:40

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 09:41 |
| 02 Std. Sol. | 0.294 | 09:41 |
| 03 Room Air | 0.000 | 09:42 |
| 04 Std. Sol. | 0.293 | 09:43 |
| 05 Room Air | 0.000 | 09:43 |
| 06 Std. Sol. | 0.293 | 09:44 |
| 07 Room Air | 0.000 | 09:45 |


08 Sim Temp = 34.0°C

Simul Ser No = 0000000

Std Sol No = 0000000

County = 63

Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Linearity 0.292 AC

Form 106-I8000

80-004197
P. 7/1/10 FM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004197
Location = BELC 8164.13.00 06/09
05/12/2010 09:31

| Test | AC | Time |
|------------------|-------|-------|
| 01 Room Air | 0.000 | 09:32 |
| 02 *Subject Test | INT* | 09:33 |
| 03 Room Air | 0.000 | 09:33 |

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Simul Ser No = 0000000
Std Sol No = 0000000
County = 63 Oper No. = 107501

Operator Signature
MARGY PEARSON

Remarks:

Interferent AC + Acetone

Form 106-I8000

*80-004197
P. 7/15 YK*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004197
Location = BELC 8164.13.00 06/09
05/12/2010 09:28

| Test | AC | Time |
|---------------|-------|-------|
| 01 Diagnostic | OK | 09:29 |
| 02 Room Air | RFI* | 09:29 |
| 03 Room Air | 0.000 | 09:30 |

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Simul Ser No = 0000000
Std Sol No = 0000000
County = 63 Oper No. = 107501

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
MARGY PEARSON

Remarks:

RFI

Form 106-I8000

80-004197
P. 73/76 *PK*


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004197
Location = BELC 8164.13.00 06/09
05/12/2010 09:46

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 09:46 |
| 02 Std. Gas | 0.082 | 09:47 |
| 03 Room Air | 0.000 | 09:47 |
| 04 Std. Gas | 0.082 | 09:48 |
| 05 Room Air | 0.000 | 09:48 |
| 06 Std. Gas | 0.081 | 09:49 |
| 07 Room Air | 0.000 | 09:49 |

Lot No = 659357
Cyl No = 38
Exp Date = 09/03/2011
County = 63 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Cal Check 0.080 AC

Form 106-I8000

80-004197
P. 74/75 KM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004197
Location = BELC 8164.13.00 06/09
05/12/2010 09:52

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 09:52 |
| 02 Std. Gas | 0.082 | 09:53 |
| 03 Room Air | 0.000 | 09:53 |
| 04 Std. Gas | 0.081 | 09:53 |
| 05 Room Air | 0.000 | 09:54 |
| 06 Std. Gas | 0.082 | 09:54 |
| 07 Room Air | 0.000 | 09:55 |

Lot No = 659357
Cyl No = 38
Exp Date = 09/03/2011
County = 63 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Cal Check 0.080 AC

Form 106-I8000

80-004197
P. 75/76 YHM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004197
Location = BELC 8164.13.00 06/09
05/12/2010 09:55

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 09:56 |
| 02 Std. Gas | 0.081 | 09:56 |
| 03 Room Air | 0.000 | 09:57 |
| 04 Std. Gas | 0.081 | 09:57 |
| 05 Room Air | 0.000 | 09:58 |
| 06 Std. Gas | 0.080 | 09:58 |
| 07 Room Air | 0.000 | 09:59 |

Lot No = 659357
Cyl No = 38
Exp Date = 09/03/2011
County = 63 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Cal Check 0.080 AC

Form 106-18000

*80-004197
P. 76/76 JH*