

**INTOXILYZER 8000 INSPECTION**

Serial No.: 80-004199

Location: Rolette Co SO

- |  | Check When Done |
|--|-----------------|
| A. Is the warm-up time less than 20 minutes?   | <u>✓</u>        |
| B. Is a three-pronged grounded outlet used?  | <u>✓</u>        |
| C. Is the breath tube heated?  | <u>✓</u>        |
| D. Is the diagnostic check complete?   | <u>✓</u>        |
| E. Is the time, date, and year correct (re-set if necessary)?  | <u>✓</u>        |
| F. Print test. (Attach test record.)   | <u>✓</u>        |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR5132</u> Lot No. <u>—</u> Rep AC <u>0.021 AC</u>   | <u>✓</u>        |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR3315</u> Lot No. <u>—</u> Rep AC <u>0.292 AC</u>  | <u>✓</u>        |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.<br>(Attach test record.)<br>Sim SN <u>DR5186</u> Lot No. <u>—</u> Rep AC <u>0.10 AC + 0.050% Acetone</u> | <u>✓</u>        |
| J. RFI check. Run CMS mode. Key radio on first room air.<br>(Attach test record.)  | <u>✓</u>        |
| K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.<br>(Attach test record.)  |                 |

|                       |                     |                             |
|-----------------------|---------------------|-----------------------------|
| Lot No. <u>659358</u> | Cyl No. <u>11</u>   | Exp. Date <u>09/03/2011</u> |
| Test 1 <u>0.080</u>   | Test 4 <u>0.081</u> | Test 7 <u>0.080</u>         |
| Test 2 <u>0.081</u>   | Test 5 <u>0.080</u> | Test 8 <u>0.081</u>         |
| Test 3 <u>0.080</u>   | Test 6 <u>0.080</u> | Test 9 <u>0.081</u>         |
| Average _____         |                     |                             |

L. Tank pressure, Level 3 Func D Sub G P: 812 psi Reg 800 psi

M. Remarks/Maintenance record: Annual - Inspection

12 May 2010  
Date

Mary Plason  
Field Inspector's Signature

15  
17 May 2010  
Date

Deb Kashner  
Reviewed by

80-004199  
P. 78 of 86  
FM

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004199  
Location = ROLL                8164.13.00 06/09  
05/12/2010                      08:10

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMNPNOPQRSTUVWXYZ!@#%&^\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-|=|  
ABCDEFGHIJKLMNPNOPQRSTUVWXYZ!@#%&^\*()\_+?

Current Instrument Setup


Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?          1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            659358  
Standard Cyl #?            11  
Standard Expiration?       09/03/2011  
Oper No?                    107501

Flow Cal. Date:            06/04/2009  
Slope                        701  
Intercept                  -580043

IR Calibration Date:        06/03/2009  
                                    3um                      9um

-----  
0th Coef(\*100):            -44451                    -27073  
1st Coef(\*100):            294251                    132474  
2nd Coef(\*100):            -1375                     1189  
H2O adj(mg/l\*10k):        315                        545

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
MARGY PEARSON

Remarks:

Print

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004199  
Location = ROLL      8164.13.00 06/09  
05/12/2010      08:20

WET CAL CHECK

| Test         | AC    | Time  |
|--------------|-------|-------|
| 01 Room Air  | 0.000 | 08:20 |
| 02 Std. Sol. | 0.017 | 08:22 |
| 03 Room Air  | 0.000 | 08:23 |
| 04 Std. Sol. | 0.018 | 08:23 |
| 05 Room Air  | 0.000 | 08:24 |
| 06 Std. Sol. | 0.017 | 08:25 |
| 07 Room Air  | 0.000 | 08:25 |

08 Sim Temp = 34.0°C

Simul Ser No = 000000

Std Sol No = 00000

County = 40

Oper No. = 107501

*Margy Pearson*

Operator Signature  
MARGY PEARSON

Remarks:

*Low 0.021 AC*

Form 106-18000

80-004199  
P. 80 of 86  
KHL

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004199  
Location = ROLL      8164.13.00 06/09  
05/12/2010      08:27

WET CAL CHECK

| Test         | AC    | Time  |
|--------------|-------|-------|
| 01 Room Air  | 0.000 | 08:28 |
| 02 Std. Sol. | 0.301 | 08:29 |
| 03 Room Air  | 0.000 | 08:29 |
| 04 Std. Sol. | 0.301 | 08:30 |
| 05 Room Air  | 0.000 | 08:30 |
| 06 Std. Sol. | 0.302 | 08:31 |
| 07 Room Air  | 0.000 | 08:32 |


08 Sim Temp = 34.0°C

Simul Ser No = 000000

Std Sol No = 00000

County = 40

Oper No. = 107501



Operator Signature  
MARGY PEARSON

Remarks:

Linearity 0.292 AC

Form 106-18000

80-004199  
P. 81 of 86  
KLN


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004199  
Location = ROLL      8164.13.00 06/09  
05/12/2010      08:12

| Test             | AC    | Time  |
|------------------|-------|-------|
| 01 Room Air      | 0.000 | 08:13 |
| 02 *Subject Test | INT*  | 08:14 |
| 03 Room Air      | 0.000 | 08:15 |

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = NA  
Dr. Lic. = ND/DIS821456  
Simul Ser No = 000000  
Std Sol No = 00000  
County = 40      Oper No. = 107501

  
\_\_\_\_\_  
Operator Signature  
MARGY PEARSON

Remarks:

*Interferent AC + Acetone*

Form 106-I8000

80-004199  
P. 82 of 86  
KTH

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004199  
Location = ROLL                    8164.13.00 06/09  
05/12/2010                            08:16

| Test          | AC    | Time  |
|---------------|-------|-------|
| 01 Diagnostic | OK    | 08:17 |
| 02 Room Air   | RFI*  | 08:17 |
| 03 Room Air   | 0.000 | 08:18 |

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = NA  
Dr. Lic. = ND/DIS821456  
Simul Ser No = 000000  
Std Sol No = 00000  
County = 40                            Oper No. = 107501

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
MARGY PEARSON

Remarks:

RFI

Form 106-I8000

80-004199  
P. 83-186  
YH


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004199  
Location = ROLL      8164.13.00 06/09  
05/12/2010      08:38

DRY CAL CHECK

| Test        | AC    | Time  |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 08:38 |
| 02 Std. Gas | 0.080 | 08:39 |
| 03 Room Air | 0.000 | 08:39 |
| 04 Std. Gas | 0.081 | 08:39 |
| 05 Room Air | 0.000 | 08:40 |
| 06 Std. Gas | 0.080 | 08:40 |
| 07 Room Air | 0.000 | 08:41 |

Lot No = 659358  
Cyl No = 11  
Exp Date = 09/03/2011  
County = 40      Oper No. = 107501

  
\_\_\_\_\_  
Operator Signature  
MARGY PEARSON

Remarks:

Cal Check 0.080 AC

Form 106-I8000

80-004199  
P. 84 of 86  
JEM


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004199  
Location = ROLL      8164.13.00 06/09  
05/12/2010      08:41

DRY CAL CHECK

| Test        | AC    | Time  |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 08:42 |
| 02 Std. Gas | 0.081 | 08:42 |
| 03 Room Air | 0.000 | 08:43 |
| 04 Std. Gas | 0.080 | 08:43 |
| 05 Room Air | 0.000 | 08:43 |
| 06 Std. Gas | 0.080 | 08:44 |
| 07 Room Air | 0.000 | 08:44 |

Lot No = 659358  
Cyl No = 11  
Exp Date = 09/03/2011  
County = 40      Oper No. = 107501

  
\_\_\_\_\_  
Operator Signature  
MARGY PEARSON

Remarks:

0.080 Cal Check

Form 106-I8000

80-004199  
P. 85 of 86  
YAM



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004199  
Location = ROLL      8164.13.00 06/09  
05/12/2010      08:44

DRY CAL CHECK

| Test        | AC    | Time  |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 08:45 |
| 02 Std. Gas | 0.080 | 08:45 |
| 03 Room Air | 0.000 | 08:46 |
| 04 Std. Gas | 0.081 | 08:46 |
| 05 Room Air | 0.000 | 08:47 |
| 06 Std. Gas | 0.081 | 08:47 |
| 07 Room Air | 0.000 | 08:48 |

Lot No = 659358  
Cyl No = 11  
Exp Date = 09/03/2011  
County = 40      Oper No. = 107501

*Margy Pearson*

Operator Signature  
MARGY PEARSON

Remarks:

*0.080 AC Cal Check*

Form 106-I8000

*80-004199  
P. 86 of 86  
VM*