

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-004942

Location: TOXL

Check When Done

- A. Is the breath tube heated? ✓
- B. Is the diagnostic check complete? ✓
- C. Is the time, date, and year correct (re-set if necessary)? ✓
- D. Print test. (Attach test record.) ✓
- E. Low AC. Use ≤ 0.03 AC in ACA mode. (Attach test record.)
Sim SN DR3804 Lot No. 201501A Rep AC 0.015 ✓
- F. High AC. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN DR3377 Lot No. 14102 Rep AC 0.400 ✓
- G. Interferent check. Use 0.05% acetone in ethanol solution in ABA mode.
(Attach test record.)
Sim SN DR5135 Lot No. ICS#2 Solution 0.10AC + 0.05% ACE ✓
- H. RFI check. Run CMS mode. Key radio during CMS test.
(Attach test record.) ✓
- I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.) 0.080 AC

Lot No. 29614080A2 Cyl No. 28 Exp. Date 12/5/16
Test 1 0.081 Test 4 0.080 Test 7 0.081
Test 2 0.080 Test 5 0.080 Test 8 0.081
Test 3 0.080 Test 6 0.081 Test 9 0.081
Average 0.081

J. Tank pressure, Level 3 Func D Sub G

P: 460 psi Reg 475 psi

K. Remarks/Maintenance record: CALIBRATION ADJUSTMENT, ANNUAL INSPECTION

12 MAY 15
Date

Charles E. Egan
Field Inspector's Signature

13 May 15 ✓
Date

Deb Kashner
Reviewed by

TOXL
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-004942
05/12/2015 14:52:49

Auto Calibration

pg 1 of 2

<<<<<		3um	>>>>>		<<<<<		9um	>>>>>	
Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1 0.000 AC									
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)					
Sample #1	0.0990	(-0.0120)	0.1700	(-0.0180)					MILLI-Q H ₂ O
Sample #2	0.0730	(0.0420)	0.1540	(-0.0070)					LOT: NA
Sample #3	0.0540	(0.0850)	0.1560	(0.0060)					DR3453
Sample #4	0.0460	(0.1080)	0.1400	(0.0290)					
Avg % Abs	0.0577	(0.0783)	0.1500	(0.0093)					
STD DEV	0.0139	(0.0335)	0.0087	(0.0182)					
REL STD DEV	24.049	(42.768)	5.812	(195.322)					

Solution = 0.050 g/210L or 0.2381 mg/l, Samples = 4, Discarded = 1 0.050 AC									
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)					
Sample #1	0.9880	(-0.0130)	1.8890	(-0.0100)					LOT: 14020
Sample #2	0.9720	(0.0100)	1.8880	(0.0110)					EXP: 1.13.16
Sample #3	0.9410	(0.0260)	1.8780	(0.0060)					DR5114
Sample #4	0.9590	(0.0160)	1.9010	(-0.0020)					
Avg % Abs	0.9573	(0.0173)	1.8890	(0.0050)					
STD DEV	0.0156	(0.0081)	0.0115	(0.0066)					
REL STD DEV	1.626	(46.632)	0.611	(131.149)					

Solution = 0.100 g/210L or 0.4762 mg/l, Samples = 4, Discarded = 1 0.100 AC									
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)					
Sample #1	1.8880	(-0.0230)	3.6590	(-0.0230)					LOT: 14030
Sample #2	1.8680	(-0.0410)	3.6480	(-0.0120)					EXP: 12+ccc
Sample #3	1.8490	(-0.0640)	3.6270	(-0.0210)					1.20.16
Sample #4	1.8440	(-0.0850)	3.6240	(-0.0480)					DR5113
Avg % Abs	1.8537	(-0.0633)	3.6330	(-0.0270)					
STD DEV	0.0127	(0.0220)	0.0131	(0.0187)					
REL STD DEV	0.683	(34.749)	0.360	(69.389)					

Solution = 0.200 g/210L or 0.9524 mg/l, Samples = 4, Discarded = 1 0.200 AC									
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)					
Sample #1	3.6030	(0.0070)	6.8390	(0.0070)					LOT: 13220
Sample #2	3.6170	(0.0160)	6.8580	(0.0100)					EXP: 8.5.15
Sample #3	3.6200	(0.0230)	6.8570	(0.0110)					DR5189
Sample #4	3.6330	(0.0390)	6.8630	(0.0170)					
Avg % Abs	3.6233	(0.0260)	6.8593	(0.0127)					
STD DEV	0.0085	(0.0118)	0.0032	(0.0038)					
REL STD DEV	0.235	(45.345)	0.047	(29.889)					

Solution = 0.400 g/210L or 1.9048 mg/l, Samples = 4, Discarded = 1 0.400 AC									
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)					
Sample #1	6.8580	(0.0040)	12.7370	(0.0000)					LOT: 13113
Sample #2	6.8610	(0.0410)	12.7570	(0.0320)					EXP: 9.25.15
Sample #3	6.9130	(0.0310)	12.7890	(0.0290)					DR5191
Sample #4	6.9140	(0.0400)	12.7960	(0.0320)					
Avg % Abs	6.8960	(0.0373)	12.7807	(0.0310)					
STD DEV	0.0303	(0.0055)	0.0208	(0.0017)					
REL STD DEV	0.440	(14.752)	0.163	(5.587)					

Charles E. Eber

TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-004942
 05/12/2015 14:52:49

Auto Calibration

pg 2 of 2

<<<<< 3um >>>>>

 Zero Order Coef -118.58
 First Order Coef 2553.87
 Second Order Coef 32.53

<<<<< 9um >>>>>

 Zero Order Coef -182.95
 First Order Coef 1311.50
 Second Order Coef 15.10

 Act Fit Residual
 (g/210L) (g/210L) (g/210L)
 0.000 0.001 -0.0006
 0.050 0.049 0.0005
 0.100 0.099 0.0007
 0.200 0.201 -0.0008
 0.400 0.400 0.0002

 Act Fit Residual
 (g/210L) (g/210L) (g/210L)
 0.000 0.000 -0.0003
 0.050 0.049 0.0007
 0.100 0.100 -0.0004
 0.200 0.200 0.0000
 0.400 0.400 0.0000

<<<<< 3um >>>>>

 Solution = 0.100 g/210L or 0.4762 mg/l, Samples = 4, Discarded = 1

<<<<< 9um >>>>>

Sample
 Sample #1 4336.00
 Sample #2 4225.00
 Sample #3 4314.00
 Sample #4 4266.00
 Avg 4268.3335
 STD DEV 44.5459
 REL STD DEV 1.044
 H2O adjust (mg/l*10k) 493

4291.00 0.100 AC
 4263.00
 4298.00 LOT: 321140CEJ
 4275.00 32114100A1
 4278.6665
 17.7858 CYL: 025
 0.416
 483 EXP: 12.5.16

Atmospheric Pressure = 957

*****CALIBRATION SUCCESSFUL*****



Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004942
Location = TOXL 8164.13.00 06/09
05/12/2015 15:35

***** Printer Test *****

abcdefghijklmnopqrstuvwxy1234567890-|=|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^&*()_+?

abcdefghijklmnopqrstuvwxy1234567890-|=|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 29614080A2
Standard Cyl #? 28
Standard Expiration? 12/05/2016
Oper No? 666666

Flow Cal. Date: 08/21/2012
Slope 638
Intercept -487461

IR Calibration Date: 05/12/2015
 3um 9um

0th Coef(*100): -11858 -18295
1st Coef(*100): 255387 131150
2nd Coef(*100): 3253 1510
H2O adj(mg/l*10k): 493 483

***** Printer Test End *****



Operator Signature
CHARLES EDER

Remarks:

PRINT TEST

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004942
Location = TOXL 8164.13.00 06/09
05/12/2015 15:40

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:41
02 Std. Sol.	0.014	15:41
03 Room Air	0.000	15:42
04 Std. Sol.	0.015	15:43
05 Room Air	0.000	15:43
06 Std. Sol.	0.014	15:44
07 Room Air	0.000	15:44

08 Sim Temp = 34.0°C

Simul Ser No = DR3804
Std Sol No = 201501A
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *LOW AC*
 0.015 AC

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004942
Location = TOXL 8164.13.00 06/09
05/12/2015 15:46

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:47
02 Std. Sol.	0.405	15:48
03 Room Air	0.000	15:48
04 Std. Sol.	0.404	15:49
05 Room Air	0.000	15:49
06 Std. Sol.	0.406	15:50
07 Room Air	0.000	15:51

08 Sim Temp = 34.0°C

Simul Ser No = DR3377
Std Sol No = 14102
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *HIGH AC*
 0.400 AC

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004942
Location = TOXL 8164.13.00 06/09
05/12/2015 15:53

Test	AC	Time
01 Room Air	0.000	15:54
02 *Subject Test	INT*	15:54
03 Room Air	0.000	15:55

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = INTERFERENT
Dr. Lic. = ND/DIS821456
Lot No = 29614080A2
Cyl No = 28
Expiration Date = 12/05/2016
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *Interferent Check*
0.10AC + 0.05% ACETONE

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004942
Location = TOXL 8164.13.00 06/09
05/12/2015 15:56

Test	AC	Time
01 Diagnostic	OK	15:57
02 Room Air	RFI*	15:57
03 Room Air	0.000	15:58

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = RFI CHECK
Dr. Lic. = ND/DIS821456
Lot No = 29614080A2
Cyl No = 28
Expiration Date = 12/05/2016
County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks:

RFI CHECK


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004942
Location = TOXL 8164.13.00 06/09
05/12/2015 16:00

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:01
02 Std. Gas	0.081	16:01
03 Room Air	0.000	16:02
04 Std. Gas	0.080	16:02
05 Room Air	0.000	16:03
06 Std. Gas	0.080	16:03
07 Room Air	0.000	16:04

Lot No = 29614080A2
Cyl No = 28
Exp Date = 12/05/2016
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *Calibration Check*
0.080AC

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004942
Location = TOXL 8164.13.00 06/09
05/12/2015 16:04

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:04
02 Std. Gas	0.080	16:05
03 Room Air	0.000	16:05
04 Std. Gas	0.080	16:06
05 Room Air	0.000	16:06
06 Std. Gas	0.081	16:06
07 Room Air	0.000	16:07

Lot No = 29614080A2
Cyl No = 28
Exp Date = 12/05/2016
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *Calibration Check*
0.080 AC

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004942
Location = TOXL 8164.13.00 06/09
05/12/2015 16:08

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:08
02 Std. Gas	0.081	16:08
03 Room Air	0.000	16:09
04 Std. Gas	0.081	16:09
05 Room Air	0.000	16:10
06 Std. Gas	0.081	16:10
07 Room Air	0.000	16:11

Lot No = 29614080A2
Cyl No = 28
Exp Date = 12/05/2016
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC

Form 106-I8000