

INTOXILYZER 8000 INSPECTION

Serial No.: 80-004939

Location: TOXL

- | | Check When Done |
|---|-----------------|
| A. Is the warm-up time less than 20 minutes? | <u>✓</u> |
| B. Is a three-pronged grounded outlet used? | <u>✓</u> |
| C. Is the breath tube heated? | <u>✓</u> |
| D. Is the diagnostic check complete? | <u>✓</u> |
| E. Is the time, date, and year correct (re-set if necessary)? | <u>✓</u> |
| F. Print test. (Attach test record.) | <u>✓</u> |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)
Sim SN <u>DR5144</u> Lot No. <u>—</u> Rep AC <u>0.021</u> | <u>✓</u> |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN <u>DR5114</u> Lot No. <u>—</u> Rep AC <u>0.298</u> | <u>✓</u> |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.
(Attach test record.)
Sim SN <u>DR5186</u> Lot No. <u>—</u> Rep AC <u>0.100</u> | <u>✓</u> |
| J. RFI check. Run CMS mode. Key radio on first room air.
(Attach test record.) | <u>✓</u> |
| K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.) | |

Lot No. <u>659358</u>	Cyl No. <u>15</u>	Exp. Date <u>09/03/2011</u>
Test 1 <u>0.081</u>	Test 4 <u>0.081</u>	Test 7 <u>0.082</u>
Test 2 <u>0.081</u>	Test 5 <u>0.081</u>	Test 8 <u>0.082</u>
Test 3 <u>0.082</u>	Test 6 <u>0.081</u>	Test 9 <u>0.082</u>
Average <u>0.081</u>		

L. Tank pressure, Level 3 Func D Sub G P: 541 psi Reg 540 psi

M. Remarks/Maintenance record: Annual inspection

14 May 2010
Date

Deb Kashur
Field Inspector's Signature

15 May 2010
Date

Mary Pearson
Reviewed by

80-004939
P. 24/34 KHL

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004939
Location = TOXL 8164.13.00 06/09
05/14/2010 15:43

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:44
02 Std. Sol.	0.020	15:45
03 Room Air	0.000	15:45
04 Std. Sol.	0.019	15:46
05 Room Air	0.000	15:46
06 Std. Sol.	0.018	15:47
07 Room Air	0.000	15:48

08 Sim Temp = 34.0°C

Simul Ser No = 00000

Std Sol No = 000

County = 08

Oper No. = 777777

Deb Kashur

Operator Signature

DEB KASHUR

Remarks:

Low AC 0.021 AC

Form 106-I8000

80-004939
P. 28/34 JH

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004939
Location = TOXL 8164.13.00 06/09
05/14/2010 15:49

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:49
02 Std. Sol.	0.296	15:50
03 Room Air	0.000	15:51
04 Std. Sol.	0.296	15:51
05 Room Air	0.000	15:52
06 Std. Sol.	0.295	15:53
07 Room Air	0.000	15:53

08 Sim Temp = 34.0°C

Simul Ser No = 00000

Std Sol No = 000

County = 08

Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Linearity test 0.298 AC

Form 106-I8000

80-004939
P. 29/34 FH

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004939
Location = TOXL 8164.13.00 06/09
05/14/2010 16:05

Test	AC	Time
01 Room Air	0.000	16:05
02 *Subject Test	INT*	16:06
03 Room Air	0.000	16:06

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 659358
Cyl No = 15
Expiration Date = 09/03/2011
County = 08 Oper No. = 777777

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

Interferent check 0.100 AC acetone

Form 106-I8000

80-004939
p. 30/34 KM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004939
Location = TOXL 8164.13.00 06/09
05/14/2010 16:07

Test	AC	Time
01 Diagnostic	OK	16:07
02 Room Air	RFI*	16:07
03 Room Air	0.000	16:08

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 659358
Cyl No = 15
Expiration Date = 09/03/2011
County = 08 Oper No. = 777777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

RFI check

Form 106-I8000

80-004939
P. 3/34 JH

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004939
Location = TOXL 8164.13.00 06/09
05/14/2010 16:08

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:09
02 Std. Gas	0.081	16:09
03 Room Air	0.000	16:09
04 Std. Gas	0.081	16:10
05 Room Air	0.000	16:10
06 Std. Gas	0.082	16:11
07 Room Air	0.000	16:11

Lot No = 659358
Cyl No = 15
Exp Date = 09/03/2011
County = 08 Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080 AC

Form 106-I8000

80-004939
P. 37/34 KR

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004939
Location = TOXL 8164.13.00 06/09
05/14/2010 16:11

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:12
02 Std. Gas	0.081	16:12
03 Room Air	0.000	16:13
04 Std. Gas	0.081	16:13
05 Room Air	0.000	16:14
06 Std. Gas	0.081	16:14
07 Room Air	0.000	16:14

Lot No = 659358

Cyl No = 15

Exp Date = 09/03/2011

County = 08

Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080 AC

Form 106-I8000

80-004939
P. 33/34 F14

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004939
Location = TOXL 8164.13.00 06/09
05/14/2010 16:15

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:15
02 Std. Gas	0.082	16:16
03 Room Air	0.000	16:16
04 Std. Gas	0.082	16:16
05 Room Air	0.000	16:17
06 Std. Gas	0.082	16:17
07 Room Air	0.000	16:18

Lot No = 659358
Cyl No = 15
Exp Date = 09/03/2011
County = 08 Oper No. = 777777

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080 AC

Form 106-18000

80-004939
P. 34/34 FM