



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59281 (5-2011)

Serial Number <i>80-004199</i>	Instrument Location <i>NAPO</i>
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level 1, Function P) ✓
 - B. ACA Test (Level 1, Function C) ✓
 - C. RFI Test (Use CMS Mode; Key Radio During First Room Air) ✓
- 7. Repair and/or Maintenance Performed (if Any): _____
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>[Signature]</i>	Date <i>5-15-12</i>
Reviewed By (Crime Laboratory Use Only) <i>Deb Kashner</i>	Date <i>15 May 2012</i>

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004199
Location = NAPO 8164.13.00 06/09
05/15/2012 11:28

***** Printer Test *****

abcdefghijklmnopqrstuvwxy1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&*()_+?

abcdefghijklmnopqrstuvwxy1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 23411080A1
Standard Cyl #? 29
Standard Expiration? 10/01/2013
Oper No? 120620

Flow Cal. Date: 07/05/2011
Slope 672
Intercept -604945

IR Calibration Date: 04/26/2012

	3um	9um
0th Coef(*100):	-20552	-25225
1st Coef(*100):	263837	131160
2nd Coef(*100):	3599	1386
H2O adj(mg/l*10k):	516	667

***** Printer Test End *****



Operator Signature
GORDON LECLAIR, JR.

Remarks: Replaced Intox - 80-004442.

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004199
Location = NAPO 8164.13.00 06/09
05/15/2012 11:42

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:42
02 Std. Gas	0.077	11:43
03 Room Air	0.000	11:43
04 Std. Gas	0.078	11:44
05 Room Air	0.000	11:44
06 Std. Gas	0.076	11:45
07 Room Air	0.000	11:45

Lot No = 23411080A1
Cyl No = 29
Exp Date = 10/01/2013
County = 24 Oper No. = 120620



Operator Signature
GORDON LECLAIR, JR.

Remarks: *ACA*

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004199
Location = NAPO 8164.13.00 06/09
05/15/2012 11:46

Test	AC	Time
01 Diagnostic	OK	11:48
02 Room Air	RFI*	11:48
03 Room Air	0.000	11:49

*Invalid Test
Inhibited - RFI

Sub Name = TEST, TEST TEST
Sub DOB = 09/10/1960
Sub Sex = Male Weight = 200
Test = DUI Cit = 0000000
Dr. Lic. = ND/NA
Lot No = 23411080A1
Cyl No = 29
Expiration Date = 10/01/2013
County = 24 Oper No. = 120620

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
GORDON LECLAIR, JR.

Remarks: ~~RFI~~ test p.l.c.