

INTOXILYZER 8000 INSPECTION

Serial No.: 80-004206

Location: WAHP (WAHPETON)

- | | Check When Done |
|---|-----------------|
| A. Is the warm-up time less than 20 minutes? | <u>✓</u> |
| B. Is a three-pronged grounded outlet used? | <u>✓</u> |
| C. Is the breath tube heated? | <u>✓</u> |
| D. Is the diagnostic check complete? | <u>✓</u> |
| E. Is the time, date, and year correct (re-set if necessary)? | <u>✓</u> |
| F. Print test. (Attach test record.) | <u>✓</u> |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)
Sim SN <u>DR3455</u> Lot No. <u>12210</u> Rep AC <u>0.020</u> | <u>✓</u> |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN <u>DR3847</u> Lot No. <u>13020</u> Rep AC <u>0.300</u> | <u>✓</u> |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.
(Attach test record.)
Sim SN <u>DR5191</u> Lot No. <u>ICS#1</u> Rep AC <u>0.100AC + 0.05% ACETONE</u> | <u>✓</u> |
| J. RFI check. Run CMS mode. Key radio on first room air.
(Attach test record.) | <u>✓</u> |

K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.)

Lot No. <u>2341108⁰AA1</u> <small>CBE</small>	Cyl No. <u>58</u>	Exp. Date <u>10/01/2013</u>
Test 1 <u>0.080</u>	Test 4 <u>0.080</u>	Test 7 <u>0.080</u>
Test 2 <u>0.081</u>	Test 5 <u>0.080</u>	Test 8 <u>0.081</u>
Test 3 <u>0.081</u>	Test 6 <u>0.081</u>	Test 9 <u>0.080</u>
Average <u>0.080</u>		

L. Tank pressure, Level 3 Func D Sub G P: 796 psi Reg 800 psi

M. Remarks/Maintenance record: REPLACED EXTERNAL BATTERY; ANNUAL INSPECTION
REPLACED O-RING SIMULATOR RETURN
15 MAY 13
Date
116 May 13
Date

Charles E. Edin
Field Inspector's Signature
Deb Kashner
Reviewed by

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004206
Location = WAHP 8164.13.00 06/09
05/15/2013 12:06

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%\$^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%\$^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 23411080A1
Standard Cyl #? 58
Standard Expiration? 10/01/2013
Oper No? 666666

Flow Cal. Date: 06/10/2009
Slope 690
Intercept -645883

IR Calibration Date: 06/07/2010
 3um 9um

0th Coef(*100): -17579 -21810
1st Coef(*100): 254057 127149
2nd Coef(*100): 4591 1650
H2O adj(mg/l*10k): 756 795

***** Printer Test End *****



Operator Signature
CHARLES EDER

Remarks:

PRINT TEST

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004206
Location = WAHP 8164.13.00 06/09
05/15/2013 12:19

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:20
02 Std. Sol.	0.018	12:21
03 Room Air	0.000	12:21
04 Std. Sol.	0.018	12:22
05 Room Air	0.000	12:23
06 Std. Sol.	0.019	12:23
07 Room Air	0.000	12:24

08 Sim Temp = 34.0°C

Simul Ser No = DR3455

Std Sol No = 12210

County = 39

Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

Low AC

0.020AC (Sim. Sol.)

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004206
Location = WAHP 8164.13.00 06/09
05/15/2013 12:25

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:26
02 Std. Sol.	0.305	12:26
03 Room Air	0.000	12:27
04 Std. Sol.	0.307	12:28
05 Room Air	0.000	12:28
06 Std. Sol.	0.308	12:29
07 Room Air	0.000	12:29

08 Sim Temp = 34.0°C


Simul Ser No = DR3847

Std Sol No = 13020

County = ~~08~~ 39

Oper No. = 666666

CEE



Operator Signature
CHARLES EDER

Remarks:

*LINEARITY TEST
0.300 AC (SIM. SOL.)*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004206
Location = WAHP 8164.13.00 06/09
05/15/2013 12:15

Test	AC	Time
01 Room Air	0.000	12:16
02 *Subject Test	INT*	12:16
03 Room Air	0.000	12:17

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = MNBVCXZLKJH0987
Dr. Lic. = ND/DIS821456
Lot No = 23411080A1
Cyl No = 58
Expiration Date = 10/01/2013
County = 39 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *INTERFERENT CHECK*
0.100 AC + 0.05% ACETONE

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004206
Location = WAHP 8164.13.00 06/09
05/15/2013 12:12

Test	AC	Time
01 Diagnostic	OK	12:13
02 Room Air	RFI*	12:13
03 Room Air	0.000	12:14

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = QWERTYUIOPASDFG
Dr. Lic. = ND/DIS821456
Lot No = 23411080A1
Cyl No = 58
Expiration Date = 10/01/2013
County = 39 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks: *RFI CHECK*


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004206
Location = WAHP 8164.13.00 06/09
05/15/2013 12:31

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:31
02 Std. Gas	0.080	12:31
03 Room Air	0.000	12:32
04 Std. Gas	0.081	12:32
05 Room Air	0.000	12:33
06 Std. Gas	0.081	12:33
07 Room Air	0.000	12:34

Lot No = 23411080A1
Cyl No = 58
Exp Date = 10/01/2013
County = 39 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC (E.G.S.)


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004206
Location = WAHP 8164.13.00 06/09
05/15/2013 12:34

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:35
02 Std. Gas	0.080	12:35
03 Room Air	0.000	12:36
04 Std. Gas	0.080	12:36
05 Room Air	0.000	12:36
06 Std. Gas	0.081	12:37
07 Room Air	0.000	12:37

Lot No = 23411080A1
Cyl No = 58
Exp Date = 10/01/2013
County = 39 Oper No. = 666666


Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC (E.G.S.)

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004206
Location = WAHP 8164.13.00 06/09
05/15/2013 12:37

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:38
02 Std. Gas	0.080	12:38
03 Room Air	0.000	12:39
04 Std. Gas	0.081	12:39
05 Room Air	0.000	12:40
06 Std. Gas	0.080	12:40
07 Room Air	0.000	12:41

Lot No = 23411080A1

Cyl No = 58

Exp Date = 10/01/2013

County = 39

Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC (E.G.S.)

Form 106-I8000