

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-005360

Location: ASHL (McINTOSH COUNTY S.O.)

- Check When Done
- A. Is the breath tube heated? ✓
 - B. Is the diagnostic check complete? ✓
 - C. Is the time, date, and year correct (re-set if necessary)? ✓
 - D. Print test. (Attach test record.) ✓
 - E. Low AC. Use ≤ 0.03 AC in ACA mode. (Attach test record.)
Sim SN DR5191 Lot No. 15160 Rep AC 0.020 ✓
 - F. High AC. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN DR3847 Lot No. 16060 Rep AC 0.300 ✓
 - G. Interferent check. Use 0.05% acetone in ethanol solution in ABA mode.
(Attach test record.)
Sim SN DR5190 Lot No. ICS 4 Solution 0.10AC + 0.05% ACETONE ✓
 - H. RFI check. Run CMS mode. Key radio during CMS test.
(Attach test record.) ✓
 - I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.) 0.080 AC

Lot No. 29614080AR Cyl No. 007 Exp. Date 12/05/2016

Test 1 <u>0.083</u>	Test 4 <u>0.082</u>	Test 7 <u>0.083</u>
Test 2 <u>0.081</u>	Test 5 <u>0.082</u>	Test 8 <u>0.083</u>
Test 3 <u>0.082</u>	Test 6 <u>0.083</u>	Test 9 <u>0.082</u>

Average 0.082

J. Tank pressure, Level 3 Func D Sub G P: 1015 psi Reg ~1000 psi

K. Remarks/Maintenance record: ANNUAL INSPECTION

5/13/2016
Date

Charles E. Eick
Field Inspector's Signature

16 May 16 ✓
Date

Deb Shanaver
Reviewed by

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005360
Location = ASHL 8164.13.00 06/09
05/13/2016 11:54

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 29614080A2
Standard Cyl #? 7
Standard Expiration? 12/05/2016
Oper No? 666666

Flow Cal. Date: 08/17/2011
Slope 639
Intercept -343114

IR Calibration Date: 08/17/2011
 3um 9um

0th Coef(*100): -17266 -25747
1st Coef(*100): 263773 142804
2nd Coef(*100): 2606 1406
H2O adj(mg/l*10k): 499 577

***** Printer Test End *****



Operator Signature
CHARLES EDER

Remarks: PRINT TEST

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005360
Location = ASHL 8164.13.00 06/09
05/13/2016 12:22

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:23
02 Std. Sol.	0.017	12:23
03 Room Air	0.000	12:24
04 Std. Sol.	0.018	12:25
05 Room Air	0.000	12:25
06 Std. Sol.	0.017	12:26
07 Room Air	0.000	12:26

08 Sim Temp = 34.0°C

Simul Ser No = DR5191

Std Sol No = 15160

County = 26

Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

Low AC
0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005360
Location = ASHL 8164.13.00 06/09
05/13/2016 12:12

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:12
02 Std. Sol.	0.298	12:13
03 Room Air	0.000	12:14
04 Std. Sol.	0.300	12:14
05 Room Air	0.000	12:15
06 Std. Sol.	0.300	12:16
07 Room Air	0.000	12:16


08 Sim Temp = 34.0°C

Simul Ser No = DR3847

Std Sol No = 16060

County = 26

Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

HIGH AC
0.300 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005360
Location = ASHL 8164.13.00 06/09
05/13/2016 12:02

Test	AC	Time
01 Room Air	0.000	12:03
02 *Subject Test	INT*	12:03
03 Room Air	0.000	12:04

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = INTERFERENT
Dr. Lic. = ND/DIS821456
Lot No = 29614080A2
Cyl No = 7
Expiration Date = 12/05/2016
County = 26 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

Interferent Check
0.10 AC + 0.05% ACETONE

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005360
Location = ASHL 8164.13.00 06/09
05/13/2016 12:00

Test	AC	Time
01 Diagnostic	OK	12:01
02 Room Air	RFI*	12:01
03 Room Air	0.000	12:02

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = RFI CHECK
Dr. Lic. = ND/DIS821456
Lot No = 29614080A2
Cyl No = 7
Expiration Date = 12/05/2016
County = 26 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks:

RFI CHECK


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005360
Location = ASHL 8164.13.00 06/09
05/13/2016 12:31

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:32
02 Std. Gas	0.082	12:32
03 Room Air	0.000	12:33
04 Std. Gas	0.082	12:33
05 Room Air	0.000	12:34
06 Std. Gas	0.083	12:34
07 Room Air	0.000	12:34

Lot No = 29614080A2
Cyl No = 7
Exp Date = 12/05/2016
County = 26 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005360
Location = ASHL 8164.13.00 06/09
05/13/2016 12:35

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:35
02 Std. Gas	0.083	12:36
03 Room Air	0.000	12:36
04 Std. Gas	0.083	12:36
05 Room Air	0.000	12:37
06 Std. Gas	0.082	12:37
07 Room Air	0.000	12:38

Lot No = 29614080A2
Cyl No = 7
Exp Date = 12/05/2016
County = 26 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC