

INTOXILYZER 8000 INSPECTION

Serial No.: 80-002666

Location: TOXL

- | | Check When Done |
|--|-----------------|
| A. Is the warm-up time less than 20 minutes? | <u>✓</u> |
| B. Is a three-pronged grounded outlet used? | <u>✓</u> |
| C. Is the breath tube heated? | <u>✓</u> |
| D. Is the diagnostic check complete? | <u>✓</u> |
| E. Is the time, date, and year correct (re-set if necessary)? | <u>✓</u> |
| F. Print test. (Attach test record.) | <u>✓</u> |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)
Sim SN <u>DR5144</u> Lot No. <u>—</u> Rep AC <u>0.021</u> | <u>✓</u> |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN <u>DR3847</u> Lot No. <u>—</u> Rep AC <u>0.294</u> | <u>✓</u> |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.
(Attach test record.)
Sim SN <u>DR3452</u> Lot No. <u>—</u> Rep AC <u>0.10AC + Acetone</u> | <u>✓</u> |
| J. RFI check. Run CMS mode. Key radio on first room air.
(Attach test record.) | <u>✓</u> |
| K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.) | |

Lot No. <u>15410080A1</u>	Cyl No. <u>21</u>	Exp. Date <u>7-1-12</u>
Test 1 <u>0.079</u>	Test 4 <u>0.079</u>	Test 7 <u>0.079</u>
Test 2 <u>0.080</u>	Test 5 <u>0.080</u>	Test 8 <u>0.079</u>
Test 3 <u>0.080</u>	Test 6 <u>0.080</u>	Test 9 <u>0.080</u>
Average <u>0.080</u>		

L. Tank pressure, Level 3 Func D Sub G P: 226 psi Reg 240 psi

M. Remarks/Maintenance record: Annual Inspection

5-18-11
Date

[Signature]
Field Inspector's Signature

19 May 2011
Date

Deb Kashur
Reviewed by

80-002666
P. 81 of 89
AM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002666
Location = TOXL 8164.13.00 06/09
05/18/2011 13:53

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 15410080A1
Standard Cyl #? 21
Standard Expiration? 07/01/2012
Oper No? 555555

Flow Cal. Date: 09/10/2007
Slope 677
Intercept -603196

IR Calibration Date: 09/10/2007
 3um 9um

0th Coef(*100): -14148 -8312
1st Coef(*100): 274152 135499
2nd Coef(*100): 1387 855
H2O adj(mg/l*10k): 430 237

***** Printer Test End *****



Operator Signature
ASHLEY AYD

Remarks: *Print test*

*80-002666
P. 82 of 89
AM*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002666
Location = TOXL 8164.13.00 06/09
05/18/2011 13:55

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:56
02 Std. Sol.	0.021	13:57
03 Room Air	0.000	13:57
04 Std. Sol.	0.021	13:58
05 Room Air	0.000	13:59
06 Std. Sol.	0.021	13:59
07 Room Air	0.000	14:00

08 Sim Temp = 34.0°C

Simul Ser No = 0

Std Sol No = 0

County = 08

Oper No. = 555555



Operator Signature
ASHLEY AYD

Remarks:

Low AC 0.021AC

Form 106-I8000

80-002666
P. 83 of 89
VM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002666
Location = TOXL 8164.13.00 06/09
05/18/2011 14:01

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:02
02 Std. Sol.	0.294	14:02
03 Room Air	0.000	14:03
04 Std. Sol.	0.294	14:04
05 Room Air	0.000	14:04
06 Std. Sol.	0.295	14:05
07 Room Air	0.000	14:05


08 Sim Temp = 34.0°C

Simul Ser No = 0

Std Sol No = 0

County = 08

Oper No. = 555555



Operator Signature
ASHLEY AYD

Remarks: *Linearity test 0.294 AC*

Form 106-I8000

80-002666
P. 84 of 89
VM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002666
Location = TOXL 8164.13.00 06/09
05/18/2011 14:07

Test	AC	Time
01 Room Air	0.000	14:08
02 *Subject Test	INT*	14:08
03 Room Air	0.000	14:09

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Simul Ser No = 0
Std Sol No = 0
County = 08 Oper No. = 555555



Operator Signature
ASHLEY AYD

Remarks: *Interferent check 0.10AC + Acetone*

Form 106-I8000

80-002666
P. 85 of 89
FM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002666
Location = TOXL 8164.13.00 06/09
05/18/2011 14:09

Test	AC	Time
01 Diagnostic	OK	14:10
02 Room Air	RFI*	14:10
03 Room Air	0.000	14:11

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Simul Ser No = 0
Std Sol No = 0
County = 08 Oper No. = 555555

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
ASHLEY AYD

Remarks: RFI

Form 106-I8000

80-002666
P. 86 of 89
AM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002666
Location = TOXL 8164.13.00 06/09
05/18/2011 ~~13:39~~ ^{13:42}
13:39

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:39
02 Std. Gas	0.079	13:40
03 Room Air	0.000	13:40
04 Std. Gas	0.080	13:41
05 Room Air	0.000	13:41
06 Std. Gas	0.080	13:42
07 Room Air	0.000	13:42

Lot No = 15410080A1
Cyl No = 21
Exp Date = 07/01/2012
County = 08 Oper No. = 555555



Operator Signature
ASHLEY AYD

Remarks: Calibration check 0.080AC

Form 106-I8000

80-002666
P. 87 of 89
FM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002666
Location = TOXL 8164.13.00 06/09
05/18/2011 13:44

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:45
02 Std. Gas	0.079	13:45
03 Room Air	0.000	13:46
04 Std. Gas	0.080	13:46
05 Room Air	0.000	13:46
06 Std. Gas	0.080	13:47
07 Room Air	0.000	13:47

Lot No = 15410080A1
Cyl No = 21
Exp Date = 07/01/2012
County = 08 Oper No. = 555555



Operator Signature
ASHLEY AYD

Remarks: Calibration check 0.080AC

Form 106-I8000

80-002666
P. 88 + 89
AM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-002666
Location = TOXL 8164.13.00 06/09
05/18/2011 13:49

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:49
02 Std. Gas	0.079	13:50
03 Room Air	0.000	13:50
04 Std. Gas	0.079	13:51
05 Room Air	0.000	13:51
06 Std. Gas	0.080	13:51
07 Room Air	0.000	13:52

Lot No = 15410080A1
Cyl No = 21
Exp Date = 07/01/2012
County = 08 Oper No. = 555555



Operator Signature
ASHLEY AYD

Remarks: Calibration check 0.080 AC

Form 106-I8000

80-002666
P. 89489
FM