



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59281 (5-2011)

Serial Number 80-004937	Instrument Location WFPO
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level 1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. RFI Test (Use CMS Mode; Key Radio During First Room Air).
- 7. Repair and/or Maintenance Performed (if Any): _____
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>[Signature]</i>	Date 5-18-12
Reviewed By (Crime Laboratory Use Only) <i>Deb Kashner</i>	Date 22 May 2012

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004937
Location = WFPD 8164.13.00 06/09
05/18/2012 11:50

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 05911080A1
Standard Cyl #? 34
Standard Expiration? 02/01/2013
Oper No? 039617

Flow Cal. Date: 01/18/2010
Slope 730
Intercept -779991

IR Calibration Date: 05/24/2011
 3um 9um

0th Coef(*100): -12537 -15884
1st Coef(*100): 262613 132554
2nd Coef(*100): 1600 951
H2O adj(mg/l*10k): 529 537

***** Printer Test End *****



Operator Signature
PETER NIELSEN

Remarks:

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004937
Location = WFPD 8164.13.00 06/09
05/18/2012 11:46

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:47
02 Std. Gas	0.079	11:47
03 Room Air	0.000	11:48
04 Std. Gas	0.079	11:48
05 Room Air	0.000	11:49
06 Std. Gas	0.079	11:49
07 Room Air	0.000	11:50

Lot No = 05911080A1
Cyl No = 34
Exp Date = 02/01/2013
County = 09 Oper No. = 039617



Operator Signature
PETER NIELSEN

Remarks:

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004937
Location = WFPD 8164.13.00 06/09
05/18/2012 11:54

Test	AC	Time
01 Diagnostic	OK	11:56
02 Room Air	RFI*	11:56
03 Room Air	0.000	11:57

*Invalid Test
Inhibited - RFI

Sub Name = TEST, TEST T
Sub DOB = 01/01/1900
Sub Sex = Male Weight = 150
Test = DUI Cit = TEST
Dr. Lic. = ND/TEST
Lot No = 05911080A1
Cyl No = 34
Expiration Date = 02/01/2013
County = 09 Oper No. = 039617

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
PETER NIELSEN

Remarks: