

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-004203

Location: RUGB

- | | Check When Done |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|
| A. Is the breath tube heated? | <u>✓</u> |
| B. Is the diagnostic check complete? | <u>✓</u> |
| C. Is the time, date, and year correct (re-set if necessary)? | <u>✓</u> |
| D. Print test. (Attach test record.) | <u>✓</u> |
| E. Low AC. Use <0.03 AC in ACA mode. (Attach test record.) Sim SN <u>DR5115</u> Lot No. <u>12111</u> Rep AC <u>0.020</u> | <u>✓</u> |
| F. High AC. Use ≥ 0.25 AC in ACA mode. (Attach test record.) Sim SN <u>DR5143</u> Lot No. <u>13080</u> Rep AC <u>0.300</u> | <u>✓</u> |
| G. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode. (Attach test record.) Sim SN <u>DR3375</u> Lot No. <u>IRS #2</u> Rep AC <u>0.100 AC + 0.05% ACETONE</u> | <u>✓</u> |
| H. RFI check. Run CMS mode. Key radio on first room air. (Attach test record.) | <u>✓</u> |
| I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode. (Attach test record.) | |

| | | |
|---------------------------|---------------------|-----------------------------|
| Lot No. <u>16512680AA</u> | Cyl No. <u>020</u> | Exp. Date <u>08/01/2014</u> |
| Test 1 <u>0.082</u> | Test 4 <u>0.083</u> | Test 7 <u>0.082</u> |
| Test 2 <u>0.083</u> | Test 5 <u>0.083</u> | Test 8 <u>0.082</u> |
| Test 3 <u>0.082</u> | Test 6 <u>0.082</u> | Test 9 <u>0.081</u> |
| Average <u>0.082</u> | | |

J. Tank pressure, Level 3 Func D Sub G P: 1011 psi Reg 1000 psi

K. Remarks/Maintenance record: REPLACED BREATH TUBE AND SIMULATOR

RETURN O-RINGS, ANNUAL INSPECTION

| | |
|--------------------|-----------------------------|
| <u>15 MAY 14</u> | <u>Charles E. Elmer</u> |
| Date | Field Inspector's Signature |
| <u>19 May 14 ✓</u> | <u>Deb Kashur</u> |
| Date | Reviewed by |

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004203
Location = RUGB 8164.13.00 06/09
05/15/2014 15:22

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopQRSTUVWXYZ!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopQRSTUVWXYZ!@#\$\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 16512080A2
Standard Cyl #? 20
Standard Expiration? 08/01/2014
Oper No? 666666

Flow Cal. Date: 01/28/2010
Slope 719
Intercept -711903

IR Calibration Date: 01/28/2010
 3um 9um

0th Coef(*100): -17469 -28710
1st Coef(*100): 264216 135022
2nd Coef(*100): 2778 1081
H2O adj(mg/l*10k): 430 605

***** Printer Test End *****



Operator Signature
CHARLES EDER

Remarks:

PRINT TEST

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004203
Location = RUGB 8164.13.00 06/09
05/15/2014 15:31

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 15:32 |
| 02 Std. Sol. | 0.017 | 15:33 |
| 03 Room Air | 0.000 | 15:34 |
| 04 Std. Sol. | 0.017 | 15:35 |
| 05 Room Air | 0.000 | 15:35 |
| 06 Std. Sol. | 0.018 | 15:36 |
| 07 Room Air | 0.000 | 15:36 |

08 Sim Temp = 34.0°C

Simul Ser No = DR5115

Std Sol No = 12111

County = 35

Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *LOW AC*

0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004203
Location = RUGB 8164.13.00 06/09
05/15/2014 15:37

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 15:38 |
| 02 Std. Sol. | 0.288 | 15:38 |
| 03 Room Air | 0.000 | 15:39 |
| 04 Std. Sol. | 0.293 | 15:40 |
| 05 Room Air | 0.000 | 15:40 |
| 06 Std. Sol. | 0.295 | 15:41 |
| 07 Room Air | 0.000 | 15:42 |

08 Sim Temp = 34.0°C

Simul Ser No = DR5143

Std Sol No = 13080

County = 35

Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

HIGH AC
0.300 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004203
Location = RUGB 8164.13.00 06/09
05/15/2014 15:25

| Test | AC | Time |
|------------------|-------|-------|
| 01 Room Air | 0.000 | 15:26 |
| 02 *Subject Test | INT* | 15:26 |
| 03 Room Air | 0.000 | 15:27 |

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = INTERFERENT
Dr. Lic. = ND/DIS821456
Lot No = 16512080A2
Cyl No = 20
Expiration Date = 08/01/2014
County = 35 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

INTERFERENT CHECK
0.100 AC + 0.05% ACETONE

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004203
Location = RUGB 8164.13.00 06/09
05/15/2014 15:22

| Test | AC | Time |
|---------------|-------|-------|
| 01 Diagnostic | OK | 15:24 |
| 02 Room Air | RFI* | 15:24 |
| 03 Room Air | 0.000 | 15:24 |

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = RFI CHECK
Dr. Lic. = ND/DIS821456
Lot No = 16512080A2
Cyl No = 20
Expiration Date = 08/01/2014
County = 35 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks:

RFI CHECK

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004203
Location = RUGB 8164.13.00 06/09
05/15/2014 15:44

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 15:44 |
| 02 Std. Gas | 0.083 | 15:45 |
| 03 Room Air | 0.000 | 15:45 |
| 04 Std. Gas | 0.083 | 15:46 |
| 05 Room Air | 0.000 | 15:46 |
| 06 Std. Gas | 0.082 | 15:46 |
| 07 Room Air | 0.000 | 15:47 |

Lot No = 16512080A2
Cyl No = 20
Exp Date = 08/01/2014
County = 35 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION
CHECK
0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004203
Location = RUGB 8164.13.00 06/09
05/15/2014 15:47

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 15:48 |
| 02 Std. Gas | 0.082 | 15:48 |
| 03 Room Air | 0.000 | 15:49 |
| 04 Std. Gas | 0.082 | 15:49 |
| 05 Room Air | 0.000 | 15:50 |
| 06 Std. Gas | 0.081 | 15:50 |
| 07 Room Air | 0.000 | 15:50 |

Lot No = 16512080A2
Cyl No = 20
Exp Date = 08/01/2014
County = 35 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION
CHECK
0.080 AC