

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-005955

Location: DUNN - Manning

- |  | Check When Done |
|--|-----------------|
| A. Is the breath tube heated?  | <u>✓</u>        |
| B. Is the diagnostic check complete?   | <u>✓</u>        |
| C. Is the time, date, and year correct (re-set if necessary)?  | <u>✓</u>        |
| D. Print test. (Attach test record.)   | <u>✓</u>        |
| E. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR5115</u> Lot No. <u>12111</u> Rep AC <u>0.020</u>  | <u>✓</u>        |
| F. High AC. Use ≥ 0.25 AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR5143</u> Lot No. <u>13080</u> Rep AC <u>0.300</u>  | <u>✓</u>        |
| G. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.<br>(Attach test record.)<br>Sim SN <u>DR3375</u> Lot No. <u>FCS 2</u> Rep AC <u>0.100 AC + 0.05% acetone</u> | <u>✓</u>        |
| H. RFI check. Run CMS mode. Key radio on first room air.<br>(Attach test record.)  | <u>✓</u>        |
| I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.<br>(Attach test record.)  |                 |

Lot No. <u>16512080A2</u>	Cyl No. <u>13</u>	Exp. Date <u>08/01/2014</u>
Test 1 <u>0.081</u>	Test 4 <u>0.082</u>	Test 7 <u>0.082</u>
Test 2 <u>0.081</u>	Test 5 <u>0.082</u>	Test 8 <u>0.082</u>
Test 3 <u>0.081</u>	Test 6 <u>0.081</u>	Test 9 <u>0.081</u>
Average <u>0.081</u>		

J. Tank pressure, Level 3 Func D Sub G P: 873 psi Reg 870 psi

K. Remarks/Maintenance record: Annual inspection

16 May 2014  
Date

19 MAY 2014  
Date

Deb Kashur  
Field Inspector's Signature

Charles E. Edr  
Reviewed by

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-005955  
Location = DUNN                    8164.13.00 06/09  
05/16/2014                            09:47

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%&^\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%&^\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?          1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            16512080A2  
Standard Cyl #?            13  
Standard Expiration?        08/01/2014  
Oper No?                    777777

Flow Cal. Date:            06/21/2013  
Slope                      757  
Intercept                  -837385

IR Calibration Date:        06/21/2013  
                                    3um                            9um

-----  
0th Coef(\*100):            -19897                    -22878  
1st Coef(\*100):            268381                    131979  
2nd Coef(\*100):            5157                      1473  
H2O adj(mg/l\*10k):        717                        600

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
DEB KASHUR

Remarks:



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005955  
Location = DUNN      8164.13.00 06/09  
05/16/2014      10:15

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:16
02 Std. Sol.	0.020	10:17
03 Room Air	0.000	10:17
04 Std. Sol.	0.020	10:18
05 Room Air	0.000	10:18
06 Std. Sol.	0.021	10:19
07 Room Air	0.000	10:20

08 Sim Temp = 34.0°C

Simul Ser No = DR5115  
Std Sol No = 12111  
County = 13      Oper No. = 777777

Deb Kashur  
Operator Signature  
DEB KASHUR

Remarks:

Low AC      0.020 AC

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005955  
Location = DUNN      8164.13.00 06/09  
05/16/2014      10:21

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:21
02 Std. Sol.	0.305	10:22
03 Room Air	0.000	10:23
04 Std. Sol.	0.305	10:23
05 Room Air	0.000	10:24
06 Std. Sol.	0.306	10:25
07 Room Air	0.000	10:25

08 Sim Temp = 34.0°C

Simul Ser No = DR5143  
Std Sol No = 13080  
County = 13      Oper No. = 777777

*Deb Kashur*

Operator Signature  
DEB KASHUR

Remarks:

*High AC 0.300 AC*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-005955  
Location = DUNN                    8164.13.00 06/09  
05/16/2014                            10:28

Test	AC	Time
01 Room Air	0.000	10:29
02 *Subject Test	INT*	10:29
03 Room Air	0.000	10:30

\*Invalid Test  
Interferent Detected

Sub Name = INTERFERENT, NA NA  
Sub DOB = 08/01/1975  
Sub Sex = Unknown                    Weight = NA  
Test = DUI                              Cit = NA  
Dr. Lic. = ND/NA  
Lot No = 16512080A2  
Cyl No = 13  
Expiration Date = 08/01/2014  
County = 13                              Oper No. = 777777



Operator Signature  
DEB KASHUR

Remarks:

*Interferent check 0.100AC + 0.05% acetone*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-005955  
Location = DUNN                    8164.13.00 06/09  
05/16/2014                            10:05

Test	AC	Time
01 Diagnostic	OK	10:06
02 Room Air	RFI*	10:06
03 Room Air	0.000	10:07

\*Invalid Test  
Inhibited - RFI

Sub Name = RFI, NA NA  
Sub DOB = 08/01/1975  
Sub Sex = Unknown                    Weight = NA  
Test = DUI                              Cit = NA  
Dr. Lic. = ND/NA  
Lot No = 16512080A2  
Cyl No = 13  
Expiration Date = 08/01/2014  
County = 13                              Oper No. = 777777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Deb Kashur

Operator Signature  
DEB KASHUR

Remarks:

RFI test

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005955  
Location = DUNN      8164.13.00 06/09  
05/16/2014      09:48

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:49
02 Std. Gas	0.081	09:49
03 Room Air	0.000	09:50
04 Std. Gas	0.081	09:50
05 Room Air	0.000	09:51
06 Std. Gas	0.081	09:51
07 Room Air	0.000	09:52

Lot No = 16512080A2  
Cyl No = 13  
Exp Date = 08/01/2014  
County = 13      Oper No. = 777777

Deb Kasher  
Operator Signature  
DEB KASHUR

Remarks:

Calibration check 0.080 Ac

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005955  
Location = DUNN      8164.13.00 06/09  
05/16/2014      09:52

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:52
02 Std. Gas	0.082	09:53
03 Room Air	0.000	09:53
04 Std. Gas	0.082	09:54
05 Room Air	0.000	09:54
06 Std. Gas	0.081	09:54
07 Room Air	0.000	09:55

Lot No = 16512080A2  
Cyl No = 13  
Exp Date = 08/01/2014  
County = 13      Oper No. = 777777

*Deb Kasher*

Operator Signature  
DEB KASHUR

Remarks:

*Calibration check 0.080 AC*

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005955  
Location = DUNN      8164.13.00 06/09  
05/16/2014      09:56

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:56
02 Std. Gas	0.082	09:56
03 Room Air	0.000	09:57
04 Std. Gas	0.082	09:57
05 Room Air	0.000	09:58
06 Std. Gas	0.081	09:58
07 Room Air	0.000	09:59

Lot No = 16512080A2  
Cyl No = 13  
Exp Date = 08/01/2014  
County = 13      Oper No. = 777777

Deb Kasher  
Operator Signature  
DEB KASHUR

Remarks:

Calibration check 0.080 AE