

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-005944

Location: WAHP

- Check When Done
- A. Is the breath tube heated? ✓
 - B. Is the diagnostic check complete? ✓
 - C. Is the time, date, and year correct (re-set if necessary)? ✓
 - D. Print test. (Attach test record.) ✓
 - E. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)
Sim SN DR5115 Lot No. 12111 Rep AC 0.020 ✓
 - F. High AC. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN DR5143 Lot No. 13080 Rep AC ~~0.080~~ 0.300
CEE ✓
 - G. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.
(Attach test record.)
Sim SN DR3375 Lot No. ICS#2 Rep AC 0.100 AC + 0.05% ACETONE ✓
 - H. RFI check. Run CMS mode. Key radio on first room air.
(Attach test record.) ✓
 - I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.)

Lot No. 16512080A2 Cyl No. 058 Exp. Date 08/01/2014

Test 1 <u>0.082</u>	Test 4 <u>0.083</u>	Test 7 <u>0.082</u>
Test 2 <u>0.082</u>	Test 5 <u>0.084</u>	Test 8 <u>0.083</u>
Test 3 <u>0.082</u>	Test 6 <u>0.082</u>	Test 9 <u>0.083</u>

Average 0.083

J. Tank pressure, Level 3 Func D Sub G P: 635 psi Reg 650 psi

K. Remarks/Maintenance record: ANNUAL INSPECTION

14 MAY 14

Date

Charles E. Eden
Field Inspector's Signature

19 May 14 ✓

Date

Deb Kashur
Reviewed by

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005944
Location = WAHP 8164.13.00 06/09
05/14/2014 15:01

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 16512080A2
Standard Cyl #? 58
Standard Expiration? 08/01/2014
Oper No? 666666

Flow Cal. Date: 06/20/2013
Slope 677
Intercept -702620

IR Calibration Date: 06/20/2013

	3um	9um
0th Coef(*100):	-20012	-21455
1st Coef(*100):	275397	137470
2nd Coef(*100):	2699	1252
H2O adj(mg/l*10k):	704	563

***** Printer Test End *****



Operator Signature
CHARLES EDER

Remarks:

PRINT TEST

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005944
Location = WAHP 8164.13.00 06/09
05/14/2014 15:17

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:18
02 Std. Sol.	0.017	15:19
03 Room Air	0.000	15:19
04 Std. Sol.	0.019	15:20
05 Room Air	0.000	15:20
06 Std. Sol.	0.020	15:21
07 Room Air	0.000	15:22

08 Sim Temp = 34.0°C

Simul Ser No = DR5115
Std Sol No = 12111
County = 39 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *Low AC*
0.020 AC

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005944
Location = WAHP 8164.13.00 06/09
05/14/2014 15:11

Test	AC	Time
01 Room Air	0.000	15:12
02 *Subject Test	INT*	15:13
03 Room Air	0.000	15:13

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = INTERFERENT
Dr. Lic. = ND/DIS821456
Lot No = 16512080A2
Cyl No = 58
Expiration Date = 08/01/2014
County = 39 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: INTERFERENT CHECK
0.100 AC + 0.05% ACETONE

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005944
Location = WAHP 8164.13.00 06/09
05/14/2014 15:09

Test	AC	Time
01 Diagnostic	OK	15:10
02 Room Air	RFI*	15:10
03 Room Air	0.000	15:11

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = RFI CHECK
Dr. Lic. = ND/DIS821456
Lot No = 16512080A2
Cyl No = 58
Expiration Date = 08/01/2014
County = 39 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks: *RFI CHECK*


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005944
Location = WAHP 8164.13.00 06/09
05/14/2014 15:14

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:14
02 Std. Gas	0.082	15:15
03 Room Air	0.000	15:15
04 Std. Gas	0.082	15:15
05 Room Air	0.000	15:16
06 Std. Gas	0.082	15:16
07 Room Air	0.000	15:17

Lot No = 16512080A2
Cyl No = 58
Exp Date = 08/01/2014
County = 39 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005944
Location = WAHP 8164.13.00 06/09
05/14/2014 15:29

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:29
02 Std. Gas	0.083	15:29
03 Room Air	0.000	15:30
04 Std. Gas	0.084	15:30
05 Room Air	0.000	15:31
06 Std. Gas	0.082	15:31
07 Room Air	0.000	15:32

Lot No = 16512080A2
Cyl No = 58
Exp Date = 08/01/2014
County = 39 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080AL

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005944
Location = WAHP 8164.13.00 06/09
05/14/2014 15:32

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:33
02 Std. Gas	0.082	15:33
03 Room Air	0.000	15:33
04 Std. Gas	0.083	15:34
05 Room Air	0.000	15:34
06 Std. Gas	0.083	15:35
07 Room Air	0.000	15:35

Lot No = 16512080A2
Cyl No = 58
Exp Date = 08/01/2014
County = 39 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC

Form 106-I8000