



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59281 (03-2014)

Serial Number <i>80-006494</i>	Instrument Location <i>VALL</i>
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level 1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. RFI Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if Any): _____
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>Dan Kuntz</i>	Date <i>5/19/2016</i>
Reviewed By (Crime Laboratory Use Only) <i>Deb Shanaver</i>	Date <i>14 Oct 16</i>

This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Certified By: *Charles E. Egan* Certified Date: *14 Nov. 16*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006494
Location = VALL 8164.13.00 06/09
05/19/2016 14:51

***** Printer Test *****

abcdefghijklmnopqrstuvwxy1234567890-|=|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%&^&*()_+?

abcdefghijklmnopqrstuvwxy1234567890-|=|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%&^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? No
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 01316080A2
Standard Cyl #? 31
Standard Expiration? 02/05/2018
Oper No? 100012

Flow Cal. Date: 08/18/2015
Slope 655
Intercept -624238

IR Calibration Date: 08/18/2015
 3um 9um

0th Coef(*100):	-24880	-16926
1st Coef(*100):	256511	140602
2nd Coef(*100):	4127	1722
H2O adj(mg/l*10k):	783	426

***** Printer Test End *****



Operator Signature
DANA RUSTEBAKKE

Remarks:

Print Test


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006494
Location = VALL 8164.13.00 06/09
05/19/2016 14:51

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:52
02 Std. Gas	0.079	14:52
03 Room Air	0.000	14:53
04 Std. Gas	0.079	14:53
05 Room Air	0.000	14:54
06 Std. Gas	0.079	14:54
07 Room Air	0.000	14:55

Lot No = 01316080A2
Cyl No = 31
Exp Date = 02/05/2018
County = 02 Oper No. = 100012



Operator Signature
DANA RUSTEBAKKE

Remarks:

ACA Test

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006494
Location = VALL 8164.13.00 06/09
05/19/2016 14:58

Test	AC	Time
01 Diagnostic	OK	15:00
02 Room Air	0.000	15:01
03 *Subject Test	RFI*	15:01
04 Room Air	0.000	15:02

*Invalid Test
Inhibited - RFI

Sub Name = TEST, RFI NA
Sub DOB = 01/01/2000
Sub Sex = Male Weight = 175
Test = DUI Cit = 123456
Dr. Lic. = ND/123456
Lot No = 01316080A2
Cyl No = 31
Expiration Date = 02/05/2018
County = 02 Oper No. = 100012

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
DANA RUSTEBAKKE

Remarks:

RFI Test

Form 106-I8000