

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-005949

Location: TOXL

- Check When Done
- A. Is the breath tube heated? ✓
  - B. Is the diagnostic check complete? ✓
  - C. Is the time, date, and year correct (re-set if necessary)? ✓
  - D. Print test. (Attach test record.) ✓
  - E. Low AC. Use  $\leq 0.03$  AC in ACA mode. (Attach test record.)  
Sim SN MP3066 Lot No. 201507F Rep AC 0.020 ✓
  - F. High AC. Use  $\geq 0.25$  AC in ACA mode. (Attach test record.)  
Sim SN MP3057 Lot No. 16102 Rep AC 0.400 ✓
  - G. Interferent check. Use 0.05% acetone in ethanol solution in ABA mode.  
(Attach test record.)  
Sim SN DR5135 Lot No. ICS4 Solution 0.10AC + 0.05% ACETONE ✓
  - H. RFI check. Run CMS mode. Key radio during CMS test.  
(Attach test record.) ✓
  - I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.  
(Attach test record.) 0.080 AC

Lot No. 13314080A3 Cyl No. 034 Exp. Date 06/05/2016

Test 1	<u>0.078</u>	Test 4	<u>0.079</u>	Test 7	<u>0.078</u>
Test 2	<u>0.078</u>	Test 5	<u>0.077</u>	Test 8	<u>0.078</u>
Test 3	<u>0.077</u>	Test 6	<u>0.078</u>	Test 9	<u>0.077</u>

Average 0.078

J. Tank pressure, Level 3 Func D Sub G P: 1007 psi Reg ~1000 psi

K. Remarks/Maintenance record: ANNUAL INSPECTION, Breath hose & Simulator return O-rings replaced by DEBS.

5/17/2016  
Date

Charles E. Edr  
Field Inspector's Signature

19 May 16 ✓  
Date

Deb Shanaver  
Reviewed by

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-005949  
Location = TOXL                    8164.13.00 06/09  
05/17/2016    08:56

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&\*()\_+?

Current Instrument Setup


Data Entry Mode:                    Enabled  
Start Test Sequence:                DABACABA  
Display Prelim Rslt?                Yes  
Display Third Digit?                Yes  
Inhib Printer(Y/N)?                No  
Display Volume?                    No  
Disable On Memfull?                Yes  
# of Print Copies?                    1  
Select Std (D/W/I)?                Dry  
Standard Value?                    0.080  
Standard Lot #?                    13314080A3  
Standard Cyl #?                    34  
Standard Expiration?                06/05/2016  
Oper No?                            666666

Flow Cal. Date:                    06/20/2013  
Slope                                692  
Intercept                            -708987

IR Calibration Date:                04/06/2015  
    3um                                9um

-----  
0th Coef(\*100):                    -24726                                -28868  
1st Coef(\*100):                    263448                                133361  
2nd Coef(\*100):                    3398                                    1497  
H2O adj (mg/l\*10k):                326                                    367

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
CHARLES EDER

Remarks:

PRINT TEST

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005949  
Location = TOXL      8164.13.00 06/09  
05/17/2016      09:12

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:12
02 Std. Sol.	0.017	09:13
03 Room Air	0.000	09:14
04 Std. Sol.	0.017	09:15
05 Room Air	0.000	09:15
06 Std. Sol.	0.018	09:16
07 Room Air	0.000	09:16

08 Sim Temp = 34.0°C

Simul Ser No = MP3066  
Std Sol No = 201507F  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

Low AC  
0.020 AC

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005949  
Location = TOXL      8164.13.00 06/09  
05/17/2016      09:28

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:28
02 Std. Sol.	0.399	09:29
03 Room Air	0.000	09:30
04 Std. Sol.	0.399	09:30
05 Room Air	0.000	09:31
06 Std. Sol.	0.400	09:32
07 Room Air	0.000	09:32

08 Sim Temp = 34.0°C

Simul Ser No = MP3057  
Std Sol No = 16102  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:      HIGH AC  
                         0.400 AC

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005949  
Location = TOXL      8164.13.00 06/09  
05/17/2016      09:36

Test	AC	Time
01 Room Air	0.000	09:37
02 *Subject Test	INT*	09:37
03 Room Air	0.000	09:38

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = INTERFERENT  
Dr. Lic. = ND/DIS821456  
Lot No = 13314080A3  
Cyl No = 34  
Expiration Date = 06/05/2016  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks:

*Interferent Check*

*0.10 AC + 0.05% ACETONE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005949  
Location = TOXL      8164.13.00 06/09  
05/17/2016      09:39

Test	AC	Time
01 Diagnostic	OK	09:40
02 Room Air	0.000	09:40
03 *Subject Test	RFI*	09:41
04 Room Air	0.000	09:41

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 01/01/1982  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = RFI CHECK  
Dr. Lic. = ND/DIS821456  
Lot No = 13314080A3  
Cyl No = 34  
Expiration Date = 06/05/2016  
County = 08      Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
CHARLES EDER

Remarks:

*RFI CHECK*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005949  
Location = TOXL      8164.13.00 06/09  
05/17/2016      09:42

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:43
02 Std. Gas	0.078	09:43
03 Room Air	0.000	09:43
04 Std. Gas	0.078	09:44
05 Room Air	0.000	09:44
06 Std. Gas	0.077	09:45
07 Room Air	0.000	09:45

Lot No = 13314080A3  
Cyl No = 34  
Exp Date = 06/05/2016  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005949  
Location = TOXL      8164.13.00 06/09  
05/17/2016      09:46

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:46
02 Std. Gas	0.079	09:47
03 Room Air	0.000	09:47
04 Std. Gas	0.077	09:48
05 Room Air	0.000	09:48
06 Std. Gas	0.078	09:48
07 Room Air	0.000	09:49

Lot No = 13314080A3  
Cyl No = 34  
Exp Date = 06/05/2016  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC

Form 106-I8000



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005949  
Location = TOXL      8164.13.00 06/09  
05/17/2016      09:49

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:50
02 Std. Gas	0.078	09:50
03 Room Air	0.000	09:51
04 Std. Gas	0.078	09:51
05 Room Air	0.000	09:52
06 Std. Gas	0.077	09:52
07 Room Air	0.000	09:53

Lot No = 13314080A3  
Cyl No = 34  
Exp Date = 06/05/2016  
County = 08      Oper No. = 666666



Operator Signature  
CHARLES EDER

Remarks: CALIBRATION CHECK  
0.080 AC