



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59281 (5-2011)

Serial Number: <u>80-004949</u>	Instrument Location: <u>HILL</u>
Reason for Install/Repair:	
<input checked="" type="checkbox"/> Install After Receiving from Crime Laboratory	<input type="checkbox"/> Install After Location Change
<input type="checkbox"/> Other (Specify)	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location). (Level 2, Function E)
- 5. Scan/Enter Gas Cylinder Information. (Level 1, Function S)
- 6. Run Tests:
 - A. Print Test (Level 1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. RFI Test (Use CMS Mode – Key Radio During First Room Air).
- 7. Repair and/or Maintenance Performed (if Any): NA
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN 50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN 50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN 59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Troy Hoiker 14275
 Field Inspector Signature

5-20-11
 Date

Deb Kashner
 Reviewed By

25 May 2011
 Date

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004949
Location = HILL 8164.13.00 06/09
05/20/2011 19:01

***** Printer Test *****

abcdefghijklmnopqrstuvwxy1234567890-=
|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

abcdefghijklmnopqrstuvwxy1234567890-=
|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 15410080A1
Standard Cyl #? 15
Standard Expiration? 07/01/2012
Oper No? 049524

Flow Cal. Date: 09/24/2010
Slope 784
Intercept -862571

IR Calibration Date: 09/24/2010
 3um 9um

0th Coef(*100): -23764 -22121
1st Coef(*100): 267615 132337
2nd Coef(*100): 1069 992
H2O adj(mg/l*10k): 531 480

***** Printer Test End *****



Operator Signature
TROY A. HISCHER

Remarks:


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004949
Location = HILL 8164.13.00 06/09
05/20/2011 19:01

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	19:02
02 Std. Gas	0.078	19:02
03 Room Air	0.000	19:03
04 Std. Gas	0.078	19:03
05 Room Air	0.000	19:04
06 Std. Gas	0.078	19:04
07 Room Air	0.000	19:05

Lot No = 15410080A1
Cyl No = 15
Exp Date = 07/01/2012
County = 49 Oper No. = 049524



Operator Signature
TROY A. HISCHER

Remarks:

Form 106-I8000

80-004949
p 115 of 116
FM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004949
Location = HILL 8164.13.00 06/09
05/20/2011 19:05

Test	AC	Time
01 Diagnostic	OK	19:07
02 Room Air	RFI*	19:07
03 Room Air	RFI*	19:07

*Invalid Test
Inhibited - RFI

Sub Name = INSTAL, INSTALL I
Sub DOB = 01/01/1999
Sub Sex = Male Weight = 180
Test = DUI Cit = 222222
Dr. Lic. = ND/INT99999
Lot No = 15410080A1
Cyl No = 15
Expiration Date = 07/01/2012
County = 49 Oper No. = 049524

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
TROY A. HISCHER

Remarks:

Form 106-I8000

80-004949
P. 116 of 116
JAM