



INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59281 (5-2011)

Serial Number <i>80-004191</i>	Instrument Location <i>Bowman Co SO</i>
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
 - A. Print Test (Level 1, Function P).
 - B. ACA Test (Level 1, Function C).
 - C. RFI Test (Use CMS Mode; Key Radio During First Room Air).
- 7. Repair and/or Maintenance Performed (if Any): *N/A*
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>[Signature]</i>	Date <i>5-23-11</i>
Reviewed By (Crime Laboratory Use Only) <i>Deb Kashner</i>	Date <i>24 May 2011</i>

*80-004191
p. 87 of 90
JFM*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004191
Location = BOWM 8164.13.00 06/09
05/23/2011 11:08

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 659358
Standard Cyl #? 51
Standard Expiration? 09/03/2011
Oper No? 020814

Flow Cal. Date: 01/28/2010
Slope 738
Intercept -777574

IR Calibration Date: 01/28/2010
 3um 9um

0th Coef(*100): -23128 -33609
1st Coef(*100): 272485 131724
2nd Coef(*100): 1490 1056
H2O adj(mg/l*10k): 432 614

***** Printer Test End *****



Operator Signature
JOHNNY Z. LAWSON

Remarks:

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004191
Location = BOWM 8164.13.00 06/09
05/23/2011 11:24

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:25
02 Std. Gas	0.078	11:25
03 Room Air	0.000	11:26
04 Std. Gas	0.078	11:26
05 Room Air	0.000	11:27
06 Std. Gas	0.077	11:27
07 Room Air	0.000	11:28

Lot No = 659358
Cyl No = 51
Exp Date = 09/03/2011
County = 06 Oper No. = 020814



Operator Signature
JOHNNY Z. LAWSON

Remarks:

Form 106-I8000

80-004191
P. 89490
JZL


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004191
Location = BOWM 8164.13.00 06/09
05/23/2011 11:29

Test	AC	Time
01 Room Air	RFI*	11:31
02 Room Air	0.000	11:31

*Invalid Test
Inhibited - RFI

Sub Name = TEST, SUBJECT Q
Sub DOB = 09/09/1987
Sub Sex = Male Weight = 175
Test = OTH Cit = 090909
Dr. Lic. = ND/TES870909
Lot No = 659358
Cyl No = 51
Expiration Date = 09/03/2011
County = 06 Oper No. = 020814



Operator Signature
JOHNNY Z. LAWSON

Remarks:

Form 106-I8000

80-004191
P. 90 of 90
FM