

INTOXILYZER® 8000 INSPECTION

(WILL)

Serial No.: 80-004950

Location: WILLIAMS Co. LEC.

- Check When Done
- A. Is the breath tube heated? ✓
 - B. Is the diagnostic check complete? ✓
 - C. Is the time, date, and year correct (re-set if necessary)? ✓
 - D. Print test. (Attach test record.) ✓
 - E. Low AC. Use ≤ 0.03 AC in ACA mode. (Attach test record.)
Sim SN DR5143 Lot No. 201501B Rep AC 0.020 ✓
 - F. High AC. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN DR3455 Lot No. 14270 Rep AC 0.300 ✓
 - G. Interferent check. Use 0.05% acetone in ethanol solution in ABA mode.
(Attach test record.)
Sim SN DR5132 Lot No. ICS #2 Solution 0.10AC+0.05% ACE ✓
 - H. RFI check. Run CMS mode. Key radio during CMS test.
(Attach test record.) ✓
 - I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.) 0.080 AC

Lot No. 13314080A2 Cyl No. 51 Exp. Date 6/5/2016
Test 1 0.082 Test 4 0.081 Test 7 0.081
Test 2 0.082 Test 5 0.081 Test 8 0.082
Test 3 0.082 Test 6 0.080 Test 9 0.082
Average 0.081

J. Tank pressure, Level 3 Func D Sub G P: 468 psi Reg 475 psi

K. Remarks/Maintenance record: ANNUAL INSPECTION

5/28/15
Date

Charles E. Edm
Field Inspector's Signature

29 May 15 ✓
Date

Deb Kashner
Reviewed by

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004950
Location = WILL 8164.13.00 06/09
05/28/2015 09:42

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 13314080A2
Standard Cyl #? 51
Standard Expiration? 06/05/2016
Oper No? 666666

Flow Cal. Date: 01/18/2010
Slope 774
Intercept -908904

IR Calibration Date: 04/06/2011
 3um 9um

0th Coef(*100): -2573 -7967
1st Coef(*100): 257525 132746
2nd Coef(*100): 5623 2169
H2O adj(mg/l*10k): 540 597

***** Printer Test End *****



Operator Signature
CHARLES EDER

Remarks:

PRINT TEST

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004950
Location = WILL 8164.13.00 06/09
05/28/2015 10:03

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 10:04 |
| 02 Std. Sol. | 0.018 | 10:05 |
| 03 Room Air | 0.000 | 10:05 |
| 04 Std. Sol. | 0.018 | 10:06 |
| 05 Room Air | 0.000 | 10:07 |
| 06 Std. Sol. | 0.018 | 10:07 |
| 07 Room Air | 0.000 | 10:08 |


08 Sim Temp = 34.0°C

Simul Ser No = DR5143

Std Sol No = 201501B

County = 53

Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

Low AC

0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004950
Location = WILL 8164.13.00 06/09
05/28/2015 10:20

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 10:20 |
| 02 Std. Sol. | 0.304 | 10:21 |
| 03 Room Air | 0.000 | 10:22 |
| 04 Std. Sol. | 0.305 | 10:22 |
| 05 Room Air | 0.000 | 10:23 |
| 06 Std. Sol. | 0.308 | 10:23 |
| 07 Room Air | 0.000 | 10:24 |

08 Sim Temp = 34.0°C

Simul Ser No = DR3455
Std Sol No = 14270
County = 53 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: HIGH AC
0.300 AC


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004950
Location = WILL 8164.13.00 06/09
05/28/2015 09:54

| Test | AC | Time |
|------------------|-------|-------|
| 01 Room Air | 0.000 | 09:55 |
| 02 *Subject Test | INT* | 09:56 |
| 03 Room Air | 0.000 | 09:56 |

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = INTERFERENT
Dr. Lic. = ND/DIS821456
Lot No = 13314080A2
Cyl No = 51
Expiration Date = 06/05/2016
County = 53 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: *Interferent Check*
0.10AC + 0.05% ACETONE

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004950
Location = WILL 8164.13.00 06/09
05/28/2015 09:47

| Test | AC | Time |
|------------------|-------|-------|
| 01 Diagnostic | OK | 09:48 |
| 02 Room Air | 0.000 | 09:48 |
| 03 *Subject Test | RFI* | 09:49 |
| 04 Room Air | 0.000 | 09:49 |

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = RFI CHECK
Dr. Lic. = ND/DIS821456
Lot No = 13314080A2
Cyl No = 51
Expiration Date = 06/05/2016
County = 53 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
CHARLES EDER

Remarks:

RFI CHECK


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004950
Location = WILL 8164.13.00 06/09
05/28/2015 09:50

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 09:51 |
| 02 Std. Gas | 0.082 | 09:51 |
| 03 Room Air | 0.000 | 09:52 |
| 04 Std. Gas | 0.082 | 09:52 |
| 05 Room Air | 0.000 | 09:52 |
| 06 Std. Gas | 0.082 | 09:53 |
| 07 Room Air | 0.000 | 09:53 |

Lot No = 13314080A2
Cyl No = 51
Exp Date = 06/05/2016
County = 53 Oper No. = 666666


Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004950
Location = WILL 8164.13.00 06/09
05/28/2015 09:58

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 09:59 |
| 02 Std. Gas | 0.081 | 09:59 |
| 03 Room Air | 0.000 | 09:59 |
| 04 Std. Gas | 0.081 | 10:00 |
| 05 Room Air | 0.000 | 10:00 |
| 06 Std. Gas | 0.080 | 10:01 |
| 07 Room Air | 0.000 | 10:01 |

Lot No = 13314080A2
Cyl No = 51
Exp Date = 06/05/2016
County = 08 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks: CALIBRATION CHECK
0.080 AC

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004950
Location = WILL 8164.13.00 06/09
05/28/2015 10:26

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 10:27 |
| 02 Std. Gas | 0.081 | 10:27 |
| 03 Room Air | 0.000 | 10:28 |
| 04 Std. Gas | 0.082 | 10:28 |
| 05 Room Air | 0.000 | 10:28 |
| 06 Std. Gas | 0.082 | 10:29 |
| 07 Room Air | 0.000 | 10:29 |

Lot No = 13314080A2
Cyl No = 51
Exp Date = 06/05/2016
County = 53 Oper No. = 666666



Operator Signature
CHARLES EDER

Remarks:

CALIBRATION CHECK
0.080 AC

Form 106-I8000