



**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION  
 Toxicology Section/Breath Alcohol Program  
 SFN 59281 (03-2014)

Serial Number <u>80-004943</u>	Instrument Location <u>Mt. of Air Base</u>
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A.  Print Test (Level 1, Function P).
  - B.  ACA Test (Level 1, Function C).
  - C.  RFI Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- N/A 7. Repair and/or Maintenance Performed (if Any): \_\_\_\_\_
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- N/A 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <u>Deb Kasher</u>	Date <u>02 June 15</u> ✓
Reviewed By (Crime Laboratory Use Only) <u>Charles E. Edin</u>	Date <u>03 JUN 15</u>

This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Certified By: Charles E. Edin Certified Date: 23 Nov. 2016

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004943  
Location = MNAF                    8164.13.00 06/09  
06/02/2015                            11:20

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890- = |  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890- = |  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

Current Instrument Setup


Data Entry Mode:                    Enabled  
Start Test Sequence:                DABACABA  
Display Prelim Rslt?                Yes  
Display Third Digit?                Yes  
Inhib Printer(Y/N)?                 No  
Display Volume?                     No  
Disable On Memfull?                 No  
# of Print Copies?                  1  
Select Std (D/W/I)?                 Dry  
Standard Value?                     0.080  
Standard Lot #?                     13314080A2  
Standard Cyl #?                     53  
Standard Expiration?                06/05/2016  
Oper No?                             777777

Flow Cal. Date:                     01/18/2010  
Slope                                693  
Intercept                            -665387

IR Calibration Date:                05/21/2015  
    3um                                    9um


-----  
0th Coef(\*100):                    -28195                                -26668  
1st Coef(\*100):                    268915                                137995  
2nd Coef(\*100):                    103                                    992  
H2O adj(mg/l\*10k):                222                                    308

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
DEB KASHUR

Remarks:




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004943  
Location = MNAF      8164.13.00 06/09  
06/02/2015      11:32

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:33
02 Std. Gas	0.077	11:33
03 Room Air	0.000	11:33
04 Std. Gas	0.077	11:34
05 Room Air	0.000	11:34
06 Std. Gas	0.077	11:35
07 Room Air	0.000	11:35

Lot No = 13314080A2  
Cyl No = 53  
Exp Date = 06/05/2016  
County = 64      Oper No. = 77777

  
\_\_\_\_\_  
Operator Signature  
DEB KASHUR

Remarks:  
Calibration check 0.080AC

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004943  
Location = MNAF      8164.13.00 06/09  
06/02/2015      11:36

Test	AC	Time
01 Diagnostic	OK	11:37
02 Room Air	RFI*	11:37
03 Room Air	0.000	11:38

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female      Weight = 150  
Test = DUI      Cit = NA  
Dr. Lic. = ND/DIS921456  
Lot No = 13314080A2  
Cyl No = 53  
Expiration Date = 06/05/2016  
County = 64      Oper No. = 777777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
DEB KASHUR

Remarks:

*RFI check*

Form 106-18000