

INTOXILYZER 8000 INSPECTION

Serial No.: 80-003059

Location: West Fargo PD

- | | Check When Done |
|--|-----------------|
| A. Is the warm-up time less than 20 minutes? | <u>✓</u> |
| B. Is a three-pronged grounded outlet used? | <u>✓</u> |
| C. Is the breath tube heated? | <u>✓</u> |
| D. Is the diagnostic check complete? | <u>✓</u> |
| E. Is the time, date, and year correct (re-set if necessary)? | <u>✓</u> |
| F. Print test. (Attach test record.) | <u>✓</u> |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.) Sim SN <u>DR5132</u> Lot No. <u>—</u> Rep AC <u>0.021 AC</u> | <u>✓</u> |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.) Sim SN <u>DR3375</u> Lot No. <u>—</u> Rep AC <u>0.292 AC</u> | <u>✓</u> |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode. (Attach test record.) Sim SN <u>DR3454</u> Lot No. <u>—</u> Rep AC <u>0.10 AC + Acetone</u> | <u>✓</u> |
| J. RFI check. Run CMS mode. Key radio on first room air. (Attach test record.) | <u>✓</u> |
| K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode. (Attach test record.) | |

| | | |
|-------------------------|---------------------|-----------------------------|
| Lot No. <u>659 357</u> | Cyl No. <u>10</u> | Exp. Date <u>09/03/2011</u> |
| Test 1 <u>0.080</u> | Test 4 <u>0.079</u> | Test 7 <u>0.081</u> |
| Test 2 <u>0.081</u> | Test 5 <u>0.080</u> | Test 8 <u>0.080</u> |
| Test 3 <u>0.080</u> | Test 6 <u>0.081</u> | Test 9 <u>0.081</u> |
| Average <u>0.080 AC</u> | | |

L. Tank pressure, Level 3 Func D Sub G P: 1001 psi Reg 1000 psi

M. Remarks/Maintenance record: Annual Inspection

30 June 2010
Date

Margy Plann
Field Inspector's Signature

14 June 2010
Date

Deb Kashur
Reviewed by

80-003059
P. 84 of 92
Jm

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003059
Location = WFPD 8164.13.00 06/09
06/03/2010 09:32

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 659357
Standard Cyl #? 10
Standard Expiration? 09/03/2011
Oper No? 107501

Flow Cal. Date: 04/07/2008
Slope 713
Intercept -692349

IR Calibration Date: 06/04/2008

| | 3um | 9um |
|--------------------|--------|--------|
| 0th Coef(*100): | -15808 | -12583 |
| 1st Coef(*100): | 279401 | 139580 |
| 2nd Coef(*100): | 1046 | 760 |
| H2O adj(mg/l*10k): | 478 | 208 |

***** Printer Test End *****



Operator Signature
MARGY PEARSON

Remarks:

Print

Form 106-I8000

80-003059
P. 85 of 92
JFM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003059
Location = WFPD 8164.13.00 06/09
06/03/2010 09:38

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 09:38 |
| 02 Std. Sol. | 0.019 | 09:39 |
| 03 Room Air | 0.000 | 09:40 |
| 04 Std. Sol. | 0.019 | 09:40 |
| 05 Room Air | 0.000 | 09:41 |
| 06 Std. Sol. | 0.019 | 09:42 |
| 07 Room Air | 0.000 | 09:42 |


08 Sim Temp = 34.0°C

Simul Ser No = 0000000

Std Sol No = 000000

County = 09

Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Low AC 0.021 AC ACA

Form 106-I8000

80-003059
P. 86 of 92
JEM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003059
Location = WFPD 8164.13.00 06/09
06/03/2010 09:43

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 09:44 |
| 02 Std. Sol. | 0.296 | 09:45 |
| 03 Room Air | 0.000 | 09:45 |
| 04 Std. Sol. | 0.295 | 09:46 |
| 05 Room Air | 0.000 | 09:47 |
| 06 Std. Sol. | 0.293 | 09:47 |
| 07 Room Air | 0.000 | 09:48 |


08 Sim Temp = 34.0°C

Simul Ser No = 0000000

Std Sol No = 000000

County = 09

Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Linearity 0.292 AC

Form 106-18000

80-003059
P. 87 of 92
FM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003059
Location = WFPD 8164.13.00 06/09
06/03/2010 09:48

| Test | AC | Time |
|------------------|-------|-------|
| 01 Room Air | 0.000 | 09:49 |
| 02 *Subject Test | INT* | 09:50 |
| 03 Room Air | 0.000 | 09:50 |

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 11/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Simul Ser No = 0000000
Std Sol No = 000000
County = 09 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Interferent Alc + Acetone

Form 106-I8000

80-003059
P. 88 of 92
RM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003059
Location = WFPD 8164.13.00 06/09
06/03/2010 09:55

| Test | AC | Time |
|---------------|-------|-------|
| 01 Diagnostic | OK | 09:56 |
| 02 Room Air | RFI* | 09:56 |
| 03 Room Air | 0.000 | 09:57 |

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 11/01/1982
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS821456
Lot No = 659357
Cyl No = 10
Expiration Date = 12/23/2011
County = 09 Oper No. = 107501

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
MARGY PEARSON

Remarks:

RFI

Form 106-I8000

80-003059
P. 87 of 92
PM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003059
Location = WFPD 8164.13.00 06/09
06/03/2010 09:57

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 09:58 |
| 02 Std. Gas | 0.080 | 09:58 |
| 03 Room Air | 0.000 | 09:59 |
| 04 Std. Gas | 0.081 | 09:59 |
| 05 Room Air | 0.000 | 09:59 |
| 06 Std. Gas | 0.080 | 10:00 |
| 07 Room Air | 0.000 | 10:00 |

Lot No = 659357
Cyl No = 10
Exp Date = 12/23/2011
County = 09 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Cal Check 0.080 AC

Form 106-I8000

80-003059
P. 90 of 92
JSM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003059
Location = WFPD 8164.13.00 06/09
06/03/2010 10:01

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 10:02 |
| 02 Std. Gas | 0.079 | 10:02 |
| 03 Room Air | 0.000 | 10:02 |
| 04 Std. Gas | 0.080 | 10:03 |
| 05 Room Air | 0.000 | 10:03 |
| 06 Std. Gas | 0.081 | 10:04 |
| 07 Room Air | 0.000 | 10:04 |

Lot No = 659357
Cyl No = 10
Exp Date = 12/23/2011
County = 09 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Cal Check 0.080 AC

Form 106-I8000

80-003059
P. 9/1/92
JFM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003059
Location = WFPD 8164.13.00 06/09
06/03/2010 10:05

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 10:05 |
| 02 Std. Gas | 0.081 | 10:06 |
| 03 Room Air | 0.000 | 10:06 |
| 04 Std. Gas | 0.080 | 10:07 |
| 05 Room Air | 0.000 | 10:07 |
| 06 Std. Gas | 0.081 | 10:08 |
| 07 Room Air | 0.000 | 10:08 |

Lot No = 659357
Cyl No = 10
Exp Date = 12/23/2011
County = 09 Oper No. = 107501



Operator Signature
MARGY PEARSON

Remarks:

Cal Check 0.080 AC

Form 106-I8000

80-003059
P. 92 of 92
VM